# ImageMaster 2D Platinum 7.0

User Manual





# Contents

1	Get	ting started	
	1.1	About this User Manual	9
	1.2	About the software	9
	1.3	System requirements	9
	1.4	Install the software	10
	1.5	eLicensing	10
		1.5.1 Access code	
		1.5.2 Find the physical address	
		1.5.3 Install the GE Healthcare eLicense server	
		1.5.4 Collect and place an eLicense file	
	1.6	Launch the software	
	1.7	Customer support	
		·	
	1.8	What's new?	14
2	Gra	phical user interface	
	2.1	About the interface	17
	2.2	Menu bar	17
		2.2.1 Keyboard shortcuts	18
	2.3	Toolbars	
		2.3.1 Default toolbars	
		2.3.2 Customize toolbars	
	2.4	Status bar	
	2.5	Display zone	
		2.5.1 Sheets	
		2.5.2 Panes	
		2.5.4 Switch order	
	2.6	Workspace	
		2.6.1 Toolbar	
		2.6.2 Navigator	
	2.7	Reports	28
		2.7.1 Dynamic content	
		2.7.2 Content based on enabled spots	
		2.7.3 Toolbars	
		2.7.4 Customize reports	
	2.0	2.7.5 Edit table cells	
	2.8	Dockable windows	
	2.9	Options	34
3	Imo	ge Pool	

	3.1	Introduction	37
	3.2	Start with good images	37
		3.2.1 Resolution	
		3.2.2 Image depth	37
		3.2.3 Calibration	
		3.2.4 Image editing	38
		3.2.5 File format	38
		3.2.6 DIGE file naming convention	38
	3.3	Preview images	
		3.3.1 Open images	
		3.3.2 Image pool sheet	39
	3.4	Process images	41
		3.4.1 Rotate	
		3.4.2 Flip	
		3.4.3 Crop	
		3.4.4 Invert gray levels	
	3.5	Calibrate images	
		3.5.1 Display calibration information about images	
		3.5.2 Prepare for a calibration of the image capture device	
		3.5.3 Create calibration of the capture device	
		3.5.4 Apply a calibration to an image	
		3.5.5 Remove a calibration	
		3.5.6 Control a calibration	48
	3.6	Describe images	
		3.6.1 Gel descriptions	
		3.6.2 Display gel descriptions	
		3.6.3 Add or edit gel descriptions	
		3.6.4 Permanent gel description categories	51
4	Pro	jects	
	4.1	Introduction	53
		4.1.1 Match hierarchies	53
		4.1.2 Classes	55
	4.2	Create a project	56
		4.2.1 The first time the software is launched	56
		4.2.2 At any time	57
		4.2.3 Add files to project	57
	4.3	Create match hierarchy	58
		4.3.1 Create a match set	58
		4.3.2 Merge a match set	
		4.3.3 Set reference	
		4.3.4 Use existing match sets	
		4.3.5 Export / import a match set	
	4.4	Create classes	
		4.4.1 Create a class	63
	4.5	Handle project items	63

		4.5.1 Display	63
		4.5.2 Remove	63
		4.5.3 Properties	64
		4.5.4 Move	
	4.6	Save projects	
		4.6.1 Project folder	
		4.6.2 Save	
		4.6.3 Share project	
	4.7	Manage projects	
		4.7.1 Add / remove projects	
		4.7.2 Add files to project	
		4.7.3 Backup / restore project	
		4.7.4 Project visibility	
		4.7.5 Project properties	66
5	Gel	S	
	5.1	Introduction	67
	5.2	Manipulate images	
		5.2.1 Move	
		5.2.2 Zoom	
		5.2.3 Same position and zoom factor	
		5.2.4 Region	
		5.2.5 Measure	
		5.2.6 Bookmarks	
	5.3	View signal intensity	
		5.3.1 Adjust contrast	
		5.3.2 3D view	
		5.3.3 Profile	
		5.3.4 Cursor information	
	5.4	Visually compare images	
		5.4.1 Sheet reference	
		5.4.2 Purpose of aligning images 5.4.3 Alian images	
		5.4.3 Align images         5.4.4 Show dual color	
		5.4.5 Spots overlapped	
		• • • • • • • • • • • • • • • • • • • •	
	5.5	Grid lines5.5.1 Display grid lines	
		5.5.2 Edit grid lines	
	5.6	Gel reports	
	5.0	5.6.1 Gel table	
	5.7	Save, export and print images	
		5.7.1 Save	
		5.7.2 Export	
		5.7.3 Print images and sheets	

# 6 Spots

	6.1	Introduction	87		
	6.2	Detect spots in non-DIGE gels	87		
		6.2.1 Procedure	87		
		6.2.2 Spot detection parameters	88		
		6.2.3 Spot quantification			
	6.3	Co-detect spots in DIGE gels	91		
		6.3.1 Procedure			
		6.3.2 Spot detection parameter			
		6.3.3 Spot quantification			
		6.3.4 Exclude spots	93		
	6.4	Select spots			
		6.4.1 Select			
		6.4.2 Spot sets			
		6.4.3 Enabled spots			
	6.5	Display spots			
		6.5.1 Spot shape			
		6.5.2 Spot color	97		
	6.6	Edit spots	97		
		6.6.1 Manual editing			
		6.6.2 Composite spots			
		6.6.3 Propagate spots	99		
	6.7	MW and pl calibration	100		
	6.8	Spot reports	102		
		6.8.1 Spot Table			
		6.8.2 Cursor information	102		
7	Mat	Matches			
	7.1	Introduction	105		
	7.2	Display a match hierarchy			
	7.3	Define landmarks			
	7.5	7.3.1 Rules for defining landmarks			
		7.3.2 Define Landmark			
	7.4	Automatic matching			
		<del>-</del>			
	7.5	Select matches			
		7.5.1 Select			
	7.0				
	7.6	Display matches			
		7.6.1 Show vectors			
		7.6.2 Show ID			
	7.7	Edit matches			
	7.8	Match reports			
		7.8.1 Match statistics table	114		
8	Dat	ta analysis			
	8.1	Introduction	117		

	8.2	Gene	ral settings	117
		8.2.1	Quantification value	117
		8.2.2	Statistics	117
	8.3	Analy	/ze gels	120
		8.3.1	Scatter plots	120
		8.3.2	Gel analysis table	123
		8.3.3	Gel analysis histograms	124
		8.3.4		
		8.3.5	DIGE histogram	133
	8.4	Analy	yze classes	135
		8.4.1	Specify classes	135
		8.4.2	Overlapping measures	135
		8.4.3	Statistical tests	
		8.4.4	Class analysis table	141
		8.4.5	Class analysis histograms	143
9	Anr	otatio	ons	
	9.1	Intro	duction	147
	9.2	Creat	te annotations and labels	148
		9.2.1	Create annotations	148
		9.2.2	Add labels to existing annotations	148
		9.2.3	Link annotations to spots	148
	9.3	Creat	te label categories	149
		9.3.1	Predefined label categories	
		9.3.2	Create new categories	
	9.4	Conn	ect to protein databases	151
		9.4.1	Introduction	151
		9.4.2	Set the database	152
		9.4.3	Query the database	
		9.4.4	Connext to an executable	154
	9.5	Creat	te specific links	
		9.5.1	Http link:	154
		9.5.2	File link:	155
		9.5.3	Text link:	
	9.6		t annotations and labels	
			Select	
			Select menu	
		9.6.3	Reports	
	9.7		ay annotations and labels	
			Annotation flag position	
			Visibility of annotations and labels	
	9.8		tation table	
	9.9		annotations and labels	
		9.9.1	Select	
		9.9.2		
		9.9.3	Annotation table	163

		9.9.4 Import annotations	164
10	Dat	a integration	
	10.1	Introduction	165
	10.2	Convert projects from earlier software versions	165
	10.3	Acquire images from Twain compatible scanners	165
		10.3.1 Select source	
		10.3.2 Acquire	166
	10.4	Export data	
		10.4.1 Reports	
		10.4.2 XML	
		10.4.3 Gel and report identifiers	167
	10.5	Export to spot excision robots	167
		10.5.1 GE Healthcare Ettan spot picker	
		10.5.2 Spot pickers from other manufacturers	
11	Und	do, redo and history	
	11.1	Undo, redo	171
		History	

# Getting started

#### 1.1 **About this User Manual**

All users must read this entire manual to fully understand the safe use of ImageMaster™ 2D Platinum 7.0. You will find detailed explanations of all the features and functionalities. The chapters in this manual are generally organized according to the logical sequence of a 2-DE gel analysis, although expert users will agree that some steps can be inverted or repeated at some point.

The user manual contains notes, as defined below.

Note: A Note is used to indicate information that is important for trouble-free and optimal use of the product.

### About the software 1.2

ImageMaster offers a flexible solution for the comprehensive visualization, exploration and analysis of 2-D gel data.

This version of ImageMaster was developed by a team from the Swiss Institute of Bioinformatics (SIB) in collaboration with Geneva Bioinformatics (GeneBio) SA and GE Healthcare.

There are two modules of ImageMaster 2D Platinum 7.0 available for purchase:

- ImageMaster 2D Platinum 7.0 **DIGE**: To be used with conventional 2-DE and DIGE gels. It is fully functional, enabling you to add and import DIGE gels directly in the workspace. You can also co-detect DIGE gels using the algorithm created by the GE Healthcare DeCyder™ software development team as well as match, report, plot histograms and perform statistical analyses on DIGE gels.
- ImageMaster 2D Platinum 7.0: To be used with conventional 2-DE gels. All menu commands related to DIGE are not functional and are grayed out in the graphical user interface.

## 1.3 System requirements

In order to install and run ImageMaster, your computer must satisfy the following requirements:

- Microsoft Windows XP or Vista operating systems.
- Administrative permission to install ImageMaster.

# 1 Getting started

# 1.4 Install the software

- At least 500 MB of RAM for ImageMaster (recommended: Intel dual-core
  processor with 1GB of RAM) and 768 MB RAM for ImageMaster DIGE
  (recommended: Intel dual-core processor with 2GB of RAM). The amount of
  memory required is determined by the number and size of image files to
  be processed simultaneously. Increased memory therefore enhances the
  performance when many and/or large images are analyzed.
- A high-quality display. To take full advantage of the software including the 3D View feature, the color resolution should be set to 24 bit (16.7 million colors). However, a color resolution of 8 bit (256 colors) is generally sufficient. It is recommended to use a screen resolution of at least 1024 x 768 pixels.
- At least 60 MB of available disk space for program files.
- Internet Explorer 6.0 (Microsoft Corporation), Netscape Navigator 7.0
  (Netscape Communications Corporation), Mozilla 1.4 (Mozilla Foundation)
  or higher versions. A browser allows you to print reports and to access
  scientific databases on the Web.

# 1.4 Install the software

You can install ImageMaster from a CD-ROM or by downloading the installation package over the Internet. When the CD-ROM is inserted into the appropriate drive on your computer, the Setup Wizard starts automatically and gives a series of on-screen instructions. Alternatively, you can double-click on the icon of the installer file (.msi or .exe file) to launch the Setup Wizard.

The ImageMaster installer creates a default directory on your hard disk called Program Files\GE Healthcare\ImageMaster 2D Platinum 7.0, in which the program files are placed. If you want to save the default directory in a different folder, then browse and select the location before continuing the installation.

Once installation is complete, it is recommended to restart your computer.

# 1.5 eLicensing

An electronic license (eLicense) file is required to enable ImageMaster to run once installation is completed.

There are two types of eLicenses:

Node-locked license (Machine license): To be used on a single computer.
 It is practical when only a few computers are used for working with ImageMaster. A node-locked license file must be placed on the computer running ImageMaster.

Floating license (Concurrent license): To be used by all computers networked to a license server. It is useful when many users, but not all at the same time, need access to the software. The number of computers that can simultaneously work with ImageMaster depends on the license, and is administered via the license server. A floating license file must be placed on the computer running the eLicense server. This server can either be installed on a computer running ImageMasteror on any other network computer (recommended: install on a network computer that is continually running).

ImageMaster checks that a valid license file is available each time the software is launched

#### 1.5.1 Access code

After ordering ImageMaster, a letter including an access code will be sent to the order's shipment address. The access code is necessary for collecting the eLicense files. Store this access code in a safe place.

### 1.5.2 Find the physical address

You must know the physical address of your computer when collecting the eLicense file. The address identifies the computer and is used by the licensing system. If purchasing a node-locked license, you need the physical address of the computer where ImageMaster is to be installed. If purchasing a floating license is used, you need the physical address of the computer where the eLicense server is to be installed

# To find the physical address:

- On the computer where the eLicense file is to be placed, choose (AII) **Programs: Accessories:: Command Prompt** in the Windows Start menu.
- In the Command Prompt window, enter the *ipconfig /all*. There must be a space between ipconfig and /all.
- Note down the **Physical Address** among the displayed information but ignore any dashes or colons. If several physical addresses are listed (dependent on the computer's network connections), any of them can be used for identification purposes. The physical address must be correctly noted for the license file to work.


#### Install the GE Healthcare eLicense server 1.5.3

This section only applies to floating licenses. If you have a node-locked license, then skip to the next section.

To install the GE Healthcare eLicense server:

# 1 Getting started1.5 eLicensing

- 1 On the computer where the eLicense server is to be installed, insert the ImageMaster 2D Platinum 7.0 installation CD-ROM.
- 2 Click the *Install GE Healthcare eLicense Server* button.
- 3 In the installation window, click **Next**.
- 4 Accept the default installation path and start the installation. By changing the default installation path, the file paths in LMTOOLS will need to be updated accordingly.
- 5 Answer **Yes** to any question about Windows Firewall (this may or may not appear) in order for the eLicense server to work properly.
- 6 Verify that the Launch GE Healthcare Software Licensing Server box is checked and click Finish. LMTOOLS automatically opens once the installation is successfully completed.
- 7 Click the **Config Services** tab and verify that all options are set as in Figure 1-1. The file paths shown (C:\Program Files\GE Healthcare\eLicense server\) are valid if the default installation was done.



**Figure 1-1.** The Config Services tab in LMTOOLS.

8 Click the **Save Service** button, even if no changes were made.

# 1.5.4 Collect and place an eLicense file

To collect and place a license file:

- Go to the web site http://www.elicensing.amershambiosciences.com/ gtlweb/login.
- 2 Enter the access code and click Continue.
- 3 Click Collect Licenses.
- Select the product (the box is checked). 4
- 5 Register your product.
- Enter the *Physical Address* (or Host-ID) of the computer where the license file is to be placed (by default the type of Host-ID is **Ethernet**). Verify that the Physical Address is correct. Click Continue.
- If the data shown is correct, then select **Collect License**. Otherwise navigate back in your Web browser to edit.
- Save the license to a file, or email the license. It is highly recommended to save the license file. If emailing, the license is delivered as raw text in the body of the message. Copy and paste the text into Notepad and save as an ASCII file with the extension .lic.
- Place the license file. For a node-locked license, the eLicense file must be placed in the ImageMaster installation directory (C:\Program Files\GE Healthcare\ImageMaster 2D Platinum, by default). For a floating license, the file must be placed in the C:\Program Files\GE Healthcare\eLicense server\Licenses on the computer running the eLicense server.
- 10 Restart the computer.

#### 1.5.5 Test the eLicense

In order to test that the license file is correctly placed, start ImageMaster on all computers where the software is installed.

For a floating license, a FLEXIm License Finder window will display asking to specify the *License Server System* or *License File*. Choose the first option and then enter the name of the computer on where the license server is installed. Click OK.

#### Launch the software 1.6

Double-click the ImageMaster icon (Figure 1-2) found on your computer's desktop to start the software. Alternatively, go to (All) Programs: ImageMaster 2D Platinum 7.0 in the Windows Start menu. A splash page appears while the software is loading. It will disappear automatically.



# 1.7 Customer support

Figure 1-2. ImageMaster icon.

# 1.7 Customer support

GE Healthcare provides technical and scientific support for ImageMaster. Please contact us if any problems arise during the installation or use of the software. Our support team is happy to help you.

# How to contact us

By email: imagemaster-support@ge.com.

You can also find a local office at www.gelifesciences.com/contact.

# Product, system and license information

Launch the software and choose *Help: About* ImageMaster to obtain product, system and license information. The product name, version number and version date that are displayed will be requested in technical support issues. You can view system and license information by clicking on the corresponding buttons.

The *Copy to clipboard* button in the License Information window can be used to paste license data into an email. In the System Information window, find facts about your computer or choose *File: Export* to save all the information in a text file.

# 1.8 What's new?

The graphical interface and user interaction modes have been entirely updated. The streamlined analysis in ImageMaster now offers enhanced usability and speed.

Features have been redesigned to help minimize manual spot editing and repetitive match editing. Ultimately, tasks are simplified and the reliability of the results are increased.

The most important new features in ImageMaster, version 7.0 are listed below:

- Fully dynamic tables, histograms, plots and 3-D views in which both content and selection are continuously updated to stay synchronized with the corresponding sheet that contains the gel images.
- Simplified import and visualization of images
- Improved population matching
- Management of multiple matches and composite spots
- Reorganized menu structure with icons
- Custom and context-related toolbars

- One-click to choose desired layout of sheets and panes
- Option to surround spot selections by boxes, for easier identification
- Single tool to select/edit spots and annotations
- Dedicated landmark tool
- Measure tool to compute pixel, pl, MW, or real-world (centimeter or inch) distances between spots
- Reviewed contrast adjustment feature
- New 3D View
- Customized report templates
- Adaptive display of histograms
- And much more...

# 1 Getting started 1.8 What's new?

# Graphical user interface

### About the interface 2.1

The graphical user interface is divided into four main parts, shown in Figure 2-1. They are the Menu Bar, the Toolbars, the Status Bar, and the Display Zone.

The Display Zone is the center of the interface. This is where gel images are arranged in sheets and panes. The Workspace and any reports are dockable windows found along the edges of the Display Zone.




Figure 2-1. The ImageMaster window. (a) Menu Bar, (b) Toolbar, (c) Status Bar, (d) Display Zone, (e) Sheet, (f) Pane, (g) Workspace, and (h) Reports.

#### 22 Menu bar

You can choose actions to be performed during your analysis from the Menu Bar. The menus are context related. This means some of the commands may not be available all of the time and either go away or are grayed out. For example,

the Select menu allows you to select spots, matches and annotations for detection and matching. Therefore, this menu is neither available nor displayed while viewing and editing images in the Image Pool.

Menu	Description
File	Close, save, import, export, print and other basic operations. You can also exit ImageMaster.
Edit	Undo/redo the last operations, show a history of operations, or edit (add, modify, delete) specific gels, spots, annotations or matches. You can also edit spot sets and enable spots.
View	Modify the settings for grid lines, profile or overview in the display, align images, show dual color or spot overlap in the current sheet, or change the way gels, spots, annotations or matches are visualized.
Select	Select specific spots, spot sets, annotations or matches.
Reports	Display tabular or graphical information about selected gels, spots, annotations or matches, and compute differences and similarities between gel images. The data analysis is based on robust statistics, factor analysis, and statistical tests.
Tools	Change display, quantification, and other options at the software level and customize the user interface (custom toolbars, keyboard shortcuts). Create and control a calibration tablet while working in the Image Pool .
Help	Access documentation and obtain license, product, version, or system information.

Most menu commands also have toolbar icons and/or keyboard shortcuts.

# 2.2.1 Keyboard shortcuts

Keyboard shortcuts, when available, are designated on the right-hand side of the corresponding menu command. A list of all shortcuts is given in the Appendix. Please note the logic behind the key combinations:

Ctrl is used to maneuver gels.

**Shift** is used to maneuver spots.

Alt is used to maneuver annotations.

Ctrl + Shift is used to maneuver matches.

You can create your own keyboard shortcuts for menu commands.

**Note:** A new keyboard shortcut must be unique. Be careful not to duplicate a keyboard shortcut that is already assigned to another menu command.

You risk destroying an existing shortcut. However, a given command can have several, different keyboard shortcuts.

# To add your own keyboard shortcut:

- Choose Tools: Customize.
- 2 Click the **Keyboard** tab.
- Select the menu in the *Category* drop-down list and then select the 3 command for which you want to create a shortcut (Figure 2-2).
- Read through the *Key assignments* list to see if the shortcut is currently assigned.
- Click in the **Press new shortcut key** box and then press the new combination on your keyboard.
- Click **Assign**. The software warns you when the shortcut is already assigned to another command and asks if you want to re-assign it. Click **No** and come up with a different shortcut.
- Click Close.




Figure 2-2. Defining keyboard shortcuts in the Customize window.

# 2.3 Toolbars

# 2.3.1 Default toolbars

The standard toolbars provided are designed so that you can quickly manipulate your gel data and apply the most frequently used features. Each of the tools are described in detail elsewhere in the User Manual.



This toolbar contains the basic tools of the software. More importantly, their functions (Move, Zoom, Region, Selection, Measure and Landmark) are not available in any menu. Although they do have corresponding keyboard shortcuts.



This toolbar contains the options Undo Zoom/Move, Redo Zoom/Move, Adjust Contrast and Cursor Information Window.



The image toolbar (Rotate, Flip, Invert Gray Levels, Crop and Add Files to Project) is contextual; it is only available when working with images in the Image Pool.

# Detect And Match Spots

This toolbar is only displayed when working with gels in a Match or Class sheet (that is, once gels have been imported into a project). You can Detect, Enable Edit, Match Gels, Add Match, Delete Match and Show Vectors.



This toolbar is only available when spot edition is enabled. Choose *Edit: Spots: Edit Enabled* to get into this mode.

#### 2.3.2 **Customize toolbars**

Toolbars can be configured according to your individual specifications.

# Toolbar position

To change the position of a toolbar, click the left edge and drag the toolbar to the position you want. You can drag a toolbar to any of the edges of the user interface. When a toolbar is dragged outside of the frame of the ImageMaster window, it becomes a floating window.

# Toolbar format

By default, the toolbar icons in the software are small (16x16 pixels). You can choose to display large icons (24x24 pixels) for better visibility. In addition, when you move the mouse over an icon, a screentip appears. Perhaps you prefer to hide all the screentips, or remove the keyboard shortcuts for the tools from the screentips.

# To modify the toolbar format:

- Choose **Tools**: **Customize**
- Click the **Options** tab in the Customize window.
- Select the desired options in the *Other* section.
- Click Close. 4

# **Custom toolbars**

You can create your own toolbars with icons for the functions you use most commonly.

# To create a custom toolbar:

- Choose Tools: Customize.
- Click the **Toolbars** tab in the Customize window. 2
- Click **New** to create a toolbar. Enter a name for the toolbar and click **OK**. An empty toolbar is created below the existing ones.
- Click the **Commands** tab in the Customize window.
- Select the menu from the *Category* list. Then select the command for which 5 you want to create an icon and drag it to the empty toolbar.
- Repeat step 5 to add icons to the toolbar.
- 7 Click Close in the Customize window.


# 2.4 Status bar

The Status Bar at the bottom of the Display Zone is an important resource. It indicates the total number of gels, spots, matches and annotations that are selected in the current sheet.

If you move your mouse over a gel image, the Status Bar also indicates the X and Y coordinates at the cursor position, as well as the image intensity. The unit of the coordinates can be changed in *Tools: Options*, under the *Display* tab.

# 2.5 Display zone

Gel images are opened in the Workspace and viewed in sheets and panes in the Display Zone (Figure 2-3). Their layouts can be arranged according to your requirements.



**Figure 2-3.** The Display Zone. There are two sheets open here, including the Image Pool. The current sheet AB is in front and contains four panes (AT1, AT2, BT1 and BT2). Only the pane AT1 is selected (green tab). By clicking the Layout icon for the pane, different options (Stacked, Tiled, One Row, One Column, Free) are displayed. Similarly, click the Layout icon in the upper right corner to re-arrange a sheet.

# 2.5.1 Sheets

When no dockable windows (Workspace, reports) are open, sheets occupy the entire Display Zone. Each sheet has a tab, with its name and an icon representing its type:

• **[Parage Pool**: This sheet contains images for viewing and basic processing.

- MatchSet: This sheet is opened by right-clicking on the name of a match set in the Workspace. Spot detection and matching must be carried out on this type of sheet.
- Class: This sheet is opened by right-clicking on the name of one or several classes in the Workspace. To carry out advanced expression analysis, you must work in this type of sheet.

When you move the mouse cursor over a sheet tab, the screentip specifies the type.

## Selection

The current sheet is always in front and its name is in bold. Select a sheet by clicking on its tab.

## Close

Click the *Close* icon in the upper right corner of the current sheet when you are finished working with it.

## Layout

The **Layout** icon to the left of the **Close** icon can be used to choose the arrangement of the panes in a sheet.

Option	Description
Stacked	Panes in a sheet are one on top of another.
Tiled	Panes in a sheet are side by side.
One Row	Panes in a sheet are in a single horizontal line.
One Column	Panes in a sheet are in a single vertical line.
Free	You specify the number of panes laid out horizontally and vertically in a sheet.

### 2.5.2 **Panes**

There can be one or more panes in a sheet (see Figure 2-3). Each pane has a tab with its name. On the left side of the tab, one or more icons describe the match or class hierarchy.

## Selection

Select a pane by clicking on its tab. Use the Shift or Ctrl keys to make multiple selections. Click the sheet tab to select all panes in a sheet. Selected panes have areen tabs.

By default, panes are laid out in Tiled mode. When working in a different mode like Stacked, bring a hidden pane to the front by clicking on its tab.

# Layout

Click the *Layout* icon on the right side of the tab to change the arrangement of the images in a pane. The options are the same as for panes in a sheet.

# 2.5.3 Images

By default, the gel name is displayed in the upper left corner of an image. The color of the name indicates if the image is selected (green) or not (gray). If an image name has a red corner, this means that it is the reference image for the matching. If an image name is a darker green or darker gray than the other images, then it is the current sheet reference. All other images in the sheet are compared to this sheet reference when using the menu commands in *View*: *Sheet*.

To hide the image names, choose *View: Global: Show Gel Names*. The image name is replaced by three dots.

Move the mouse over an image name to view a screentip specifying the match hierarchy. In a Class sheet, the class and path to which it belongs are also given. If the image is the sheet reference, then this fact is included in the screentip.

## Selection

Select an image by clicking on its name. Use the Shift or Ctrl keys to make multiple selections. Click the pane tab to select all images in a pane.

When images are hidden like in Stacked mode, bring an image to the front by clicking on its tab at the bottom of the pane. Quickly sift through images using the Page Up and Page Down keys on your keyboard or click on the navigation triangles in the lower right corner of a pane.

## **Scrollbars**

The scrollbars on the right and bottom edges of each image can be used to change the zoom factor of an image by dragging one of the ends of the scrollbar. You can also move the visible area of the gel horizontally or vertically by clicking in the middle of the scrollbar and dragging left/right or up/down.

Click on the gray square at the intersection of two scrollbars to reset the image to its full image size.

If you want to view a specific area of your gel image, choose *View: Global: Scrollbars: Adjust*. In the Adjust Visible Area window, set the exact horizontal and vertical start and end coordinates for the area to be displayed. You can do this in terms of different units: Pixel coordinate, Percentage and Real coordinate (pl and MW). Please note that pl\_MW annotations must be defined in order to use this function

To hide the scrollbars, choose *View: Global: Scrollbars: Show.* In that case, You can still move and zoom images by using the Move and Zoom tools in the toolbar.

#### 2.5.4 Switch order


You can change the order in which images are displayed by dragging the gel name onto another image. It is then inserted before this image. Similarly, you can re-order panes by dragging their tabs to a new position.

Swap panes or images by choosing *View: Sheet: (Navigation:) Switch* or the Ctrl+F shortcut. This reverses the last re-ordering operation applied to a pane or image.

## Workspace 2.6

The Workspace plays an important role in the software. It allows you to organize your gels into projects, to specify how the gels are to be matched together, and to define your classes (or groups) for statistical analysis. A brief description of the Workspace window is presented here. The utility of the Workspace is discussed in another chapter.

The Workspace window is a dockable window and has two main parts: the Workspace toolbar and the Navigator (Figure 2-4).



**Figure 2-4.** The expanded view of the Workspace window. The Workspace toolbar and Navigator (at the left) are always visible. The file details in the expanded view (at the right) are only displayed when desired.

## 2.6.1 Toolbar

Commands to create new projects, and files to projects, remove, backup and restore projects are found under the **Project** icon in the Workspace toolbar.

Click the **Expanded View** icon to enlarge the Workspace window to include details of all files selected in the Navigator. The files can be sorted in ascending or descending order by clicking a column header. Click the **Expanded View** icon again to hide the file information.

# 2.6.2 Navigator

The Navigator displays all files and folders you can view (Image Pool) and analyze (projects). The look and feel is similar to Windows Explorer. There is a hierarchical structure of folders, subfolders and files that can be expanded or collapsed, dragged and dropped in a new location, copied and pasted, etc.

When launching the software for the first time, there is an empty Image Pool folder. Images are opened through the Image Pool. To further analyze images (detect, match, carry out statistical analysis, etc.), they must be transferred to a project. You can work on many projects in the Workspace, each of which will contain one or more root structures with Match and Classes subfolders.

# **Current sheet**

Any item in bold in the Navigator corresponds to items displayed in the current sheet.

## Contextual menus

A contextual menu containing relevant options appears when you right-click an item in the Navigator.

# **Navigator icons**

Each folder type in the Navigator has a specific icon (Image Pool, Project, Match, Class). In addition, images and match set folders have icons that indicate their status. You can distinguish DIGE gels from non-DIGE gels, know which gels have been detected and matched, and recognize gels used as references for matching.

Icon	Meaning
<b>6</b>	Image Pool folder.
	Project.
Ģ.	Match folder.
Ę.	Classes folder

Table 2-1. The folder icons in the Navigator.

Icon	Meaning
	Undetected image / reference image.
<b>**</b>	Detected image / reference image.
<b>*</b>	Matched image / reference image.

<b>A</b>	DIGE image (in class).
& A	DIGE gel / reference gel.
<b>№</b> ♣	Matched DIGE gel / reference gel.
<b>№</b> №	Match set / reference match set.
<b>№</b> ♣	Matched match set / reference match set.
<b>■</b> 13	Class.

**Table 2-2.** The Navigator icons inform you about the type and status of the image, gel and match set.

# 2.7 Reports

Reports are highly practical for organizing and describing your gel data. They make it much easier to process all of the information.

Reports are not necessarily in table format. Graphical representations of data such as histograms, scatter plots and 3D views are treated as reports as well. All reports are dockable windows.

Report	Description
3D View	A three-dimensional view of selected gel regions or areas around selected spots.
Analyze Gels	Information about each selected spot match such as its Match ID, value for each spot in the match, and chosen statistical measures calculated on all spots in the match.
	Scatter Plots also provide information (slope, offset, correlation coefficient, fitting error) that compare the spot values for two gels.
	The Match Statistics Table displays the number of matches and percentage of matches for each gel at the selected level of the match hierarchy.
	The DIGE Histogram displays the frequency distribution of the DIGE volume ratios.

Report	Description
Analyze Classes	Central tendency, dispersion, and overlapping measures for classes of gels computed for all selected matches.  Differences between the spot values in several classes can also be quantified with Statistical Tests such as ANOVA, Mann-Whitney U test and Kolmogorov-Smirnov test.
Gel Table	Summarized information about the selected gels, such as Gel ID, file name and path, gel calibration data, gel resolution and size, number of detected spots, user defined properties such as sample type, staining or date of the experiment, and much more. Predefined templates (Properties, Files, Descriptions, Calibration) allow you to quickly display a specific subset of fields. You can also create custom templates.
Spot Table	Specific information about selected spots such as Spot ID and coordinates, quantification values, attached labels, etc.
Annotation Table	Information about annotations, including the label content for each category, the annotation coordinates and Spot ID (if the annotation is linked to a spot).

#### 2.7.1 **Dynamic content**

Report content is continuously updated. A report selection is synchronized to reflect the most current data from the corresponding sheet that contains the gel images.

**Note:** Reports are attached to their corresponding sheet. If the sheet is closed, the reports will be closed as well.

## 2.7.2 Content based on enabled spots

By default, all spots are enabled and therefore represented in the reports. But the content of reports can be limited to a subset of spots by deactivating spots that are not of interest.

#### 2.7.3 **Toolbars**

Most reports have a toolbar with the following icons in addition to some reportspecific icons and functionalities that are described in later chapters of the User Manual:



# Suspend synchronization

By default, the selection in the report is continuously synchronized with the corresponding sheet that contains the images and with other reports. Click the **Suspend Synchronization** icon to stop the synchronization and render the selection in the report independent of the selection on the gel, and vice versa.



Enter a name and select the desired format for the file to be saved. Tables can be saved in tab-delimited text format (.txt), as a Microsoft Excel Workbook (.xls), or in XML format (.xml). Graphics can be saved in PNG, TIFF or BMP formats.



Normal print options are available when printing graphical reports. When printing tabular reports, the table is first displayed in your default Web browser. This is because the XSL stylesheet located in the Template\Reports folder of the ImageMaster installation directory is used to transform the XML report into an attractive table. You can then use the print option in your browser to get a printout.

# Copy to clipboard

Export your data directly into another application. First select the desired lines in a table, or the desired graphics in a window, using the Shift or Ctrl keys. Then copy the selection to the clipboard by choosing the **Copy to Clipboard** icon. Paste directly into the preferred software.



# Previous in selection

Click the **Previous in Selection** icon to skip to the first selected item encountered when scrolling towards the top of your table.



# Next in selection

Click the Next in Selection icon to skip to the first selected item encountered when scrolling towards the bottom of your table.



# **Settings**

Click this icon to change the display settings for your report. In most cases, this means setting the visibility (hidden or shown) of the columns. In some cases (Histograms and Match Statistics Table), you can define other display options.



# Y Select by value

This feature allows you to select items in the report based on a numerical search criterion. Click the Select by Value icon, then select the measure (that is, column) you want to use for refinement, and finally set the lower and/or upper limits of your search interval.

## 2.7.4 **Customize reports**

# Sorting data in columns

Data in tabular reports can be sorted by the column content. If you click once on the header of a specific column, a triangle is displayed indicating that the column's numerical or textual data is sorted in ascending order. When you click once more on the header, the triangle inverts indicating that the data is sorted

in the opposite order. Note that ascending order means that numbers are sorted from 0 to 9 and text is sorted from A to Z.

# Column visibility

Load a predefined template or save your own template indicating what columns should be hidden or shown using the Settings icon in a report toolbar. This is particularly useful when you only want the essential information to appear for clarity or for printing purposes.

# To apply an existing report template:


- Click the **Settings** icon in a report toolbar.
- In the template window (Figure 2-5), select a template name with the **Load** icon. By default, only the Gel Table has more than one predefined template (Properties, Files, Decriptions, Calibration).
- Click OK.

# To create and save your own report template:

- Click the **Settings** icon in a report toolbar.
- In the template window, select (box is checked) the attributes that you want to show in the report and deselect (box is unchecked) the attributes to hide columns.
- To save the template for later use, click the **Save** icon. Choose a name from the list or select < New...: to create one. Click OK.
- Click **OK** in the template window.

# To delete an existing report template:

- Click the **Settings** icon in a report toolbar.
- 2 In the template window, click the *Remove* icon and choose the template to be deleted. Click **OK**.



3 Click **Yes** when asked for confirmation.

Figure 2-5. Report template window.

## Column order

You can reorganize table columns directly in a report window. To do so, drag the column header to its new position. You will see red arrows at the insertion point.

## Column size

You can enlarge or reduce the column size. Drag the right edge of the column header until the column is the width you want. To resize columns so that their whole content is displayed, double-click on their separator. The column to the left of the cursor is resized.

## 2.7.5 Edit table cells

In the Spot table, cells containing annotation labels are editable. Double-click in a cell to start typing or editing your label. When finished, a single click in any cell quits the editing mode.

Please note that annotation categories of the data type Set are displayed as check boxes. A checked box means that the item belongs to the Set. An empty box means that it does not belong to the Set.

The same is true for spot sets represented in any of the tables: a checked box means that the spot belongs to the set. An empty box means that the spot does not belong to the set.

# To check several boxes simultaneously:

- 1 Use the Shift or Ctrl keys to select the rows in which you want to check or uncheck the box.
- 2 Click in the box of one of the selected rows. All rows are now either checked or unchecked.

#### Dockable windows 2.8

The Workspace, reports, and Adjust Contrast window are dockable windows. Dockable windows make it easier for you to work with numerous windows at the same time. These windows can be docked (that is, fixed in place) against the left, right, top or bottom edges of the ImageMaster window and always lie on top when visible. A visible window can be in one of three modes:

## Pinned

The pinned mode 4 enables dockable windows to be locked into position around the edges of the ImageMaster window. Once a window is pinned, you can move it to a different location by dragging the title bar (Figure 2-6). Guides indicate where the window may be docked. By moving the cursor over the quides, a shaded blue box appears showing where the window will reside if the left mouse button is released. If you move the dockable window to a nonpredefined location, it becomes a floating window. Moving a docked window may affect the location and size of other docked windows.

## **Un-pinned**

A visible window in un-pinned mode 🔁 automatically collapses when not in use to become a tab at the edge of the ImageMaster window (in Figure 1-7, the Workspace is in this mode). When you click on a docked tab, the window slides back into view and is ready for use. You can also click on Minimize \_ to minimize a window in un-pinned mode.

## Floating


A dockable window in floating mode will always appear on top. It can be dragged to any position within the software or even outside the ImageMaster window. You can switch in and out of floating mode by double-clicking on the title bar of a dockable window.

# Tabbed groups

Dockable windows can be organized into tabbed groups. This feature extends your ability to maximize the use of limited screen space by combining multiple dockable windows into one window. In order to form a tabbed group, drag the title bar of a dockable window into the center of another. You will see the nested tabs at the bottom of the docked window. In order to separate a tabbed group,

# 2.9 Options

drag a tab away from the docked window or double-click on the tab. Please note that tabular reports and graphical reports cannot be grouped together.



**Figure 2-6.** The Class Analysis Histograms window previously docked at the right edge will now be docked to the right of the Class Analysis Table window, in the lower left corner of the ImageMaster window.

# 2.9 Options

You can set various parameters that influence your work in ImageMaster. These settings are accessed by choosing *Tools: Options* in the menu. More detailed information about the options are provided in the related chapters of the User manual. However, an overview of the settings, per tab, is given below.

# Display

• Indicate the default spot colors (for enabled, disabled, selected and overlapped spots) and the color for the match vectors.

- Choose the units to express coordinates in the Status Bar and Measure tool.
- Indicate if spot selections should be surrounded by boxes, for easier localization

# **Annotation**

Define the annotation categories that should always be available in the software, and set their attributes and display properties.

# Gel descriptions

Define the gel descriptions that should always be available in the software.

# Quantification

- Specify the spot quantification value to be used for non-DIGE (Value) and DIGE (DIGE Value) experiments. This spot quantification value is used in Gel and Class Analysis Reports.
- Compute the spot areas in mm<sup>2</sup>, based on gel resolution (default), or pixels.

- 2 Graphical user interface2.9 Options

# 3 Image Pool

#### 3.1 Introduction

This chapter is about opening and processing images in the Image Pool. You must be familiar with the concepts described in the last chapter (the Workspace, sheets, panes and images) in order to get the most out of it.

### 3.2 Start with good images

Your gels must first be converted into an image file by an appropriate imaging device. During digitization, a gel is resolved into a two-dimensional matrix of squares, or pixels. Each pixel in the generated image file is characterized by its X and Y coordinates, and its signal intensity, or gray value.

To make any analysis meaningful, it is important to start with good quality image files. The following paragraphs give some helpful tips on what resolution, depth and image formats should be used to obtain the best possible results.

#### 3.2.1 Resolution

The scanning resolution of a gel is critical as it influences the amount of visible detail in the image. A low resolution corresponds to a large pixel size or a small number of pixels (or dpi, dots per inch). When the image resolution is too low, individual spots cannot be distinguished. On the other hand, when the scan resolution is too high, the image file becomes very large. This slows down the gel analysis significantly. A resolution between 150 and 300 dpi is generally sufficient for gel analysis.

#### 3.2.2 Image depth

The range of potential gray levels in an image varies according to the image depth. In the case of an 8-bit image, one pixel has  $256 \, (2^8)$  possible gray values (0 to 255). Images scanned with a higher image depth contain more information. A 16-bit image ( $2^{16}$  = 65536 gray levels) will reveal more subtleties. We strongly recommend an image depth of at least 12 bits for gel analysis. 16 bits is preferred.

Please make sure to use gray scale images for your analysis, and not color images. The extra color information (one intensity value for each color channel) is of no value.

#### 3.2.3 Calibration

Some imaging devices can measure more gray levels (for example, 100000 or more) than can be stored in the availabe image formats. These instruments encode the gray values using a nonlinear calibration curve (generally encoding low intensity values with higher accuracy than the high intensity values), to conserve as much information as possible in the saved files. ImageMaster uses the calibration curve contained in these image files to recalculate the measured gray values and use them for display and quantification.

Other devices convert raw pixel values into real-world units (generally optical density or OD). This means that the range of gray values is the same no matter what the original image depth. There are image capture devices, such as the GE Healthcare ImageScanner™ in conjunction with the LabScan™ software, that even allow you to perform this type of calibration. When this is the case, ImageMaster takes into account the conversion tables or calibration formulas stored in the files. You can also calibrate such an image capture device within ImageMaster using calibration step wedges or calibration strips.

#### 3.2.4 Image editing

General purpose graphics software such as Adobe Photoshop ignore or even remove calibration information. Therefore, you should not use them to flip, rotate, crop or invert your images. Instead, use the software that came with your scanner or the dedicated tools in ImageMaster.

#### 3.2.5 File format

The supported input formats are GEL (Molecular Dynamics), MEL (ImageMaster), TIFF (Tag Image File Format), IMG (Fuji), GSC and 1SC (Bio-Rad). Please note that the TIFF format does not include calibration information.

#### 3.2.6 DIGE file naming convention

To facilitate the import of DIGE images, it is recommended that the file names for the group of two or three images contain a common string and their respective dye names (Cy2, Cy3, Cy5).

## 3.3 Preview images

#### 3.3.1 Open images

You can open gel images that are in any of the above-mentioned input formats (.mel, .gel, .tif, .img, .gsc, .1sc).

#### To open gel images:

1 Do one of the following:

- Choose File: Open.
- In the Workspace window, click the **Project** icon. In the Add Files window, search for Files of Type "All Image Files (\*.mel; \*.tif; \*.gel; \*.img)".
- In the Navigator of the Workspace window, right-click on the *Image* **Pool** folder and select **Add**
- In the Open Files or Add Files window, browse the directory where the image file is located, select its name and click **Open**. Use the Shift or Ctrl keys to make multiple selections. Note that you can select files from multiple DIGE gels.
- For the different types of images, do the following:
  - For non-DIGE images: Specify the staining. If the staining cannot be found in the list, you can type a new one. Click **OK**.
  - For DIGE images: Select the images that are part of the same DIGE gel (Figure 3-1). By default, the software proposes image combinations based on the file names. Click **Add** to confirm one DIGE gel at a time, or **Add All** if the proposed combinations are all correct. The suggested names for the created DIGE gels must be edited and confirmed individually.



Figure 3-1. Create DIGE Gel window. The first list contains the Cy2 images to be opened, the second the Cy3 images and the third the Cy5 images. The software proposes the combination of images based on the file names.

**Note:** The software checks the image resolution before opening a file. If the resolution is too low, you get a warning message. If it is higher than 300 dpi, you are able to scale the image in order to reduce its size.

#### 3.3.2 Image pool sheet

The gel images appear in the Image Pool folder in the Workspace and are automatically opened in the Image Pool sheet (Figure 3-2). Non-DIGE images all appear in a single pane with the name Files. A separate pane is displayed for each DIGE gel, containing the images belonging to that gel.

#### Remove

Gels remain in the Image Pool until they are added to a Project. When you are finished working with an image, right-click on it and select *Remove*. To remove all images in the Image Pool, right-click on the Image Pool folder and select *Remove All*. Please note that this action never affects the original image files (which always remain unchanged in their original location).

#### Hide / display

If you do not want to display all gels in the Image Pool sheet, you can right-click on certain items in the Image Pool folder and select *Hide*. To show a hidden image, right-click on it and select *Display*.

To open the Image Pool sheet at any time, right-click on the Image Pool folder in the Workspace and select *Display*.

#### **Properties**

To rename a gel or image, or to see the file path, right-click on a gel and select *Properties*.




Figure 3-2. Image Pool folder in the Workspace and Image Pool sheet.

#### 3.4 **Process images**

The Image toolbar provides the basic tools to process images: Rotate, Flip, Crop and Invert Gray Levels, which can also be found in the *Edit*: *Gels* menu. These tools are only available when the Image Pool is the current sheet. They are also optional and should only be applied when needed.

Processing your images in ImageMaster does not affect your original image files. Copies are saved in a temporary folder until the images are added to a project and saved in the corresponding project folder.


Please note that these operations are simultaneously carried out on the two or three images in a DIGE gel.

#### A SN Rotate 3.4.1

Selected images can be rotated by 90° CCW (counterclockwise) or 90° CW (clockwise). Free rotation (the third icon) can also be applied, although it should be avoided since it modifies the original data.

#### To freely rotate an image:

- Select the gel image to be rotated in the Image Pool sheet.
- Click the **Free Rotate** icon.
- A red grid appears on the image. The bold horizontal grid line plays the role of landmark to help you visualize the rotation. It becomes the new horizontal in your rotated image.
- Click anywhere in the image and rotate the grid while holding the left mouse button. Release the button when the bold line is parallel with what should be the new horizontal in your image (Figure 3-3). You can also manually enter a rotation angle in the Rotation Tool window.



**Figure 3-3.** Rotation Tool window. The grid is rotated until its bold line is parallel with what should be the new horizontal reference. When the mouse button is released, the gel image is rotated.

#### 

When gel images are scanned in the wrong direction, you can Flip Horizontally or Flip Vertically to produce their correct mirror image.

# 3.4.3 **4** Crop

You can crop your gel images with the *Crop* tool. This creates new gels that only contain the selected area and removes the outer area.

When you crop a gel, you get the choice to create a new image (a number is appended to the existing name) or to overwrite the image in the Image Pool.

**Note:** Gels in the Image Pool are copies of the original image file.

# Crop area

The crop area is a region that can have an anchor attached to it. You can position the anchor on an easily recognizable protein spot. As the region moves with the anchor, and vice versa, you can easily crop a similar part of each gel by correctly positioning the crop areas (of the same size) in the gels (Figure 3-4).

#### To define a crop area:

- Click the **Region** tool and place the cursor at the top left position of the area you want to crop. Hold down the left mouse button and move the cursor to the bottom right position (a dashed box is displayed). Release the mouse button at the end point.
- Move a crop area by clicking inside the box and dragging it. Change the size of the area by dragging a corner or edge. If the box is reduced to its minimum size, its appearance changes to a blue circle. To remove an area, double-click on the image.
- If you want to attach an anchor to the crop area, hold the Alt key while clicking on an easily recognizable protein spot. A dark blue circle will be centered on the spot. Note that this anchor may be located inside or outside the crop area.
- To change the position of the anchor, hold the Alt key and click another spot. To remove the anchor, hold the Alt key and click on the anchor.
- Propagate the crop area to the other images in your sheet by holding the Shift key while clicking in the crop area.
- Adjust the position of the crop area in each gel by moving it so that the anchor is centered on the appropriate spot.




Figure 3-4. Identical crop areas in three gel images.

Click the Crop icon to crop the image according to your area.



Crop areas can be exported and imported to ensure that the final size of all your cropped gels is identical, even between work sessions.

#### To export a crop area:

- 1 With the **Region** tool, define a crop area as described above.
- 2 Select the gel in which you defined the crop area.
- 3 Choose *Edit*: *Gels*: *Crop Area*: *Export* to save the crop area to a file with the extension .cpt (Crop Tool).
- 4 Browse to locate a folder and enter a file name. Click **Save**.

#### To import a crop area:

- 1 Select the gels to be cropped to a previously defined crop area.
- 2 Choose Edit: Gels: Crop Area: Import and select the previously saved file. Click Open.
- 3 The crop area appears on the selected images.
- 4 By clicking inside the area with the **Region** tool and dragging it, you can move the crop area to superimpose the anchor on a characteristic spot.

# 3.4.4 Invert gray levels

You can invert the gray levels of selected gels. This means that if your image shows white spots on a black background, the inversion displays black spots on a white background (the required mode for analysis in ImageMaster).

**Note:** When your images open in the software with white spots on a black background this may indicate that you used incorrect scanner settings or that your files are not properly imported. Please verify your image acquisition parameters. If you think the software does not correctly support your files (that is, saved in one of the recommended input formats), please contact our technical support service.

# 3.5 Calibrate images

#### 3.5.1 Display calibration information about images

It is recommended to use calibrated images for your gel analysis. To verify if your images were calibrated, and possibly view the calibration information:

Choose Reports: Gel Table, click the Settings icon and select the
Calibration template from the Load drop-down list. The Calibration Unit,
Formula, Name, Creator and Date are displayed in the report. If nothing
appears in these columns, the gels were not calibrated. Create a
calibration according to the following sections.

If your image was calibrated with LabScan 5.0 or 6.0, you can select the image and choose Edit: Gels: Show Calibration Plot to view the calibration curve. See below for details about the calibration curve.

If you are digitizing your images using a flatbed document scanner and do not have calibrated images yet, you can prepare and apply a calibration at this point, before adding your images to a project. It is generally recommended to perform this type of intensity calibration on the image capture device once every month.

#### 3.5.2 Prepare for a calibration of the image capture device

#### Scan a calibration step tablet

To calibrate the image capture device, you need to scan a calibration step tablet or calibration strip along with your gels. These step tablets have known intensity values (expressed in optical density, OD, or diffuse density, DD) published by the manufacturer of the step tablet.

Please note that for the purpose of 2D gel analysis, it is only useful to calibrate the image capture device when working in transparent mode. No calibration needs to be done when you do reflective scanning. When both a transparent and a reflective calibration strip are provided, be sure to use the appropriate calibration step tablet.

#### Calibration tablet file

The OD values for the step tablet have to be specified in a Calibration Tablet File, together with other information such as the height and width of the tablet, and the number of steps.

An example of such a Calibration Tablet File (Kodak2.tab) can be found in the Template\Tablet folder of the ImageMaster installation directory. This Calibration Tablet File is made for use with the Kodak Step Tablet no. 2. If you do not use this specific step tablet, you can copy the file and edit the data to make your own Calibration Tablet File. You can edit the file with tools such as Windows Notepad.


As the intensity values supplied with your step tablet are generally expressed in diffuse density (DD), you have to convert them to OD values. For this purpose, the manufacturer of the step tablet should provide the appropriate relationship. For the Kodak Step Tablets no. 2 and 3, for example, this is OD = 1.4 DD.

#### 3.5.3 Create calibration of the capture device

You can create a calibration once you have scanned the step tablet and have a correct Calibration Tablet File.

#### To create a calibration:

- 1 Open the step tablet image file. If necessary, rotate the image so that the light steps are displayed at the top.
- 2 Choose **Tools**: **Calibration Tablet**: **Create**.
- 3 In the Load Step Tablet Definition window, browse to the folder where you saved the Calibration Tablet File (.tab) specifically tailored to your step tablet (see above), select the file and click *Open*.
- 4 A red calibration step overlay appears on the image and the Create Calibration window is displayed (Figure 3-5).
- 5 Adjust the position of the steps by dragging the overlay while holding the left mouse button. The size of a step on the red overlay (the distance between two short horizontal lines) should correspond exactly to the size of a step on the image. If this is not the case, you must adjust the height of the tablet in the Calibration Tablet File.



**Figure 3-5.** Image of the step tablet with the red calibration step overlay. The Create Calibration window shows the calibration curve and the OD values for the different steps (on the left). When selecting a step in the list, the corresponding step becomes

automatically highlighted in green on the step tablet overlay and in the calibration graph.

- At the left of the Create Calibration window, you see the theoretical optical density (OD) values of the different steps in the tablet (the values you entered in the Calibration Tablet File). Some steps are automatically deselected (grayed out) because of their unreliable values. You can deselect additional ones if you estimate that they should be excluded from the calibration process. At the right of the Create Calibration window you see the calibration curve between the logarithmic transmittance values on the X-axis, and the OD intensities on the Y-axis. Note that the measured intensity values for each step are calculated as median intensities over all the pixels in the small rectangle area for each step (on the step overlay). The horizontal dispersion intervals in blue (or gray for deselected spots) represent the intensity ranges when 10% of the less intense and 10% of the most intense pixel values are removed. The calibration formula and error are given below the graph.
- You can display some reports to judge the quality of your scanner calibration (see below for more details).
- Once you are satisfied with the calibration, close the Create Calibration window. The software asks if you want to apply this new calibration. If you answer Yes, ImageMaster applies the calibration to the image.

The icons in the toolbar of the Create Calibration window and their use are:



**Open** another step tablet definition.




Save the calibration (with the extension .cal).



Print the calibration graph.



**Copy** the calibration graph to the clipboard.



- The *Fitting Table* displays the calibration formula.
- The *Calibration Table* displays for each step: the step number, the measured average gray level, the theoretical intensity value, and the fitting error (difference between the curve and the point).

### 3.5.4 Apply a calibration to an image

Once you have created a calibration using a step tablet, and saved this calibration in the step tablet image file or a calibration file (.cal), you can apply the calibration to newly-scanned image files.

#### To apply a calibration to a gel:

- 1 Display and select the gels the calibration should be applied to.
- 2 Choose Edit: Gels: Apply Calibration.
- 3 Select the source of the calibration information. This can be an open gel (or step tablet image) that was already calibrated or you can select a file (.cal or .mel) from the hard disk.

All pixel and spot values subsequently calculated and displayed in the software correspond to the calibrated values.

#### 3.5.5 Remove a calibration

You can remove a calibration from a gel. Note that this can also be done for gels that were already calibrated when you imported them into the software.

#### To remove a calibration from a gel:

- 1 Display and select the gels from which you want to remove the calibration.
- 2 Choose **Edit**: **Gels**: **Reset Calibration**.
- 3 Confirm your choice.

#### 3.5.6 Control a calibration

The Control calibration mode allows you to verify if you are using the correct calibration. It requires a different, specially calibrated step tablet (for example, Kodak Step Tablet no. 3), which you compare to your previous calibration results. In other words, you have a calibration step tablet for everyday use and a specially calibrated control step tablet to verify your calibration periodically.

#### To control a calibration:

- Scan your control step tablet (for example, Kodak Step Tablet no. 3) and open the image in ImageMaster. If necessary, rotate the image such that the light steps are displayed at the top.
- 2 Choose Tools: Calibration Tablet: Control.
- 3 You are asked to load the calibration to be controlled. This calibration could have previously been saved using the *Save* icon in the Create Calibration window (.cal) or can simply come from a calibrated image file (.mel) such as the calibrated step tablet image obtained in the section above.

- 5 The red calibration step overlay appears on the image of the control step tablet and the Control Calibration window is displayed.
- 6 Select the **Spot** or **Annotation** tool in the ImageMaster toolbar and adjust the position of the steps.
- You should now verify that the calibration curve is passing through the data points correctly and with minimum dispersion intervals. If this is not the case, try to find out why your current calibration does not seem to work properly.

### 3.6 Describe images

### 3.6.1 Gel descriptions

A gel description is information you can enter about the gel image, which can be used for later reference by yourself or any colleagues. This information can include sample type, gel running protocol, date of the experiment, operator name, pH range, SDS gel percentage. All of this data is entered as Gel Descriptions.

When opening images in the software, you are asked to specify the staining. Staining is a standard Gel Description. This information is simply informative, that is, it is not used by the software except when you specify the Cy<sup>TM</sup> dye for DIGE images that do not have this information in their file name.

Each image can only have one description of a given category. For example, a gel description category Treatment could contain the description "Drug 1" for

certain images and "Drug 2" for other images (Figure 3-6). You can define any number of categories for an image or a given set of images.




Figure 3-6. Gel descriptions in a Gel Table (using the Descriptions template).

#### 3.6.2 Display gel descriptions

Gel descriptions can be displayed in a Gel Table. Select the predefined Descriptions template in this table to display only the columns for the gel description categories.

#### To display gel descriptions:

- 1 Choose **Reports**: **Gel Table**.
- 2 Click the **Settings** icon.
- 3 In the Gel Table window, select the predefined Descriptions template from the *Load* icon drop-down list. The boxes for the FileName and all gel description categories (by default, only Staining and Comment) are checked.
- 4 Click OK.

### 3.6.3 Add or edit gel descriptions

You can edit a gel description for one or more selected gels.

#### To add or edit gel descriptions:

- 1 Select the images for which you want to add or edit descriptions.
- 2 Choose Edit: Gels: Edit Description.
- 3 In the Add Gel Description window, select an existing gel description category or create one. Click **OK**.
- 4 Enter the gel description that applies to the selected images. Click **OK**.
- 5 The new gel descriptions are displayed in the Gel Table.

Alternatively, you can define gel descriptions by clicking on the *Add Description* icon in the Gel Table toolbar. It's behaviour is identical to that of the *Edit: Gels: Edit Description* menu.

You can delete all gel descriptions of a certain category for the selected gels by choosing *Edit*: *Gels*: *Delete Description*.

#### 3.6.4 Permanent gel description categories

When gel description categories are created as described above, they are only used for the gels that were selected during the creation of the categories. To make gel description categories permanent in the software (always availabe from the category list), you must define them in the Options.

#### To create permanent gel description categories:.

- 1 Choose **Tools: Options**.
- 2 Click the **Add in the Gel Descriptions** tab.
- 3 Enter a category name in the *Add Category* box and click *OK*.
- 4 The category is displayed in the permanent category list.

To remove a permanent category, select it from the list and click **Delete**.

- 3 Image Pool3.6 Describe images

# 4 Projects

#### 4.1 Introduction

A project includes all gels, spots, matches, annotations, spot sets, and other information produced and analyzed during the course of a specific gel study. You can create or add many projects in the Workspace.

A project in the Workspace can include one or more match hierarchies (for example, AB in Figure 4-1), each of which contains a Match folder and a Classes folder. The Match folder describes how gels or populations of gels, called match sets, should be matched together. The Classes folder is where the biological question is stated, through the definintion of classes of gels to be compared.





Figure 4-1. Project structure.

The next sections further explain the principles of match hierarchies and classes, and how they can be used.

#### 4.1.1 Match hierarchies

All images in an experiment are not equally easy to compare, even when the gels are run in a highly controlled way. Typically, gels belonging to the same biological group are easier to match (that is, corresponding spots are easier to find) than images from different biological populations.

It is therefore recommended to use hierarchical match structures to create more efficient match designs. Figure 4-2 and Figure 4-3 show an example of a match hierarchy, or root match set, AB and submatch sets A, B, AT1, AT2, BT1 and BT2. Gels or match sets with a red marker are used as the reference in the matching and always appear first in the list.



**Figure 4-2.** The Match folder. In this example experiment, a set of samples from bacteria were cultivated with either substrate A or with substrate B. Under both growing conditions, two treatments were tested. Therefore, 4 different populations exist. A gel was run for each of the 3 samples in a population, giving a total of 12 gel images.




Figure 4-3. A tree view of the example match hierarchy described above.

An important advantage of this type of match hierarchy is that it minimizes the number of difficult match combinations. Instead of having to match 9 images (those belonging to AT2, BT1 and BT2) to a global reference (for example, A\_T1\_Gel1) that is not belonging to the same population, only 3 slightly more complex match combinations must be performed (AT1 versus AT2, BT1 versus BT2, and A versus B). This significantly reduces time spent on match editing.

There is another important reason why it is useful to adopt hierarchical population matching instead of matching all images against a unique arbitrary reference image. Only spots matched with a spot in another gel are included in Gel and Class Analysis Tables (all spots are of course presented in the Spot Table). The likelihood of a spot being matched is much higher when matching with a gel from the same biological population. Spots that are represented in a single population (submatch set) are therefore included in the analysis, even if they are not in the global match reference. This considerably reduces the number of spots missed in the analysis.

**Note:** DIGE gels are inherent match sets. A DIGE gel (the entity with its composing images) is treated as any other non-DIGE gel when setting up a match hierarchy.

Once created in the Workspace, you can display a complete match hierarchy (for example, AB in Figure 4-2) in a sheet and carry out spot detection on all the included gels. After defining one or two landmarks, the entire experiment is matched in a matter of seconds, and matches are automatically propagated at each level of the match hierarchy.

#### 4.1.2 Classes

In the Classes folder, you state your biological questions. This means that you define a class for each set of gels that you want to compare with other such entities. Your goal, by comparing classes and therefore the gels within those

classes, is to find the protein expression variations between different biological states.

The classes created in the Classes folder of a given match hierarchy can contain any of the images in that hierarchy. One image can be in different classes. In Figure 4-4, for example, the same 12 images can be compared as part of the 4 classes AT1, AT2, BT1 and BT2, or as part of the two classes T1 and T2.

You can define classes at any time, even in the very beginning of your gel analysis study when no spots are detected. To carry out statistical analysis, however, your gels must be detected and matched.




Figure 4-4. The Classes folder.

## 4.2 Create a project

#### 4.2.1 The first time the software is launched

As long as no projects have been created or added to the Workspace (for example, the first time you open the software), you are prompted to create a

project (see Figure 4-5). Enter a project name, browse to where the project folder should be saved, and possibly add a comment.




Figure 4-5. New Project window.

### 4.2.2 At any time

You can create a new project at any time by selecting **New** from the **Project** icon drop-down list in the Workspace toolbar (Figure 4-6). Again, enter the **Project Name**, **Location** and **Comment** (Figure 4-5).





Figure 4-6. Project icon drop down menu in the workspace.

### 4.2.3 Add files to project

Use the *Add Files to Project* icon in the Image toolbar to add gels from the Image Pool to a project (Figure 4-7). Add an existing project or create a new in

#### 4.3 Create match hierarchy

the Add Files to Project window. An existing or new match set name must be entered as well.



**Figure 4-7.** Add Files to Project window. A project can be created by choosing <New...: from the Project list.


### 4.3 Create match hierarchy

#### 4.3.1 Create a match set

The easiest way to create a match hierarchy is by adding gels from the Image Pool to a project. The idea is to select only the gels that should be added to a particular match set (for example, gels from the same population or same experimental batch) and then to create this match set. Different options exist:

- Select gels in the Image Pool folder and drag them onto a project name. Enter a name for the new match set.
- Select gels in the Image Pool sheet and drag them onto a project name in the Workspace. Enter a name for the new match set.
- Select gels in the Image Pool sheet and click Add Files to Project in the Image toolbar. Set or create the destination project, and then click <New...: in the MatchSets field (Figure 4-8). Enter a name for the new match set.

Alternatively, create an empty match set by right-clicking on the project name and selecting *Create MatchSet* in the contextual menu. Drag images from the Image Pool (folder or sheet) into the new match set. Gels selected in the Image Pool sheet can also be added to an existing match set by clicking the *Add Files to Project* icon.




**Figure 4-8.** Add Files to Project. The three images of population A\_T1 are selected in the Image Pool sheet. After having clicked the Add Files to Project icon, the destination project and new name for the MatchSet can be set.

#### 4.3.2 Merge a match set

Once match sets have been created, they can be merged into higher level match structures for further matching. Select the match sets, right-click on one

of them, choose *Merge MatchSet* from the contextual window, and give a name for the new match set (Figure 4-9). A hierarchy (Figure 4-2) can be created.



**Figure 4-9.** Match sets BT1 and BT2 are selected for merging. Match sets AT1 and AT2 were already merged this way into a match set A.

#### 4.3.3 Set reference

The reference for each match set must be carefully chosen. This is because automatic matching compares the spots in the reference to those in the other images. If a spot is absent from the reference, it cannot be matched automatically (although it can be matched manually with spots in other gels). The rule of thumb is to choose the gel or match set with the most and the best quality spots as the reference.


Within each match set, the gel or match set that has a red marker and appears first in the list is used as the reference in the matching process. To change the match reference, drag the desired gel or match set onto the name of its parent match set so that it moves into the first position.

**Note:** You can change the match reference as long as the images in your match set have not been matched. Once they are, the reference image can no


longer be changed.

### 4.3.4 Use existing match sets

You may wish to carry out several analyses, using different matching schemes. You can copy match sets to use them in another configuration. In Figure 4-10, for example, the match sets A\_T1, A\_T2, B\_T1 and B\_T2, were copied to be used in populations per treatment (T1, T2) rather than growing substrate (A, B). Because existing matches are conserved when copying match sets, this can save a lot of work.



**Figure 4-10.** Match sets AT1, AT2, BT1 and BT2 were copied to be used in the match hierarchy T1T2.



**Figure 4-11.** A tree view of the match hierarchy in Figure 4-10. The match sets A\_T1, A\_T2, B\_T1, and B\_T2 have been reused in a different configuration compared to the example in Figure 4-3.

You can simply drag one or more match sets to a new destination, or copy (right-click on the match set) and paste them in the new destination (right-click on the destination project or match set). The destination can be a match set or open project. Copies of the match sets, including the existing matches, are created. A number is appended to the original name. You can then rename the copied match sets (see below).

To copy a match set in the same project, hold down the Ctrl key while dragging the match set onto the project name.

### 4.3.5 Export / import a match set

The easiest way to provide others with access to the data is to export a match hierarchy (the match set that is the parent of the project, containing the folders Match and Classes). Simply right-click on the parent match set and choose *Export MatchSet*. The entire match set (including images, matches, spots, annotations, spot sets) is compressed into a single .exp file.

A .exp file can be imported into a project by right-clicking on the project name and choosing *Import MatchSet*.

If the project data are saved in a folder on a shared network, colleagues having access to this folder can open and work with the project, and therefore the match sets.

#### 4.4 Create classes

#### 4.4.1 Create a class

In the Match folder of a hierarchy, select all images that should be added to a single class (group of biologically-related images). You can select match sets as well. However, only one type of item (images or match sets) should be selected at a time.

- Option 1: Drag the selection onto the Classes folder of the same hierarchy.
   Enter a name for the new class and click OK. You can also use the Copy and Paste options in the contextual menus.
- Option 2: Right-click on one of the selected items and choose Add In Class. Enter a name for the new class and click OK.

Alternatively, create an empty class by right-clicking on the Classes folder and selecting *Create Class* in the contextual menu. Drag images or match sets into the new class.

# 4.5 Handle project items

Several operations are available for most items in a project. These are often found in a contextual menu by right-clicking an item.

### 4.5.1 Display

To display one or more images, match sets or classes in a new sheet, right-click on a selected item and choose Display. Please note that only complete match sets or classes are displayed. This means that if you select one image in a match set, all images in the match set are shown.

In the resulting sheet, only the highest level and lowest level items (match sets or classes) are specifically visualized, in the sheet and panes, respectively. To select an intermediary level, click its icon in the pane tab of the match reference for that level.

Once the images are displayed in a sheet, you can start working with them. You can change the layout settings to focus on certain images and hide others.

#### 4.5.2 Remove

To permanently remove an item from the project, select **Remove** in its contextual menu. This will delete the item from your hard disk.

#### 4.5.3 Properties

When choosing **Properties** for an item, some of its attributes are shown such as its name, creator, file path, etc. Here is where you can change the name of an item or enter a comment to describe the item.

#### 4.5.4 Move

You can rearrange your images, match sets and classes to change their position in the list or move them into another match set or class. Drag your image, match set or class to the desired position. It is generally inserted after the item you drop it on. Whenever there is a possibility to insert it inside or after an entity, you are asked to specify.

### 4.6 Save projects

#### 4.6.1 Project folder

When you create a project, you are asked to specify a name and a location on your hard disk. All the data related to the project (images, spots, matches, annotations, spot sets) are saved in this location, in a folder with the name of the project. All users that have access to this folder are able to open, view and work with the project.

The project folder contains the following files and folders:

- **Projectname.prj**: This project file is the link between all the other data files in the project folder. If you want to add an existing project to your workspace, you must search for, and open, its .prj file.
- Raw Images: This folder contains the raw image files in .mel format.
- MatchSet: This folder contains the gel (.gda) and match (.mda) data in subfolders. The .gda files are used to store the spot information. The .mda files contain the match information.

#### 4.6.2 Save

The software automatically saves your work when you close a sheet or exit the software. You can also save your work by choosing *File: Save*. This saves all your data in the corresponding project folder on the hard disk.

#### 4.6.3 Share project

If the project data are saved in a folder on a shared network, colleagues having access to this folder can open and work with the project. The project must be added to the colleague's Workspace using the procedure described below.

### 4.7 Manage projects

#### 4.7.1 Add / remove projects

#### Remove

You can remove a project from your workspace by selecting *Remove* in the *Project* icon drop-down list in the Workspace toolbar. Select the project(s) to be removed and confirm. This does not delete the project folder from your hard disk.

To permanently remove a project folder, you must delete it from your hard disk. Make sure that this folder is not accessed by colleagues in your network before you delete it.

#### Add

A previously removed project can be re-inserted to your workspace at any time. Another user can also add it to his/her workspace.

#### To insert an existing project:

- 1 Select **Add** in the **Project** icon drop-down list in the Workspace toolbar.
- In the Add Files window, browse the directory where the project file (.prj) is located, select its name and click Open.

### 4.7.2 Add files to project

**Add** in the **Project** icon drop-down list in the Workspace toolbar also enables you to add project elements instead of entire projects. You can open the following file types:

- Project (\*.prj)
- Project Backup (\*.bkp)
- Export File (\*.exp)
- All Image Files (\*.mel; \*.tif; \*.gel; \*.img)
- MatchSet Data (\*.mda)
- Gel Data (\*.gda)

#### 4.7.3 Backup / restore project

It is good practice to do regular backups so that you can recover your work at any time. With the Backup function in the Workspace, one or more project(s) can be archived by writing all project-related data (images, spots, matches, annotations, spot sets, enabled spots) into a single compressed file with the extension .bkp. This file can then be restored when needed.

#### To backup a project:

- 1 In the Workspace, select **Backup** in the **Project** icon drop-down list.
- 2 Select one or more projects to backup, and confirm when prompted.
- 3 In the Backup Project window, browse to the directory you want, enter a file name and click **Save**.
- 4 The backup file (with the extension .bkp) is archived.

#### To restore a project:

- 1 In the Workspace, select **Restore** in the **Project** icon drop-down list.
- 2 In the Restore Project window, browse the directory where the backup file is located, select its name and click *Open*.
- 3 Select the project(s) to restore from the list. Use the Ctrl or Shift keys to make multiple selections. Click *Restore*.
- If the project already exists in the Workspace, there is the possibility to replace it with the archived project. If the project is not present in your workspace, you are given two options for restoring a project. By restoring to the original file path, the files are put back in the original project folder. By restoring to a new file path, the project is replicated in a new location, which must be specified. Click *Restore*.

### 4.7.4 Project visibility

When there are several projects in the Workspace, it can be useful to temporarily hide some of them. Select *Display* in the *Project* icon drop-down list and check (or uncheck) the Visibility box in front of the projects that should be hidden (or shown).

#### 4.7.5 Project properties

To quickly view details (creator, comments, etc.) about a project, right-click on the project and select *Properties* in the contextual menu.

You can also modify the Project Name and Comment.

#### 5.1 Introduction

This chapter presents operations that specifically relate to gel images. You will learn how to manipulate images, use different tools to view signal intensities, and discover ways to visually compare images. The last section explains how to save, export and print gel images.

### 5.2 Manipulate images

The main toolbar provides the following tools to deal with images:



Select this tool, click in the image and hold down the left mouse button while moving the cursor. The image changes position. Release the button at the position you want.

#### **Alternatives**

There are other ways to change positions in gels:

- Use the View: Gels: Navigation: Move menu.
- Use the shortcut keys for the above-mentioned menu commands.
- Use the scrollbars.

#### Move all gels

To move all gels in the current sheet in the same way, hold down the Shift key while changing the position in one of the gels.



Select this tool, click repeatedly in the area of the gel where you want to zoom in. Right-click repeatedly on the gel to zoom out.

You can also define a zoom area: place the cursor at the top left corner of the area, hold down the left mouse button, and move to the bottom right position (a red box is displayed). Release the mouse button at the end point.

#### Apply to all gels

To move all gels in the current sheet to the same position with the same zoom factor, hold down the Shift key while zooming in or out on one of the gels.

#### **Alternatives**


There are other ways to zoom gels:

- Use the mouse scroll wheel to zoom in or out.
- Use the **View**: **Gels**: **Navigation**: **Zoom** menu.
- Use the shortcut keys for the above-mentioned menu commands.
- Use the scrollbars.
- To temporarily enlarge an area in a gel image, hold down the Ctrl key while the **Zoom** tool is activated. The area under the cursor is enlarged as if you were looking through a magnifying glass.

# Overview option

When you zoom in on a gel, it can be helpful to have an overview of where the region is localized on the full image (Figure 5-1). This overview enables you to easily locate and move to any region you want on your gel.

Choose *View:* Global: Show Overview to show or hide the overview of each image in its lower right corner. The green rectangle in the overview corresponds to the current view of the gel. You can drag the green rectangle to another position to display a new region.



**Figure 5-1.** Overview option activated. In the lower right corner of all images in the Display Zone is a small overview of the entire gel with a green rectangle corresponding to the visible gel area.

### 5.2.3 Same position and zoom factor

To move all gels in the current sheet to the same position with the same zoom factor, double-click on one of the gels. The corresponding position in the different gels is estimated by interpolating between the surrounding matches, or if no matches exist, between the two nearest landmarks. Finally, when no landmarks exist, the gels are aligned at the same location using the X and Y coordinates.


**Note:** For more information about landmarks, see 7.3.

# 5.2.4 Region

A region is a rectangular area in an image that is of interest for displaying a 3D View, for previewing spot detection parameters or adjust contrast settings, for cropping, etc.

#### 5.2 Manipulate images

Select this tool, place the cursor at the top left corner of the area you want to define, hold down the left mouse button, and move to the bottom right position (a dashed box is displayed). Release the mouse button at the end point (Figure 5-2).



**Figure 5-2.** Region tool. The left gel shows a region after selection, the middle gel during selection, and the right gel shows when the region has been reduced to its minimum size.

#### Edit region

You can move a region by clicking inside the box and dragging it. You can also change the size of the box by dragging a corner or edge. To remove a region, double-click on the gel. If the box is reduced to its minimum size, then its appearance changes to a blue circle.

#### Apply to all gels

To define the same region on all gels in the current sheet, hold down the Shift key while drawing the box on one of the gels or while clicking in an existing region.

## 5.2.5 Measure

Select this tool to measure pixel, pI, MW, or real world (centimeter or inch) distances between two pixels in an image. Click on a pixel (for example, center of first spot), hold down the left mouse button and move to the next pixel (for example, center of second spot). The horizontal and vertical distances between the start and end points are displayed (Figure 5-3).

The units of the displayed coordinates can be changed by choosing **Tools**: **Options** and looking under the **Display** tab.




Figure 5-3. Measured distance between pixels on an image. In the example, the horizontal and vertical distances between the starting (left) and end (right) points are 2.51 cm and 1.11 cm, respectively.

#### 5.2.6 **Bookmarks**

You can bookmark a view of your images in a sheet. This saves the currently visible areas. Bookmarks allow you to return to the same area on your gels at a later time or to show a similar area on the same images opened in a different sheet.



### 🔀 Create

Choose View: Sheet: Bookmarks: Create and enter a name.



To go back to previously visible areas, choose View: Sheet: Bookmarks: Load and select the name of the corresponding bookmark. Select the name of another sheet in the list to reproduce the visible areas of that sheet in the current sheet.



To remove a bookmark from the list, choose View: Sheet: Bookmarks: Delete.

#### 5.3 View signal intensity

The digitized image is composed of individual pixels, each of which is characterized by its X and Y coordinates, and its signal intensity (raw pixel value). This section describes different approaches to explore the signal intensity and adapt the way it is visualized.

#### • Adjust contrast 5.3.1

Sometimes the gray levels displayed by default are so low that small spots are hardly visible. To emphasize these very faint spots, you can adjust the contrast of the image and/or display images using pseudo colors.

#### 5.3 View signal intensity

Choose *View: Gels: Adjust Contrast* and select a region of interest on one or more images using the *Region* tool. These regions will allow you to preview the contrast and pseudo color settings before applying them to the selected images. The size and position of the preview regions can be adjusted at any time. After having adjusted the settings as described below, click *Apply*. The new contrast settings are saved with the image file.

**Note:** Any Adjust Contrast changes only influence how the image is displayed on your screen and do not affect the underlying data, spot detection, and quantitation.

### **Gray level mapping**

Some scanners are able to scan 2-DE images with 100000 gray levels or more. Because common computer screens are only able to display 256 gray levels, mapping must be undertaken between the image gray levels and the 256 screen gray levels. By default, the software uses a linear mapping function where the lightest point in the image is mapped to 0 (white) and the darkest point is mapped to 255 (black).

This is illustrated in Figure 5-4. The histogram displays the frequency with which each gray level (from 3524 to 28530) occurs in the first image of the sheet. The low gray levels (on the left of the histogram) corresponding to the background occur very often, whereas the high gray levels (on the right of the histogram) corresponding to the darkest spots are much less frequent. The red line indicates that the minimum and maximum gray levels (3524 and 28530, respectively) are remapped by the default linear mapping. The vertical axis for the red remapping graph corresponds to the 256 screen gray levels.





Figure 5-4. Adjust contrast function.

#### Change the default mapping

The software offers two complementary ways to change the default mapping function in order to improve the visual display of the gels (Figure 5-5):

You can define the minimum and maximum gray levels (that is, look at only the light or the dark regions in the images). Do this by decreasing the size of the gray level range that is to be remapped linearly to the screen gray levels. Move the left or right borders of the slider that is found below the histogram function. Once the size (interval) of the slider is decreased, you can also move the interval to the left or right. Alternatively, you can type valid numbers in the boxes at the lower left and right corners of the histogram.



**Figure 5-5.** Adjustments to the gray level mapping for the selected image are immediately reflected in the preview region. In this example, the maximum gray level was set to 13300 by moving the right side of the slider to the left. All pixels with a gray level higher than 13300 appear as black. By setting the Bending slider to -3, the image becomes darker.

You can additionally use a non-linear mapping function. The Bending
parameter, that is, slider to the left of the histogram, expands or
compresses the contrast at the dark or light ends of the range. When the
bending parameter is positive, the image is lighter. When the bending
parameter is negative, the image is darker.

#### **Contrast**

Select the image for which you want to see the current minimum and maximum gray level settings. When you make changes to these settings, the list will display *Modified* to reflect this fact.

#### Unit

Two different units can be used for displaying the gray level minimum and maximum:

- Intensity uses the raw pixel values as displayed in all the reports. When a
  calibration is done, these correspond to the calibrated pixel values.
- %: Chooses the scale as a percentage of the total gray level range in the histogram.

#### Region

You can display the histogram for the gray levels in the:

- Whole Image. This is the default option.
- Selected Region. The software only considers the gray levels that are
  present in the selected region. This option is practical in combination with
  the % unit and a relatively small region. In this case, you enter an adaptive
  mode that allows you to adjust to the local gray levels. The effect is a
  regional increase in contrast that is useful for viewing very faint spots.

#### Colors

The software offers various color palettes to visualize the intensities in your image. The Gray+Saturation palette is helpful when you want to visualize saturated spots or the background. It corresponds to the default Gray option, except for the maximum gray value (saturation) that is displayed in red, and the minimum value (background) that appears in blue. To decrease the stringency on what is considered saturation or background, you can modify the minimum and maximum gray levels.

#### Invert

You can inverse the gray levels by checking the *Invert* box.

## 5.3.2 A 3D view

Another way to examine the intensity variations in gel images is by looking at the three-dimensional view of gel regions (Figure 5-6). In this type of view, the X and Y axes represent the pl and MW values, whereas the pixel intensity is plotted along the third dimension (Z axis). The resulting image shows a peak for each protein spot, with a peak height that is proportional to the spot intensity. 3D views can be rotated in any direction to look at the interesting spots from all sides, thus facilitating spot editing or matching decisions.

The 3D View window can be displayed by choosing *Reports: 3D View*. If spots are selected (see 6.4.1), the area around these spots is visualized for all images in the sheet (Figure 5-7). If only regions were defined on images in a sheet, these regions are shown in the 3D View (Figure 5-6). If both spots are selected and regions are defined, the area around the spots is shown by default but you can choose to display the regions.

The layout of the images in the 3D View reflects as much as possible the layout of the images in the sheet. Bold black lines separate images that are in different


### 5.3 View signal intensity

panes. The pixel gray levels are overlaid on 3D views; the top of the peaks are therefore always darker. You can restrict the number of images displayed in the 3D View by setting their *Visibility* in the *Display Options* (see below).

By default, 3D views can be rotated by dragging the mouse while holding down the left button. If your mouse has a scroll wheel, it can be used to zoom in or out on the view. You can change to one of the other image manipulation tools avaible in the *Tools* icon drop-down list (see below for more details).

**Note:** All views in the 3D View window are manipulated simultaneously.

Right-clicking is reserved for the contextual menu. It allows you to select the spot where your cursor is positioned and to quickly access one of the image manipulation tools described below.



**Figure 5-6.** 3D views of the regions defined in A\_T1\_Gel1 and A\_T2\_Gel2.




Figure 5-7. 3D views for the selected spots.

The following tools are available in the 3D View toolbar:



To manipulate 3D views with your left mouse button, you can switch to one of the following options:

- **Rotation**: Rotate the images using your mouse.
- **Zoom/Contrast**: Zoom in or out by dragging the mouse horizontally. Adjust the contrast, or height of the spots, by dragging the mouse vertically on the view.
- *Translation*: Choose one of the suboptions (Z, Y, X, XY) to move the views along the desired axis.
- **Auto Rotate**: Let the images rotate automatically around the Z axis so that you can view the spots from all sides.
- Show Default View: Return to the original view.



- Stack: Display the images one on top of the other, instead of using the
  default tiled view.
- Next/Previous: Move to the next or previous image in the stack.
- Animate: Switch automatically between the different images in the stack.
   This option is useful for visualizing the expression variations in a set of images.
- Transparent Mode: Display one of the images (Reference) with a transparent surface. This option is useful for visualizing expression variations in a single, static view. To move one view with respect to the other, choose Translation in the Tools icon. Hold down the Ctrl key (gray) or the Shift key (green for the Reference) while dragging the view.
- **Set Reference**: Select the image to be used as the Reference for the Transparent Mode.
- *Transparency Settings*: Choose the color for the display of the Reference, and increase or decrease the transparency.


## Display options

- Visibility: Select the images to be displayed in the 3D View window.
- Spot Shape: Set the way the spots are displayed in the 3D View (Crossed, Outlined, Filled, None).
- Display Options: Display the X (pink), Y (purple) and Z (blue) axes, the
  coordinates of the point on which these axes are centered, and change
  the way the surfaces of the 3D Views are visualized (Grid, Wireframe,
  Smooth).
- Color Palette: Combine the 3D View with one of the color palettes available.

## 5.3.3 Profile

It is sometimes difficult to judge if spots should be split, are saturated or not, or have other problems such as so-called donut structures (low intensities in the center compared to the borders). It is important to identify such problems as they will lead to incorrect spot quantitation.

The Profile function, activated or deactivated by choosing *View: Global: Show Profile* can help in such cases. Red curves represent the intensity variations of the gel in the vertical and horizontal directions at the position of the mouse cursor (Figure 5-8). The Profile can clarify the intensity changes in a gel and assist in making editing decisions, without the need to open an additional window, as is the case with the *3D View* tool.



**Figure 5-8.** Profile feature. Red curves represent the intensity variations of the gel in the vertical (right) and horizontal (bottom) directions at the position of the mouse cursor. Green lines indicate the exact position of the cursor, whereas the numbers indicate the minimum and maximum gray levels in a specific profile view.

#### 

At any given time, the Cursor Information window can be used to display information on pixels such as the pixel intensity and the X and Y coordinates expressed in pixels, in pl and MW units (if available), and in cm or inches. Note that if spots are detected, the Cursor Information window also displays spot information

The Cursor Information window is available by choosing *View: Global: Cursor Information*, or by clicking the corresponding icon in the Display toolbar. Place the cursor over the pixel for which you want to display information.

## 5.4 Visually compare images

When working in MatchSet or Classes sheets, the software provides different tools to visually compare images to the sheet reference.

# 5.4.1 Sheet reference

You can set the sheet reference by choosing *View: Sheet: Set Reference*. All other images in the sheet are compared to this reference image for the options in the *View: Sheet* menu (Align Images, Show Dual Color, Spot Overlap).

### 5.4.2 Purpose of aligning images

The Align Images feature facilitates the visual comparison of images that demonstrate large variations in protein migration. It is especially helpful when the images are in Stacked mode. It modifies (Figure 5-9) the images in the sheet so that they are better superimposed with the sheet reference and therefore with each other. This allows easy identification of corresponding spots.

The Align Images feature is a purely visual tool. It is not used in, and will not improve, the matching process. The only reason to align your images is to ease their visual comparison (possibly in combination with the Dual Color mode).


# 5.4.3 Align images

To align images, the software needs to know which positions in the different images correspond to each other, that is, represent the same protein form. This is done by defining landmarks as described in 7.3. Landmarks should be added gradually so that you can monitor their individual influence on the alignment (automatically updated after each landmark addition). Immediately remove landmarks that decrease the alignment quality. The alignment algorithm then deforms the images to superimpose the landmarks.

**Note:** As opposed to matching, landmarks do not need to be linked with spots to carry out image alignment. Therefore, no spots need to be detected.

#### To align images in a sheet:

- 1 Set the desired image as the sheet reference, with View: Sheet: Set Reference.
- 2 Select the *Landmark* tool and define a few landmarks, bearing in mind the rules listed in 7.3.
- 3 Choose View: Sheet: Align Images.
- 4 If parts of your images are still not sufficiently aligned, you can add extra landmarks. The alignment is automatically updated.
- To view the original images, choose *View: Sheet: Align Images*. The original images replace the aligned ones.



**Figure 5-9.** Image (a) before and (b) after alignment. The match vectors and landmarks are displayed.

Alternatively, double-click in an image (with the **Move** tool selected) to locally superimpose the different images:

- If the images are matched, the image positions are synchronized based on the surrounding matches.
- If the images are not matched but have landmarks, the image positions are synchronized based on a simple interpolation between the two nearest landmarks.
- If no matches or landmarks are present, the image positions are synchronized based on the same image location (X,Y coordinates).

## 5.4.4 Show dual color

Choose *View:* Sheet: Show Dual Color to display each of the images (the sheet reference and the current image) in one of two colors: red and cyan (Figure 5-10). When the pixel colors of the two superimposed gels are added:

- Overlapping spots appear as shades of gray.
- Cyan spots are present only in the current gel.
- Red spots are present only in the sheet reference.
- Halos of cyan or red around dark spots indicate that the protein is over or under expressed, respectively, compared to the sheet reference.

The less color you see in the dual color mode, the more similar the gels. Of course, this is only true if the gels are correctly aligned (see above) and superimposed.

Note: It is recommended to perform any operations related to gel alignment

(especially the definition of landmarks), before entering the Dual Color mode to avoid a slow down due to the recalculation of the overlaid images.




Figure 5-10. Dual color view.

## 5.4.5 Spots overlapped

Chooose *View: Sheet: Spot Overlap* once spots have been detected and the visible spots on the sheet reference will be shown in blue on the image (Figure 5-11). Thus, you can easily compare the position and size of the red spots in the current gels with the blue spots on the sheet reference (shown in Crossed or Outlined mode).

You can change the default color to be used for overlapped spots by selecting the *Display* tab in *Tools : Options*, and clicking the *Overlapped* box.




Figure 5-11. Spots overlapped with (a) outlined and (b) crossed spot shapes.

### 5.5 Grid lines

# 5.5.1 Display grid lines

Choose *View: Global: Grid Lines: Show* to display grid lines over your images (Figure 5-12). Grid lines can be used to evaluate distances between spots, in terms of pixel coordinates, pI/MW units, or length units (cm or inch). Grid lines are also a helpful way to visualize deformations in aligned gels because the grid lines are warped in the same way as the image.

# 5.5.2 Edit grid lines

Choose *View: Global: Grid Lines: Edit* to change the grid properties. The software partitions the visible area into the *Number of subdivisions* entered by the user (Figure 5-12). The graduations can be attached to the *Gel* or the *Screen*. This means that the software divides the gel width/height or the visible screen area by the number of subdivisions. The *Coordinate Units* can be Centimeters, Inches, Pixels, or pl/MW units, provided data is available in the annotation category pl\_MW. The pl/MW grid can also be displayed over gels that do not contain this information, but that were matched to a gel having pl and MW values.




Figure 5-12. Grid lines in cm attached to the gel image.

## 5.6 Gel reports

# 5.6.1 Gel table

Choose *Reports: Gel Table* to display a table (Figure 5-13) with summarized information about the images in the current sheet. The available columns are:

- File name, ID, path, size, modification and creation dates.
- Image height and width expressed in pixels (Rows and Columns), pixel size (PixSize).
- Minimum and maximum gray levels before (MinGray and MaxGray) and after calibration (MinValue and MaxValue).
- Calibration information (Calibration Formula, Unit, Name, Creator and Date).

- Number of spots and annotations.
- Staining method and user-defined gel descriptions.
- MatchSet to which the image belongs and the Class to which it was assigned.




Figure 5-13. Gel Table.

You can customize the Gel Table to display only the columns of interest. Four predefined report templates are already available from the *Load* icon in the *Settings* of the Gel Table toolbar:

- Properties: Shows various properties of the gel images: File Name, image height and width in pixels, maximum gray levels before and after calibration, Modification Date, as well as the MatchSet and Class to which it belongs.
- *Files*: Shows information on the image files: File Name, ID, Path, Size, Creation and Modification Dates, as well as the MatchSet and Class to which they belong.
- Descriptions: Shows the user-defined gel descriptions.
- *Calibration*: Provides calibration-related information: File Name, ID and Path, as well as Calibration Formula, Name, Creator and Date.

## 5.7 Save, export and print images

#### 5.7.1 Save

The software automatically saves your images and all associated data as part of the project when you close a sheet or exit the software. You can also save your work by choosing *File: Save*.

#### 5.7.2 Export

You can either export images to a file or to the clipboard for direct pasting into another software. The procedure is very similar to exporting images to files.

#### 5.7 Save, export and print images

However, you can only export one image at a time to the clipboard. Finally, you can also export a view of the current sheet.

The export options are available in the *File: Export* menu:

- 🕅 Image to Clipboard
- 🔯 Image to File
- Sheet to Clipboard
- Sheet to File

Rather than saving your gel images in the ImageMaster file format, you may want to export them to a different file format (TIFF, BMP or PNG). The gel images are exported as 8-bit, flat, rasterized images without any structure. This means that gel components such as spots and annotations are saved exactly as they appear on the screen, but are no longer recognizable as ImageMaster objects and therefore become part of the image. Consequently, exported gel images should only be used for presentation purposes and not for further analysis with any software package.

Adapt the size of the exported image by zooming in or out. If a region is selected on the image, you are given the choice to export the entire image or the selected region only.

#### 5.7.3 Print images and sheets

#### **Print options**

The software provides various printing options in the *File: Print* menu:

- Images: Print selected images or image areas (if defined with the Region tool). This option prints one image per page.
- Sheet: Print the current sheet.

Whatever your choice, the image is printed as it is displayed on the screen retaining objects and properties such as spots, annotations, contrast mapping and pseudo colors, alignment, zoom, grid, etc.

**Note:** With a zoom factor of 1, the printed gel image takes the full paper width. You can adapt the zoom factor to decrease the size of the printed image.

#### Page Setup

You can change print parameters such as printer name, paper size, paper orientation, etc. Choose *File: Page Setup*. This command opens the standard print window where printer-related settings can be modified.

#### 6.1 Introduction

Once gels have been added to a project and you have taken a good look at them, you are ready to detect spots. A spot delineates a small region in the gel where protein is present. This shape is automatically differentiated by a spot detection algorithm and quantified; its intensity, area and volume are computed.

There are two different spot detection algorithms implemented:

- Non-DIGE images are detected with the ImageMaster algorithm.
- DIGE images are co-detected with GE Healthcare's DeCyder 2D algorithm.

## 6.2 Detect spots in non-DIGE gels

#### 6.2.1 Procedure

The ImageMaster spot detection algorithm is optimized to give relevant biological results with minimum user interaction. You can preview spot detection and fine-tune a few parameters before automatically locating all the spots in your image.

#### To detect spots automatically:

- Display the gels to be detected by selecting a match set in the Workspace. Right-click and select **Display** in the contextual menu.
- 2 Choose *View*: *Spots*: *Outlined* to visualize the spot borders.
- 3 Click the *Region* tool. Draw preview regions on one or more gels. Note that you can still draw, resize or move regions while setting the detection parameters.
- 4 Select the gel images for spot detection.
- 5 Choose *Edit*: *Spots*: *Detect*.
- The Detect Spots window appears on-screen and the spots in the preview regions of the selected images are detected with the default parameters. If you do not want to recalculate the spots in the preview regions for each parameter change, turn the **Auto Preview** option off. To manually refresh the preview regions, simply click the **Preview** button.
- 7 Adjust the detection parameters (Figure 6-1). In particular optimize **Smooth** to detect all real spots and split the overlapping ones. Subsequently, filter

- out the noise by changing the **Saliency** and **Min Area** values. See below for more details on spot detection parameters.
- 8 When you are satisfied with the preview, click **OK** to detect all spots in the selected gels using the specified parameter values. Note that you can still change your gel selection at this point.
- 9 The spot shapes are displayed on the images.




Figure 6-1. Adjusting spot detection parameters in real time.

#### 6.2.2 Spot detection parameters

Spot detection parameters are best adjusted in the following order:

- Smooth: Set this parameter first. It fixes the number of times ImageMaster smooths the image before detecting spots, using a smooth-by-diffusion algorithm. The Smooth parameter should be optimized to detect all real spots and split as many overlapping spots as possible without being concerned about noise spots (these can be filtered out with the Saliency and Min Area parameters).
- **Saliency**: This parameter is a measure based on the spot curvature. It indicates how far a spot stands out with respect to its environment. Real spots generally have high saliency values whereas artifacts and

background noise have small saliencies. Although the Saliency is an efficient quantity for filtering spots, it is also highly dependent on the images (for example, image resolution and depth). Some gels need a saliency value of 10 for correct filtering. Others may necessitate a value of 5000. To estimate the saliency range to use with your images, you can display the Cursor Information window or Spot Table and look at the saliency value given for a spot that you would like to suppress. Enter this value in the Saliency field in the Detect Spots window. The spot detection algorithm then discards all spots with saliencies smaller than the specified threshold

 Min Area: After setting an appropriate Saliency to filter out all noise spots, there may still be noise in your gel that cannot be eliminated without suppressing real spots. This often happens with dust particles that consist of a few very dark pixels. Get rid of these artifacts by using the Min Area parameter. It eliminates spots that have an area smaller than the specified threshold (expressed in number of pixels).

## 6.2.3 Spot quantification

The software automatically computes the amount of protein present in each spot. Figure 6-2 illustrates the principles of spot quantification in the ImageMaster algorithm. Measuring the protein quantification values in this way

Intensity

Gels: 1 Spots: 0 Annotations: 0

protein expression variations (relative quantification).

File View Edit Show Select Analyze Reports Tools Window Help

Gals (Gals)

Gals


Gals (Gals)

Gals

Gals (Gals)

Gal

has the advantage of being more robust and reproducible when calculating protein expression variations (relative quantification).



- Intensity: The software first calculates the intensity of a spot. The intensity is based on the highest calibrated pixel intensities in the spot from which the background has been withdrawn. The background is defined as the minimum pixel value in the spot neighborhood.
- Area: The area of a spot is not determined at the spot base because the
  base is often arbitrary and difficult to determine. ImageMaster computes
  the area at 75% of the spot intensity, as measured from the peak of the
  spot. The spot outlines displayed in ImageMaster exactly encircle this
  computed spot area (expressed in mm²).

B

- Vol: The volume of a spot is calculated as the volume above the spot outline, which is situated at 75% of the spot height (as measured from the peak of the spot). In Figure 6-2, the measured volume of the spot is striped. Please note that the volume values, like the intensities, depend on pixel intensity calibration.
- %Vol: The relative volume of a spot is calculated as indicated below. It is a
  normalized value that remains relatively independent of variations due to
  protein loading and staining by considering the total volume over all the
  spots in the image. This means that in an image where, globally, the spots
  are darker than in another image, the majority of spot volumes is higher.
  However, the bulk of %Vol should be similar to those in the compared
  image, at least for gels with similar spot patterns.

%Vol = 
$$\frac{\text{Vol}}{\sum_{S=1}^{n} \text{Vol}_{S}} \times 100$$
  
where Vol<sub>S</sub> is the volume of spot S in a gel containing n spots.

**Note:** When you look at a detected spot in the 3D View, you will notice that the borders are not localized at the base of the spot, especially for intense spots. This is normal because the spot contours displayed effectively reflect the quantified parts of the spot, and they do not correspond to the 'whole' spot, which is difficult to define.

## 6.3 Co-detect spots in DIGE gels

The co-detection algorithm is designed to simultaneously process one, two or three images derived from a single gel.

- Single detection: one image.
- Double detection: two images.
- Triple detection: three images.

Single detection is performed on images of fluorescently post-stained gels used for picking, a case where there is a single image associated with the gel.

Double and triple detection takes advantage of the inherent co-migration benefits of the CyDye™ DIGE Fluor dyes. A set of co-run images are merged together to include all spot features in a single image. Spot detection and spot boundary definition is then performed using pixel data from all the individual raw images and the merged image. The resultant spot map is overlaid back onto the original image files. Since the spot boundaries and the detection areas are identical for all images, the spots are effectively already matched. This process results in highly accurate volume ratio calculations.

#### 6.3.1 Procedure

#### To perform spot detection on DIGE images:

- Display the gels to be detected by selecting a match set in the Workspace. Right-click and select *Display* in the contextual menu.
- 2 If desired, select a subset of gels to be detected.
- 3 Choose *Edit*: *Spots*: *Detect*.
- 4 In the DIGE Spot Detection window, enter an estimation (see below) of the *Number of Spots* present in the images. Click *OK*.

#### 6.3.2 Spot detection parameter

When detecting DIGE images, you must enter an estimation of the Number of Spots present in the images. It is recommended that this value be overestimated to compensate for the detection of non-protein spots on the image, for example, dust particles which can subsequently be excluded from the analysis using spot filtering.

If all the spots are not identified, the spot detection process can be repeated with a higher number of estimated spots.

For example, for a mammalian lysate run on an 24 cm pH 4-7 Immobiline<sup>™</sup> DryStrip and a large format gel, such as the Ettan<sup>™</sup> DALT Gel (20 cm x 26 cm), a value of 2500 for the estimated Number of Spots should be satisfactory.

#### 6.3.3 Spot quantification

Volume, area, intensity, slope, and volume ratio for individual spots are automatically calculated and included in the Spot Table and Cursor Information window.

- **Vol**: Spot volumes (sum of pixel intensities within the spot boundary) are always expressed with background subtracted. Background is subtracted on a spot specific basis by excluding the lowest 10th percentile pixel value on the spot boundary from all other pixel values within the spot boundary. The spot volume is the sum of these corrected values.
- Vol ratio: Volume ratios (volume of current image spot/volume of DIGE reference image spot) indicate the change in spot volume between two images.


**Note:** When using single detection the volume ratio value is 1.0 for all spots since there is no second image.

### 6.3.4 Exclude spots

To exclude small, non-protein spots on DIGE images, it is recommended to filter the spots based on the Volume. More precisely, spots that have a maximum volume in the two or three images that is lower than a given threshold can be deleted.

#### To exclude DIGE spots based on maximum volume:

- 1 Choose **Reports: Spot Table**.
- 2 Click the **Select by Value** icon in the toolbar of the table.
- 3 Select Vol in the displayed list. Make sure the >= box is checked and enter a value to be used as the cutoff for the volume. Deselect the <= box and click OK (Figure 6-3).
- 4 All spots with a volume higher than the given threshold are selected, together with their matched spots (which may have volumes that are smaller than the threshold). These are the spots that are kept.
- 5 Choose **Select: Spots: Inverse Selection**.
- The current selection includes all spots for which the maximum volume in any of the two or three DIGE images is lower than the given threshold.
- 7 Choose *Edit*: *Spots*: *Delete* to permanently delete the selected spots from your DIGE gels.



**Figure 6-3.** Select by Value can be used to filter spots based on their volume.

**Note:** Instead of deleting spots, you can exclude them from the analysis by disabling them, as described in 6.4.3.

## 6.4 Select spots

## 6.4.1 Select

Spots can be selected with the **Select** tool. Once selected, they are highlighted in green, unless the default spot colors were modified. All matched spots are selected as well.

If an annotation is attached to a spot, the annotation is also selected. Similarly, if you select an annotation or label with the **Select** tool, the linked spot is also selected.

To select more than one spot, select the first one and then hold down the Shift key while clicking on additional spots.

To select all spots in a region, place the cursor at the top left position of the desired region, hold down the left mouse button, and then drag the cursor to the bottom right position. All spots in the designated region are selected.

To deselect all spots, click in the gel (not on a spot).

#### **Surrounding box**

By default, spot selections are surrounded by green boxes. This makes it is easier to localize selected spots, especially when working at low zoom factors. To deactivate this option, choose *Tools: Options*, and uncheck the corresponding box in the *Display* tab.

#### **Alternatives**

There are other ways to select spots:

- Choose a command in the Select menu.
- Select spots in reports.

#### 6.4.2 Spot sets

It is possible to focus your analysis on particular spots by creating and saving spot sets for later selection or combination.



#### Create

Select spots you want to include in a spot set, either manually or by selection in a report. Then choose *Edit: Spot Sets: Create*, enter a name for the new spot set and click *OK*.

Spot sets can be visualized as columns in various reports (Figure 6-4). If they are not displayed by default, you can add them by checking the corresponding box in the **Settings** of the report window (only for tabular reports). Once the column is displayed, you will see a checked box for spots that belong to the set, or an empty box for spots that do not belong to the set. Click in a box to change its

state. If several spots are selected, you can change the state by clicking in one box while holding down the Shift key.



**Figure 6-4.** Spot sets Anova p<0.001 and Ratio: 2 in a Class Analysis Table. New spot sets can be created by clicking on the Create Spot Set icon.

Spot sets can also be created by clicking the *Create Spot Set* icon in a report toolbar. All currently selected spots are automatically included in the newly created spot set.



Delete a spot set by choosing *Edit: Spot Sets: Delete* and selecting the spot set(s) to be deleted.



You can combine two spot sets using logical operators, in order to create a new spot set or change the current selection (Figure 6-5). Four operators are available:

- And: Keeps spots that belong to both spot sets.
- **Or**: Keeps spots that belong to either one or both spot sets.
- **Not equal**: Keeps spots that belong to only one of the two spot sets.
- **Exclude**: Keeps spots that belong to the first spot set and do not belong to the second spot set.




Figure 6-5. Combine Spot Sets window.

#### To combine two spot sets:

- Choose Edit: Spot Sets: Combine. 1
- 2 From the *Inputs* drop-down lists, select the spot sets to be combined. The current selection can also be used as input.
- 3 Specify the operator (And, Or, Not Equal or Exclude).
- Enter a new name for the resulting spot set. Alternatively, the output can be immediately reflected in the current selection. Click **OK**.

#### Select

To select spot sets, choose Select: Spot Sets.

#### 6.4.3 **Enabled spots**

By default, all spots are enabled and therefore represented in the reports. To exclude a specific subset of spots, you can disable protein spots that are not of interest, or specifically define a set of spots to be enabled. Only enabled spots appear in reports.

The options related to the creation, saving or loading of sets of enabled spots can be chosen from the **Edit**: **Enabled Spots** menu:



Select the spots you want to focus your analysis on and choose this option to enable the selected spots. After deselection, excluded spots are disabled and appear in yellow.



Use this option to add selected spots to the current set of enabled spots.



## Remove

Use this option to remove selected spots from the current set of enabled spots.



Use this option to save the currently enabled spots as a new spot set. Enter a name and click **OK**.



Use this option to enable spots belonging to an existing spot set.

#### Select

To select the enabled spots, choose **Select: All Enabled Spots**.

## 6.5 Display spots

#### 6.5.1





## Spot shape

Once spots are detected, you can choose how to display their shapes (outlined, crossed, filled, outline/filled) on the gels from the *View: Spots* menu.

#### 6.5.2 Spot color

By default, enabled spots are displayed in red, disabled spots in yellow, selected spots in green and overlapped spots in blue. You can change these default colors in the *Display* tab in the Options window (accessible by choosing *Tools: Options*). Click the colored box you want to change and the Color window opens. Choose the preferred color from the spectrum and click *OK*.

## 6.6 Edit spots

**Note:** Except for deleting spots, spot editing is not allowed on DIGE gels.

Quantitative protein data, and in particular the spot volume, are highly dependent on an optimal and reproducible definition of the spot borders and a correct splitting of partially overlapped spots. To guarantee reproducibility of quantitative work it is therefore recommended to create spots by using the automatic spot detection algorithm in the software and to avoid manual editing as much as possible.

However, spot detection differences can still occur. In particular, some spots are differently split in gels to be compared. The software offers the following solutions to deal with detection variations between gels without calling for spot editing:

- Create multiple matches. In practice, this means that you can match "composite spots" that are treated as unique entities in the quantitation.
- Propagate all or selected spots from one image to the other images.

**Note:** Both solutions require prior matching. Therefore, before doing any spot editing, first match your images and only then consider to use one of the options described below

# 6.6.1 Manual editing

You must enter the special spot-editing mode to manually edit spots. Choose *Edit: Spots: Edit Enabled*, or click the corresponding icon in the Detect And Match Spots toolbar. The Edit Spots toolbar displays. Selecting the Edit Enabled option again disables spot editing.

**Note:** For edited spots, the Saliency value becomes zero. This can be used to auickly check which spots have been edited.

The following spot editing tools are available:



Click this icon and draw the outline of the new spot. Alternatively, double-click the desired location in your image, and set the disc radius for the circular spot to be drawn in the *Create Spot* window.



Click this icon, select the spot to be deleted, and confirm.



Click this icon, select a spot to be split, and draw a line through the spot at the position where the separation should occur. Make sure you start and finish outside the spot.



Click this icon, select two or more spots to be merged, and draw a trajectory through the selected spots. Make sure to start and finish in the same spot.



Click this icon, select a spot to be grown, and outline the area you would like to add. Make sure to start and finish within the selected spot.



Click this icon, select a spot to be shrunk, and outline the area you would like to suppress. Make sure to start and finish outside the selected spot.

#### 6.6.2 Composite spots

As shown in Figure 6-6, the software allows you to match several spots in one gel with multiple spots in other gels. In the figure, the three selected spots in gel  $A_T1_Gel1$  are matched to the single green spot in  $A_T2_Gel1$  and the two selected spots in  $B_T1_Gel1$ , and so on.

Once the match has been effectively created, the three spots in A\_T1\_Gel1, for example, are treated as a single entity in the quantification. The quantification value for A\_T1\_Gel1 displayed in the different reports is the one obtained after combination of the values from the three individual spots.

This is an efficient solution for dealing with spot detection differences without subjective and time-consuming spot editing. When selecting a spot on a gel, any matched spots on the other images are automatically selected as well. Select all

spots to be matched and click the **Add Match** icon in the Detect And Match Spots toolbar.

Figure 6-6. Composite spots can be defined through multiple matches

### 6.6.3 Propagate spots

You can propagate all or selected spots from one image to the other matched images (Figure 6-7). This allows you to quantify identical areas on all gels.

- **For matched spots**: the spot in the destination image is replaced with the shape of the spot in the source image.
- For non-matched spots: the spot from the source image is copied to the
  equivalent location in the destination image. This position is extrapolated
  from the surrounding match vectors.

Matches are automatically created between the original and the propagated spots.




Figure 6-7. All spots from gel A\_T1\_Gel1 have been propagated to the other images.

#### To propagate spots from one matched image to another:

- 1 Select spots on one image.
- 2 Choose Edit: Spots: Propagate.
- 3 Select one or more images you want to copy the spots to.
- 4 The new spots are added to the selected images.

# 6.7 MW and pl calibration

If you have a gel with pI/MW standards, the software can compute approximate pI and MW values for all the spots/pixels in this image, as well as any other images matched to it.

Define pl\_MW annotations for a certain number of spots/pixels in the gel. The calculated pl and MW values for all spots in this gel and any matched gels are

automatically available in the Spot Report, DIGE Report or Cursor Information window. In addition, pl and MW grid lines can be displayed over the images.

#### To define pl\_MW annotations on an image:

- Select the image for which you know the pl and MW values for several protein spots. Spots may or may not have already been detected in this image.
- 2 Click **Select** in the toolbar.
- 3 Double-click on a spot (pixel) for which you know the pl and/or MW values.
- 4 Select the **pI\_MW** category in the Create Annotation by Click window.
- 5 Enter the known pl and MW values, respectively, separated by a space. Replacing one of the values with -1 means that no value is set.
- 6 Do this for a sufficient number of protein spots that are distributed over the whole image area. Obviously, the more spots and annotations, the better the approximated pI and MW values will be.

ImageMaster does the following to calculate approximate pI and MW values. In the case of pI, it looks up the two closest annotations to the left and to the right of the spot for which the pI will be determined and then interpolates between these two points. Since ImageMaster does not have any information about the experimental (possibly non-linear) pI scale, the calculated values are only approximate. In the case of MW of the spots, the procedure is similar, except that ImageMaster searches for the closest spots above and below the spot for which the MW will be determined and it makes a logarithmic interpolation.

Extrapolating pI and MW values is more complicated. For example, if the pI of a spot on the extreme right side of your gel is to be determined, the software looks for the two closest spots to the left of the spot in question. If these two spots are sufficiently distant from each other (in order to decrease the error), the value for the spot in question can then be extrapolated.

Normally, the pI and MW values in the Spot Report or Cursor Information window should be the same as in the defined pI\_MW annotations. However, this is only the case if the annotations are attached to actual spots and not to pixel positions in the image. If an annotation is attached to a pixel, the pI and MW values for the spot that lies closest to it will be slightly different from that of the pixel (to which the annotation is attached). You can solve this ambiguity by linking the annotation to the spot.

## 6.8 Spot reports

# 6.8.1 Spot Table

The Spot Table (Figure 6-8), obtained by choosing *Reports: Spot Table*, displays summarized information about enabled spots:

- Name of the image on which they were detected.
- Spot ID (see below), and Match ID if the spot was matched.
- Coordinates of the spot's center of gravity (X and Y).
- Quantification values: Intensity, Area, Vol, and %Vol. Depending on the spot detection algorithm used, the Saliency, Vol Ratio, and Slope will also be given.
- Calculated pl and MW values, if pl\_MW annotations were defined on the image, or a matched image.
- All linked labels and spot sets.




Figure 6-8. Spot Table.

#### Spot ID

Each spot in a gel has a unique identifier, called the Spot ID. Spot IDs of deleted spots are not reused. ImageMaster attributes a new ID to each new spot. When a spot is split, the child spot for which the coordinates are closest to the parent spot keeps the existing spot ID, the other child spot gets a new ID. When two spots are merged, the resulting spot is attributed the ID of the initial spot that was closest to the new center of gravity.

#### 

The Cursor Information window is available from the menu *View: Global: Cursor Information Window* or by clicking the corresponding icon in the Display toolbar. It can be used at any time to display information on pixels and spots located at the position of your mouse cursor (Figure 6-9).

Information on the pixel under the cursor:

- Name of the image.
- Calibrated pixel intensity.
- X and Y coordinates, expressed in pixels, in pl and MW units (if available), and in cm or inches.

Information on the spot under the cursor:

- Spot ID.
- Quantification values: Intensity, Area, Vol, and %Vol. Depending on the spot detection algorithm used, the Saliency, Vol Ratio, and Slope are also provided.




Figure 6-9. Cursor Information window.

6 Spots6.8 Spot reports

#### 7.1 Introduction

Matching is a key operation in 2-DE image analysis. Basically, the image matching algorithm compares gel images to find matches between related spots, that is, spots representing the same protein in the gels. A match is therefore composed of spot n-tuples (S1, S2, ..., Sn) where S1 is a Spot ID in the first gel, ..., and Sn a Spot ID in the last gel.

The matching algorithm always starts from the reference image or match set, and looks for each spot in this reference, if corresponding spots in the other images are found. If a spot is absent from the reference, it cannot be matched automatically. However, if you have several match sets in your hierarchy, there is a good chance that the spot is on the reference of at least one of them. If so, the spot will turn up in the analysis. Subsequently, additional spots can be matched to it manually (Figure 7-1).

Matches are propagated at each level of the hierarchy. This means that once all match sets are effectively matched, spots from one gel can be directly compared with those in any of the other gels.

**Note:** A DIGE gel is an inherent match set for which the co-run images are automatically matched. To subsequently match different DIGE gels, proceed like any other match sets.



**Figure 7-1.** Reference and matching. The selected spot is absent from groups AT1 and AT2, including the sheet reference A\_T1\_Gel1. Since the spot is present in the references B\_T1\_Gel1 and B\_T2\_Gel1, it is matched to corresponding spots and turns up in the analysis (Gel and Class Analysis Tables).

## 7.2 Display a match hierarchy

In the example below, the match hierarchy AB was displayed by right-clicking on its name in the Workspace and selecting *Display*. In the resulting sheet, only the lowest level (AT1, AT2, BT1, BT2) and highest level (AB) items are specifically visualized, in the panes and sheet, respectively. To select intermediary levels (A, B), you can use the corresponding icon in the pane tab of the reference (for example, AT1) for that level (A).

At the same time, the red markers in these icons indicate what item in a given match level is used as the reference in the matching.




Figure 7-2. Match hierarchy AB displayed in a sheet. The match level A was selected by clicking on the corresponding icon in the pane tab of match set AT1.

#### 7.3 **Define landmarks**

ImageMaster is designed to match gels with minimum user intervention. Nevertheless, when the gels are very distorted or different, you may need to help the matching process by specifying a few landmarks. Landmarks are points that relate corresponding spots in each of the gels to be matched.

In some cases, no landmarks are required. More often, a single landmark is sufficient for quick and efficient matching. If the matching results are not satisfactory, you can repeat the automatic matching procedure using additional landmarks.

### 7.3.1 Rules for defining landmarks

The addition of suspect or badly positioned landmarks can worsen the alignment results. The following rules must be considered when defining landmarks:

- The number of landmarks should be kept to a minimum. There is no point in putting a landmark on each spot.
- Landmarks should be well distributed over the whole images (covering both the X and Y directions). To correct for local distortions, it is recommended to define landmarks around the distorted regions rather than within those regions.
- Landmarks should only be defined on spots that clearly represent the same protein form. Protein variants definitely should not be used as landmarks
- Landmarks should be placed on small, sharp spots (of similar area), rather than on large diffuse ones (which may differ considerably in size) to reduce the error in the position.
- When a spot is missing on a gel (sometimes happens to border spots), you should not put a landmark (that is, validate a landmark) in a hypothetical spot position. Missing landmarks are not an issue.

## 7.3.2 Define 😉 Landmark

Landmarks can be defined using the dedicated tool as described below.

#### To define landmarks:


- 1 Click the **Landmark** icon in the toolbar.
- 2 Place the mouse cursor over a known, well-defined spot in the first reference gel and click. A validated landmark symbol (bold orange circle) appears on the spot (Figure 7-3).
- In the other images, drag the non-validated symbols (green circle with orange plus sign) onto the corresponding spots. If the symbol is already on a good spot, double-click to validate it.
- 4 Repeat steps 2 and 3 to add more landmarks.

In certain gels, the symbols only become visible once the landmark has been validated in the reference. In the example, the landmark must be validated in the image B\_T2\_Gel1 before any symbols appear in the images B\_T2\_Gel2 and B\_T2\_Gel3.

Sometimes, you may want to move or zoom your images during the landmarking process. When you click the Move or Zoom tools, the orange landmark symbols disappear, and only labels with the landmark numbers are

left. When you click the *Landmark* tool again, the symbols are reactivated and you can continue defining landmarks.

To delete landmarks, you must delete the corresponding annotations. Choose **Select: Annotations: By Category** and select **Landmarks** in the list. Then choose **Edit: Annotations: Labels: Delete**.



**Figure 7-3.** Placing of landmarks. Note the difference between the validated and non-validated landmark symbols.


In general, it is good to validate one or two landmarks on all the images in the hierarchy, so that they are used for matching at all levels. If it turns out after the matching that the gels within the lowest hierarchies (for example, AT1, AT2, BT1 and BT2) have been properly matched, but that the higher level matches (for example, A and B) are not satisfactory, you can add additional landmarks to the higher levels only. Do this by validating the landmark only on the reference images.

## 7.4 Automatic matching

An entire hierarchy can be matched automatically.

## To match the gels in a match hierarchy:

- 1 Select the gels to be matched in the sheet. All the gels, even in a multilevel hierarchy, can be selected.
- 2 Choose **Edit**: **Matches**: **Match Gels**. Or click **Match Gels** in the toolbar.
- 3 If several match sets are selected, you are asked to choose the ones to be matched in the Match Gels window (Figure 7-4). Items preceded by a ~ sign still require matching. Use the Ctrl or Shift keys to make multiple selections. All match sets can be selected at the same time. Click **OK**.



**Figure 7-4.** Match Gels window. In this example, all seven match sets were selected for matching. The items preceded by a  $\sim$  sign are not yet matched.

Match vectors are displayed in blue. The vector pattern is proof of consistency. If there is a mismatch, the vector has a different length and/or orientation.

If the matching did not work properly, you can rematch a particular match level after having added additional landmarks to the appropriate images.

## 7.5 Select matches

### **7.5.1** Select

When you select a spot with the **Select** tool, the matched spots are automatically selected.

### **Alternatives**

There are other ways to select matches:

 Choose an option in the Select: Matches menu (All, Inverse Selection, Multiple Matches).

### 7.5.2 Match count

The software allows you to select spots present:

- In a certain number of gels: the *Match Count* in a *Gel Analysis Table* gives the number of gels in which the spot is present (detected and matched).
- In a certain number of classes: The Match Count in a Class Analysis Table gives the number of classes in which the spot is present (detected and matched)

Click the **Select By Value** icon in these tables to refine your selection. For example, select only spots that are present in at least X out of Y gels, or X out of Y classes.

## 7.6 Display matches

You can visualize matches by selecting them. The software also provides the following tools to display matches.


## 7.6.1 Show vectors

Right after matching, the software automatically displays the match vectors in blue. Vectors link the spots in a gel with the corresponding spots in the sheet reference. This sheet reference has a darker green gel name and should not be mistaken for the match reference, which has a red component. The sheet reference can be changed by choosing *View: Sheet: Set Reference*.

A blue upside down triangle on a spot indicates that the spot was matched to one or more spots in other gels, but not to a spot in the sheet reference (Figure 7-5).

A spot with a triangle means that the corresponding position in the sheet reference lies outside the visible area (Figure 7-5).

To hide the match vectors, or on the contrary, to display them when they are not visible, choose *View: Matches: Show Vectors*.



**Figure 7-5.** Different representations of matches. The selected spots in match set AT2 with a blue upside down triangle are matched, but not with a spot in the sheet reference. Spots with blue triangles, as the one indicated with a green arrow, are matched, but the corresponding position in the sheet reference lies outside the visible area.

To minimize the match vectors (after having moved or zoomed an image), select the *Move* tool and double-click in the image so that it is synchronized with the other images and the sheet reference.

You can change the default color for the match vectors by choosing **Tools**: **Options** and going to the **Display** tab.

### 7.6.2 Show ID

Choose *View: Matches: Show ID* to display Match IDs for selected spots on selected gels. Previously, the gels must be matched. To hide the Match IDs again, possibly in selected gels only, choose *View: Matches: Hide All ID*.

Sometimes it may be an advantage to manually add or delete matches after the automatic matching procedure. Use the options in the *Edit: Matches* menu, or the corresponding icons in the Detect and Match Spots toolbar.

Note that matches were created in a hierarchical manner. This means that spots may be properly matched in one match set level, but not in another. To delete your matches, it is therefore important to select the appropriate images/panes.



Select spots that were incorrectly matched, make sure the appropriate match level (gels/panes) is selected, and click the *Delete Match* icon to remove the spots from the match.




Select all spots to be matched and click the **Add Match** icon to add the spots to the match

While selecting spots, you may see that some spots were already matched. When selecting one of them, the others are automatically selected. Make sure that the existing matches are correct. If this is not the case, delete them before proceeding.

### Multiple matches

ImageMaster enables the creation of multiple matches (Figure 7-6). In contrast to a single match, where only a single spot is selected per gel, a multiple match implies that one or more spots in one image can be matched to several spots in other images.

All the spots from such a multiple match on a given gel image are considered as a single spot in the subsequent data analysis. The calculated quantification values for this composite spot reflect the size, intensity, and abundance of the combined spots. Therefore, this is a solution to avoid spot editing.



**Figure 7-6.** Multiple matches. The selected spot was correctly matched throughout match set A and within BT2. But due to spot splitting in B\_T1\_Gel1, it was not matched within match sets BT1 and B. Rather than merging the two spots in B\_T1\_Gel1, the spots can both be included in the match and treated as a single entity in the spot quantification.

### **Review matches**

To review matching, choose **Select: Matches: All** and specify the hierarchical level at which you want to select the matches. Matched spots are highlighted in green. The matching of any red spots should be examined. However, this can also be done during data analysis.

## 7.8 Match reports

### 7.8.1 Match statistics table

Choose *Reports*: *Analyze Gels*: *Match Statistics Table* to display the number and percent of matches found for each of the images. By default, the numbers are calculated based on *All MatchSets* in the hierarchy. But you can select a

particular match set by clicking the **Choose MatchSet** icon in the toolbar of the




Figure 7-7. Match Statistics Table.

Match Statistics Table (Figure 7-7).

# 7 Matches7.8 Match reports

## 8.1 Introduction

To study the variations in protein expression among a series of gels, the gels should be matched together (be part of the same match hierarchy).

Data analysis can be carried out at two different levels:

- Analyze Gels: Study protein expression changes within a set of gels, without taking populations into consideration. This type of analysis can be carried out on both MatchSet and Classes sheets. The analytical methods used include scatter plots, descriptive statistics, histograms, and factor analysis.
- Analyze Classes: Find significant protein expression changes between classes of gels. For this type of analysis, images must be placed in classes and opened in a Classes sheet. The analytical methods used include descriptive statistics per class, histograms, overlapping measures, and statistical tests.

## 8.2 General settings

The following concepts and settings are used in the different analytical methods described later in this chapter.

### 8.2.1 Ougntification value

The quantification value found in various tables, graphs and plots is a software option. Choose *Tools: Options* and select the *Quantification* tab to set:

- Value: To be used for the analysis of conventional 2-D gels. You can choose between Intensity, Area, Vol, %Vol, or Saliency. The default unit is %Vol.
- **DIGE Value:** To be used for the analysis of DIGE gels. You can choose between Intensity, Area, Vol, %Vol, Vol Ratio, or Slope. The default unit is Vol Ratio.




Central tendency and dispersion are the most frequently-used descriptive statistics. They are calculated in Gel or Class Analysis Tables and Histograms.

## 8.2 General settings

These statistics summarize spot values from a match. The central tendency allows you to localize a center for the data, whereas the dispersion indicates how closely the data points fall around the center.

**Note:** Absent spots, with zero values, are also considered in the calculation of statistics (for both central tendencies and dispersion).

Specify and display the statistics by clicking the *Statistics* icon in a report toolbar (Figure 8-1):



**Figure 8-1.** Statistics common to Gel or Class Analysis Tables and Histograms. The third section is specific to different report types and is described in the corresponding sections of the User Manual

### Central tendency

The central tendency gives the general location of a variable. This is commonly calculated by the arithmetic mean (also known as the average or center of gravity of a distribution), the median (the middle value which divides the sample in two equal parts) or the midrange (middle location between the two sample extremes)

Mean and midrange values are very sensitive to extreme values (outliers) and can be seriously affected by a single observation. On the other hand, the median is highly resistant to outliers. A compromise is given by the trimmed mean (or trimmed midrange) where a predefined number of outliers are removed from the sample. The trimmed measures are more robust than the mean (or midrange) but are more sensitive than the median.

The percentage slider in the Statistics window allows you to remove outliers and obtain the different central tendencies. A 100% value means that all the spot

values available in a match are used to calculate the statistics (no outliers are suppressed). With a value of 80%, for example, 10% of the minimum values and 10% of the maximum values are discarded from the sample and the trimmed measure is calculated with the remaining values.

- Select *Mean* and the arithmetic mean is calculated, that is, the sum of all the sample values divided by the sample size.
- Select *Midrange* and the midpoint of the sample value is calculated, that is, the middle location between the two sample extremes.
- Obtain a Trimmed Mean (or Trimmed Midrange) by selecting Mean (or Midrange) and discarding the desired percentage of outliers with the percentage slider.
- Obtain the **Median** by selecting **Mean** and discarding 50% of the values at each of the extremities, that is, select 0% in the percentage slider.
- Select Reference... and the value from the spot in a specified reference gel is taken as the central tendency. This option is only available for Gel Analysis Reports.

### **Dispersion**

The dispersion measures the variability of the sample data as indicated by how clustered or scattered the data points are around their center value. There are numerous measures of variability: standard deviation, range, interquartile range, and so on.

Like the statistics for central tendency, these measures make use of all the available sample data and can be heavily influenced by outliers. Therefore, you can also restrict the sample to the central values by trimming out the extreme values with the percentage slider.


- Select Mean Squared Deviation (M.S.D.) and the square root of the average squared difference of each sample value to the center location is calculated.
- Select *Mean Absolute Deviation (M.A.D.)* and the mean of the absolute difference between each value and the central value is calculated. It is not affected as much by outliers as the Mean Squared Deviation because the differences are not squared.
- Select *Half-range Size* and the difference between the largest and the smallest values divided by 2 is calculated.

### **Examples**


• The Mean 100% and the Mean Squared Deviation 100% are the most commonly-used statistics (Figure 8-2, a). Note that the standard deviation is the Mean Squared Deviation multiplied by N, where N is the sample size.

This difference comes from the fact that the standard deviation should be an unbiased estimator

- The Median (Mean 0%) and Mean Absolute Deviation 100% are much more robust to outliers than the statistics above (Figure 8-2, b).
- The Midrange 100% and Half-range 100% define an interval that includes all sample values (Figure 8-2, c).
- The Midrange 50% and Half-range 50% are known as order statistics and interquartile ranges (Figure 8-3).



**Figure 8-2.** Histograms showing the sensitivities of central and dispersion values. (a) Mean 100% and M.S.D. 100%, (b) Median and M.A.D. (c) Midrange 100% and Half-range 100%.



**Figure 8-3.** Histograms showing the effect of suppressing outliers. Midrange and half-range values are given for (a) 100%, (b) 80%, (c) 50% and (d) 33%.

**Note:** Histograms are described further in 8.3.3.

## 8.3 Analyze gels

## 8.3.1 Scatter plots

To analyze gel similarities or experimental variations such as disparities in stain intensity or sample loading, you can produce Scatter Plots for matched spots (Figure 8-4) by choosing *Reports: Analyse Gels: Scatter Plots*.




Figure 8-4. Scatter plot for matched spots.

Scatter plots give an idea of the relationship between the spot values from two gels by searching for the linear dependence between the spot values of one gel (variable X) and the corresponding values in the sheet reference (variable Y). Remember that you can change the sheet reference (darker gel name) by choosing *View: Sheet: Set Reference*.

The linear dependence is defined as the best-fit line through the data points. The best-fit line is described by a slope and its offset from the equation  $y = \text{slope} \times x + \text{offset}$ .

The goodness-of-fit for this approximation is given by the correlation coefficient *Corr.* This coefficient can vary between -1 and 1, where an absolute value near 1 indicates a good fit. The spot values of one gel can be predicted, to some

extent, by the values of the other gel. On the other hand, a low value indicates that the data could not be approximated by a straight line.

The types of conclusions that can be drawn from the regression line equations and the correlation coefficients are:

1.0 × × + 0	and Corr = 1	indicates that the spot values for all matched spots are the same in the two gels.		
1.2 × × + 0	and Corr = 0.95	indicates that almost all spot values are approximately 20% higher in the sheet reference.		
1.0 × × + 0.2	and Corr = 0.95	indicates that almost all spot values are shifted by +0.2 with respect to the sheet reference.		

In general, when the data are highly correlated (Corr close to 1) but the best-fit line is far from identity ( $1.0 \times x + 0$ ), you should search for possible reasons to explain why your values are systematically biased. Stain intensity variations, differences in protein loading, or image acquisition problems are some typical causes.

In the Scatter Plots window you can visualize a scatter plot for each gel in the sheet versus the sheet reference, together with the best-fit line, correlation coefficient and the number of matches displayed.

Scatter plots are interactive. You can click on the points representing the matched spots. This in turn selects the spots in the gels and other reports.



Move the slider in the toolbar of the Scatter Plots window to view the scatter plots for the other images in the sheet.

## Scatter table

Click the *Scatter Table* icon in the toolbar to show or hide the table below the plot. The Scatter Table displays, for each pair of matched spots in the scatter plot, the corresponding spot values in the gels and the error in relation to the regression line.

## Copy formula to clipboard

Select the *Copy Formula to Clipboard* option from the *Save* icon in the drop-down list to copy the regression formula and correlation coefficient to the clipboard for use in other applications.

### Gel analysis table 8.3.2

The Gel Analysis Table and Histograms provide valuable tools for looking at unusual matches within a set of gel images. Analyzing protein expression changes, checking spot detection or verifying matching operations are a few of the potential uses. Choose **Reports**: **Analyse Gels**: **Table** to display the Gel Analysis Table (Figure 8-5). This report describes, for each match:

- The spot values from each gel.
- The Max value, that is, the highest of all these spot values.
- The **Central Tendency** and **Dispersion** over all the gels in the sheet, regardless of if they belong to different populations or not.
- The Coefficient of Variation (Coef. Variation), which is the dispersion divided by the central tendency. It measures the relative variability of the spots in a match by correcting for the magnitude of the data values, thus giving a measure that has no units. When you choose the Median and Mean Absolute Deviation statistics, this measure is also known as the Coefficient of Dispersion.
- The Range Ratio, which is the maximum value divided by the smallest value in the sample specified. To specify the sample, click the **Statistics** icon in the Gel Analysis Table toolbar, and suppress outliers by setting the percentage slider in the *Range Ratio* section to the desired value.
- The Separability, which is the highest difference between two consecutively sorted values in the whole sample. It measures the greatest gap that you can have if you want to split the spot values in a match into two separate classes.
- The *Match Count*, which is the number of gels in which the spot is present and matched.




Figure 8-5. Gel Analysis Table.

In addition to the standard functionalities for saving, printing, copying content to the clipboard and navigating in the report, the following tools in the Gel Analysis Report are particularly useful:



## Y Select by value

Select items in the table based on a numerical search criterion. Choose the measure (that is, column) to be used for refinement, and set the lower and/or upper limits of your search interval.



## Create spot set

Create a spot set to annotate spots of interest for use at a later stage. All currently selected spots are automatically included in the newly created spot set, which appears as a column in the table. You see a checked box for spots that belong to the set, or an empty box for spots that do not belong to the set.



## Factor analysis table

Carry out a factor analysis. For more information, see 8.3.4.



Set the statistics to be used in the report. These settings are common to the Gel Analysis Table and Gel Analysis Histograms.



## **i** Settings

Some of the above-mentioned columns may not be displayed by default. Click the **Settings** icon to show or hide columns in the table.



## Normalization

In order to simplify comparisons across matches, the spot values in the Gel Analysis Table can be normalized relative to their gel analysis statistics. Select the desired type of normalization in the drop-down list in the toolbar. The following options are available:


Value	Raw spot value.	
Relative	Spot value – Central tendency	
Ratio	Spot value Central tendency	
Normalized	Spot value – Central tendency    Dispersion	



## Gel analysis histograms

Histograms are a way to look at matched spots. The Gel Analysis Histograms window (Figure 8-6) is displayed by choosing Reports: Analyze Gels: Histograms.

In the histograms, the vertical orange bars correspond to the spot values, the blue horizontal line represents the chosen central tendency and the red lines



delimit the range defined by [Central value - Dispersion, Central value + Dispersion]. The Match ID is displayed at the bottom of each histogram.

Figure 8-6. Gel Analysis Histograms.

The following tools in the Gel Analysis Histograms are very useful:



To see the matches displayed on additional pages, use the slider in the toolbar. When a match is selected on an image or in another report, it is automatically highlighted in a histograms window.




Set the statistics to be used in the histograms. These settings are common to the Gel Analysis Table and Gel Analysis Histograms.



## **Settings**


Click the **Settings** icon to select one of the following options:

Sorted Values, to sort the spot values in ascending order (Figure 8-7).



**Figure 8-7.** Histograms on matches with (a) unsorted values and (b) values sorted in ascending order.

 Adaptive Gradations, to adjust the histogram gradations according to the spot values in each match. Deselect this option to display an identical gradation in all histograms (Figure 8-8).



**Figure 8-8.** Histograms on matches with adaptive gradations set (a) individually for each histogram and (b) according to the minimum and maximum values in all histograms.

• **Show Labels**, to display a table with the gel names. Alternatively, when you place the mouse over a letter in the histograms, a screentip displays the full image name and the spot value in the image.

Value	v	Normalization
-------	---	---------------

You can display normalized spot values to simplify the comparisons across matches. These normalized values are particularly useful in combination with the histograms:

Value	Raw spot value.				
Relative					
	Spot value – Central tendency				
	This normalization sets the central tendency values to 0, and if the Adaptive Gradations option (see above) is deactivated, you have a good overview of the dispersion and therefore of the homogeneity of the matches. This normalization is sensitive to high spot values.				
Ratio					
	Spot value Central tendency				
	Central tendency				
	This normalization divides all values by the central tendency and thus gives a ratio for all data. If you deactivate the Adaptive Gradations option (see above), all histograms have the same scale and thus it becomes easier to detect matches that do not have homogenous values. This normalization is more sensitive to low spot values.				

## 8.3.4 Factor analysis

The visual task of comparing gels is rather difficult when dealing with a large number of gels that consist of thousands of spots. It can be hard to determine if different sample populations exist and to characterize their different protein expression profiles. Factor analysis helps here by condensing the information contained in such huge data sets into a smaller number of factors, or dimensions, that explain most of the variance observed. The factor analysis tool is used to examine the interrelationships between large numbers of variables (that is,, spot values for a series of gels) and to explain these relationships (for example, gel populations) in terms of common underlying factors (or associations with specific spot patterns).

Factor analysis is a complex statistical technique, whose comprehensive description is beyond the scope of this manual. For more information, references are given in the Appendix.

### Perform a factor analysis

A factor analysis is carried out on all or selected matches in the Gel Analysis Table. You must judge which of the options described below is most applicable to your analysis.

## To carry out a factor analysis:

- 1 Click the *Factor Analysis Table* icon in the Gel Analysis Table toolbar.
- 2 If any matches are selected, the software asks if you want to use all or only the selected rows
- The Factor Analysis Table (Figure 8-9) is displayed with the lines corresponding to the axes that can be drawn in the Factor Projection Plot. Select two axes to be displayed in this plot (generally the first two ones).


4 Click the *Factor Projection Report* icon in the Factor Analysis Table toolbar.

Fa	Factor Analysis Table [Mean(100%) and M.S.D]   ■ □ ▼ ×								
6	<b>□</b> ▼ <b> </b>								
	Axe	%Variance	A_T1_Gel1	A_T1_Gel2	A_T1_Gel3	A_T2_Gel1	A_T2_Gel2	A_T2_Gel3	
1	1	91.8391	0.329719	0.427523	0.441135	0.400617	0.417108	0.423599	
2	2	5.41506	-0.282403	-0.435063	-0.495152	0.432026	0.357638	0.413814	
3	3	1.01623	-0.3449	0.269066	-0.044982	0.178557	0.561399	-0.677918	
4	4	0.821133	-0.669511	0.589342	-0.0975621	-0.104748	-0.235228	0.358619	
5	5	0.59571	0.46165	0.444102	-0.688504	0.243621	-0.205315	-0.11878	
6	6	0.31278	-0.176669	-0.104868	0.27321	0.742061	-0.534285	-0.21687	

Figure 8-9. Factor Analysis Table.

The Factor Projection Plot (Figure 8-10) displays the projection of each match (cross) and each gel (blue vector) on the two factorial axes. The blue curve represents a part of the correlation circle; its form is linked to the scale of the axes. When matches are selected on the plot, they are automatically highlighted on the gels and/or any open reports.

The Factor Projection Table (Figure 8-10) displays the contribution of each match to the two axes displayed in the Factor Projection Plot. The *Quality* measures if the projection of the match is well represented on the factorial subspace.



**Figure 8-10.** Factor Projection Report includes the Factor Projection Plot (top) and the Factor Projection Table (bottom).

You can use the following tools in the Factor Projection Report:



Display the Factor Analysis Table from which the two axes for the current Factor Projection Report were selected.

## Factor projection table

Show or hide the Factor Projection Table corresponding to the displayed plot.

## Factor projection plot

Show or hide the Factor Projection Plot corresponding to the displayed table.

Choose the items to be displayed from the drop-down list in the toolbar of the plot. You can choose matches (crosses), gels (blue vectors), or both.

## Matches displayed

Enter the number of matches to be displayed in the plot.

### Interpret a factor analysis

How to interpret a factor analysis is explained with an example of six gels run to compare the effect of two treatments T1 (A\_T1\_Gel1, A\_T1\_Gel2, A\_T1\_Gel3) and T2 (A\_T2\_Gel1, A\_T2\_Gel2, A\_T2\_Gel3) on bacteria cultivated on substrate A. The gels were detected and matched, and a Gel Analysis Table displayed.

The Factor Analysis Table (Figure 8-9) summarizes the variance accounted for by successive axes (or factors), expressed as a percent of the total variance. Thus, factor 1 accounts for 91.8% of the variance, factor 2 for 5.4%, and so on. The coordinates for each gel projected on these axes is also listed.

The number of factors equals the number of gels being analyzed. Factor analysis cannot be performed with less than two gels. The more gels you use, the more reliable the results are likely to be, and the more factors are calculated. Since the first factors are generally the best ones for characterizing gels and matches that behave similarly, the factors are ranked in order of importance.

Figure 8-11 shows the Factor Projection Plot obtained when the first two axes in Figure 8-9 were selected. In the example, only the 20 most significant matches are displayed on the projection plot. The further away a spot is from the origin, the more important it is likely to be in terms of characterizing the gels. If all matches were shown, one would find that many of them cluster around the origin of the graph. This illustrates that the majority of matches, that is, protein spots, are not meaningful in classifying the gels.

**Figure 8-11.** Factor Projection Report. The ten matches with highest contribution to the second axis are selected.

To find a possible meaning for a given factor, one should first identify the matches that largely contribute to this factor. The Factor Projection Table is used for this purpose. When sorting the matches according to their contribution to the first axis, one discovers that the matches at the top of the table are those with the highest relative volumes, as found in the Gel Analysis Table. In fact, the first axis is generally correlated with protein abundance.

The Factor Projection Table also contains the Quality for each match. This number tells you how close the distance of the projection is to reality. Matches with very similar behavior (similar expression profiles across gels) are close in space. However, when projected onto a two-dimensional subspace, matches that are actually far apart may appear close together. It is therefore important to look at the Quality to judge if matches are effectively close. If the values are high for both matches, the chance is great that they are indeed nearby and have

a similar behavior. If one of the matches has a low value, any interpretation becomes tentative.


Gels that are adjacent on the graph are likely to be similar to each other. They may correspond to the same population. This is clearly the case in the example above. The gels from T1 cluster together below the horizontal axis, whereas the gels corresponding to T2 lie above.

The closer a match is to a given set of gels, the more characteristic it is likely to be of those gels. That is, the more important the match is in determining why those gels are different from other gels. In Figure 8-11, one can observe that the matches 48, 123, 180, 181 and 248 are close to the gels belonging to T1. The histograms in Figure 8-12 show that these matches have high spot values for the T1 gels and low values for the T2 gels. Matches 237 and 338, on the other hand, are characteristic of the T2 population, with higher spot values in the T2 gels. In our example, the second axis appears very important for separating the gels into two classes. It is related to the ratio between the mean spot values in each population of gels. The matches in the upper part of the graph have ratios that favor the T2 population, while those below the horizontal axis have ratios that favor the T1 gels.

### Comments on factor analysis

Factor analysis is used to examine the protein expression pattern within each match. The quality of the output depends on the quality of spot matching. Hence, it may be useful to exclude spots that are not well matched across all gels (using *Match Count* in the Gel Analysis Table). Nevertheless, in cases where a majority of spots were properly matched, the inclusion of all matches in the factor analysis can yield good results with no preliminary match filtering necessary.

This statistical method, based on data variation and their standard deviations, highlights the natural formation of populations among the gels and allows identification of matches (that is, matched spots) that are characteristic of these classes. However, one should be very critical when analyzing factor analysis plots since the results can be greatly influenced by outliers, bad matches, and so on. Factor analysis can provide valuable indications in some cases, but not in others.



**Figure 8-12.** Gel Analysis Histograms for the ten matches ranked highest for their contribution to the second axis in the factor analysis example. Matches 48, 123, 180, 181 and 248 are characteristic of population T1, matches 237 and 338 are characteristic of population T2.

## 8.3.5 DIGE histogram

A DIGE Histogram can be displayed when detected DIGE gels are selected. It shows, for each non-reference image, data associated with detected spots in the image, plotted against Log Volume Ratio on the X-axis (Figure 8-13). It has two different Y-axes:

- The left Y-axis displays the spot frequency. The blue curve represents the frequency distribution of the log volume ratios.
- The right Y-axis represents the Measure parameter (see below) selected from the *Measure* tool in the DIGE Histogram window. A plotted single data point on the histogram represents an individual protein spot.

The name of the DIGE reference for the current image, used for the calculation of the Log Volume Ratio, appears at the lower right corner of the DIGE Histogram window.

The DIGE Histogram is dynamic. You can click on the data points representing spots to select them on the gels and any other open reports.




Figure 8-13. DIGE Histogram.

Two report-specific tools are available:



Move the slider in the toolbar of the DIGE Histogram window to view the histograms for the other images in the sheet.



Select one of the following options from the *Measure* icon:

- *Max Slope*: Largest gradient associated with the co-detected spots.
- Area: Number of pixels within the spot boundary.
- Max Intensity: Largest pixel value associated with the co-detected spots.
- Max Volume: Volume of the largest co-detected spot.

### 8.4 Analyze classes

### 8.4.1 **Specify classes**

To identify protein expression variations between populations of gels, you need to specify what gels belong to which population by creating classes. Classes are created in the Workspace.

Often, you already know the populations for your set of gels. This is the case when you are comparing gels from healthy tissue extracts with those from disease-associated samples. When you have no preliminary knowledge of the populations in the set of gels, you can draw possible conclusions from factor analysis.


### 8.4.2 Overlapping measures

Spot values for a given match within each class can be summarized by the central tendency and dispersion. In addition to these descriptors, the software computes overlapping measures between the class intervals, where class intervals are defined by the ranges [Central value - Dispersion, Central value + Dispersion].

The overlapping measures quantify the overlap between these intervals, and thus evaluate how different the protein expression changes between the classes really are. They take into account both the difference between the central tendency in each population and the dispersion. The following overlapping measures are available:

- **Gap**: Maximum difference between the range of the current class and the range of one of the other classes (Figure 8-14, c-b in the case of Class A). Negative values indicate overlap whereas positive values are nonoverlapping class ranges.
- Ratio: Maximum ratio between the lower limit of one of the other classes and the upper limit of the current class (Figure 8-14, c/b in the case of Class A). Absolute values smaller than 1 indicate overlap, whereas absolute values higher than 1 show that there is no overlap. In order to easily distinguish matched spots that are over or under expressed in one of the classes, the ratio value is preceded by a minus sign when the protein spot is under expressed for the class in question, compared to the other class. Positive values are attributed to the Ratio value in the overexpressed class.
- Normalized: Maximum percentage of the current class range not overlapping with the range of one of the other classes (Figure 8-14, (c-a)/ (b-a) in the case of Class A). A value smaller than 1 indicates overlap. For example, 0.25 implies that 25% of the current class range is not recovered by one of the other classes. In the same way, a value of 1.5 indicates that

there is a gap equivalent to 50% of the current class range to the furthest other class. The normalized overlapping is not symmetrical, the value from Class A compared to Class B is not the same as the value from Class B compared to Class A.



**Figure 8-14.** Scheme demonstrating how the Gap, Ratio, and Normalized values are calculated for two non-overlapping classes (upper part) and two overlapping classes (lower part). Arrows above or below each class range illustrate how the Normalized measure relates to this class range.

**Note:** The formulas mentioned above only apply to Class A in the current example. Their purpose here is to illustrate the principles of overlapping measures. Many different cases (and therefore formulas) exist.

A value of 1e6 (1000000) in the Class Analysis Histograms characterizes the case where the protein is completely absent from a class. The software cannot compute ranges. A value of 0 for the Ratio or Normalized measures indicates that the particular class is entirely covered by another one.

**Note:** As indicated in their definition, the Gap, Ratio, and Normalized values always calculate the MAXIMUM difference, ratio or percentage with respect to ANY of the other classes. This means that when the Ratio values for three populations (for example, A, B, C) are compared, the software calculates the ratio of A with respect to B and of A with respect

to C, but only displays the highest value in the column for class A. The number shown does not indicate with respect to which class the value was obtained. The idea is to quickly enable you to find a match (that is, protein marker) that distinguishes the current class from any of the other classes. Once such protein markers are found, the Class Analysis Histograms can be used to study the match in more detail.

### 8.4.3 Statistical tests

The software provides three statistical tests: ANOVA, Mann-Whitney U test, and the Kolmogorov-Smirnov test. These tests are used to analyze differences in protein expression between classes of gels. The idea is to draw conclusions about the significance of the protein expression changes by extrapolating information from the data you collected. For example, when you have two samples (classes) with different means (that is, different means for the spot values of a particular match), you might want to know if the data were sampled from populations with different means or if the populations have the same mean with the observed difference being a coincidence of random sampling.

You can calculate the probability of observing a certain difference (or larger) between sample means in an experiment of this size, for populations that in reality have the same mean. If the probability is small, you can conclude that the difference is not likely to be caused by random sampling and assume instead that the populations have different means.

Note: You can display the desired statistical values for each match in the Class Analysis Table. These values should be considered as qualitative indications of the variations in protein expression between two populations. To draw quantitative conclusions, you must verify that the restrictive assumptions of the various tests are met. In addition, one should always check the results by visual inspection of the spots, since the conclusions may be erroneous due to inaccuracies in detection or matching.

**Note:** The given statistical values are useless when the samples (classes) do not consist of more than two gels. Your objective should always be to work with the largest possible sample sizes.

### One-way ANOVA

Analysis of Variance (ANOVA) is one of the most important statistical tests available for biologists. It is essentially an extension of the logic of Student's t-tests to those situations where the comparison of the means of several groups is required. Thus, when comparing two means, ANOVA gives the same results as the t-test for independent samples.

One-way ANOVA tests the null hypothesis that all populations have identical means. It generates a P value that answers this question: If the null hypothesis

is true, what is the probability that randomly selected samples vary as much (or more) than actually occurred?

It is based on the same assumptions as the t-test:

- The samples are randomly selected from, or at least representative of, the larger populations.
- The two samples are independent. There is no relationship between the individuals in one sample as compared to the other.
- The data are sampled from populations that approximate a Gaussian distribution.

If you are not able to assume that your data are sampled from Gaussian populations, then non-parametric tests like the Mann-Whitney or Kolmogorov-Smirnov tests can provide a better analysis. However, these test only allow you to compare two samples at the same time.

### Mann-Whitney or Wilcoxon test

The Mann-Whitney U test or rank sum test is the non-parametric substitute for the two-sample t-test when the assumption of normality (Gaussian bell-shaped distribution) is not valid. It is equivalent to the Wilcoxon rank sum test. It should only be used for comparing two unpaired samples. The assumptions of the Mann-Whitney U test are:

- The variable of interest is continuous (not discrete) and the measurement scale is at least ordinal. This means that repeated values (ties) are not acceptable. When ties are present in your data, there is an approximation provided in the calculations, but the exact results no longer hold.
- The distributions of the two samples are identical (although not necessarily normal) and differ only in location (that is, central tendency).
- The two samples are independent.

To perform the Mann-Whitney test, the software first ranks all the spot values from low to high, paying no attention to which of the two classes (for example, X and Y) each value belongs. Then each value is given a rank number. The smallest number gets a rank of 1. The largest number gets a rank of N, where N is the total number of spot values in the two classes. If two values are the same, then they both get the average of the two ranks for which they tie. Finally, the ranks in each class are summed, thus giving WX and WY, which are used to calculate the Mann-Whitney test statistic, U. The formula for UX is as follows (the formula for UY is obtained by replacing X by Y):

$$U_{X} = W_{X} - \frac{n_{X}(n_{X} + 1)}{2}$$

The smaller of the two calculated U values corresponds to the number of shifts needed in order that the spot values from the two populations do not overlap. For the first example in Figure 8-15, no shifts are necessary since the spot values of classes X and Y are already separated. On the contrary, for the second example the Mann-Whitney test indicates that the first two values from Y (or the three last values from X) have to be swapped four times in order to separate the samples.

The software displays the smaller of the two calculated U values in the Class Analysis Table. The lower the number, the higher the probability is that the means of the two samples are different. Knowing this value, and the sample size, you can easily look up the probabilities in a Mann-Whitney table.

Attention should be payed when analyzing the results of a Mann-Whitney test. First, the assumptions are often violated. This is the case when spots are completely absent in one of the classes (in this case you have repeated values of 0). Moreover, if you have small samples, the Mann-Whitney test is meaningless. If the total sample size is seven or less, the test always gives a probability (of finding different means, in the case of identical populations) greater than 0.05, no matter how little the samples differ.




Figure 8-15. Two examples to illustrate the Mann-Whitney and Kolmogorov-Smirnov

tests. In (a) the two classes do not overlap at all, whereas in (b) 4 shifts are needed to completely separate the spot values of the two classes. The bold values correspond to the spot values and ranks from class Y, the others belong to class X.

## Kolmogorov-Smirnov test

The Kolmogorov-Smirnov test tries to determine if two data sets differ significantly. It is used to test if two samples may reasonably be assumed to come from the same distribution. It has the advantage of making no assumption about the distribution of the data and is frequently preferred over the Mann-Whitney rank sum test where there are many ties (repeated values). Other tests (for example, the t-test) may be more sensitive if the data meet the requirements needed for that test. The assumptions of the Kolmogorov-Smirnov test are:


- The probability distributions are continuous.
- The measurement scale is at least ordinal.
- The two samples are mutually independent.

In the Kolmogorov-Smirnov test, the data points in each sample (spot values for a particular match in a class) are sorted in ascending order and converted into an empirical distribution function (EDF). This function gives the fraction of data points to the left of a given value z. For the second example in Figure 8-15, the ordered data points from class X are: 0.034, 0.045, 0.056, 0.064, 0.069, and 0.078. The fraction of data points to the left of each of these z values can easily be calculated and plotted (full line) in an Empirical Distribution Plot (Figure 8-16):

It is clear that no data lie strictly below 0.034, 1/6 = 17% of the data is strictly smaller than 0.045, 2/6 = 33% of the data is strictly smaller than 0.056, 3/6 = 50% of the data is strictly smaller than 0.064, and so on.

The same procedure can be followed for the second sample (class Y in our example, the dashed line in Figure 8-16). The Kolmogorov-Smirnov test statistic D is then defined as the maximum distance between the empirical distribution functions (EDF) for the two samples. In the example, D is 0.63 (0.83-0.20). If D is greater than a particular decision limit (critical value found in a Kolmogorov-

Smirnov table), there is a statistically significant difference between the two samples. However, the test provides no insight as to what causes the difference.



**Figure 8-16.** The Empirical Distribution Plot for the spot values of match 613 (Figure 8-15), in Classes X and Y. The Kolmogorov-Smirnov statistic D corresponds to the maximum distance between the two empirical distribution functions.

### 8.4.4 Class analysis table

The Class Analysis Table provides valuable data for finding significant protein expression changes between populations of gels. Using this data, it is possible to differentiate one class from the others based on a few matches/spots. Choose *Reports: Analyse Classes: Table* to display the Class Analysis Table (Figure 8-17). For each match, there is a description of:

- The Center (central tendency), Dispersion, Gap, Ratio or Normalized values for each class, depending on the selection made in the drop-down list (see below). By default, the Center values are displayed.
- The *Max* value, that is, the highest from all these class values.
- The *Match Count*, which is the number of classes in which the spot is present and matched.
- The results from the statistical tests: the ANOVA probability P, the Wilcoxon U statistic, and the Kolmogorov D statistic.

The columns corresponding to spot sets.




Figure 8-17. Class Analysis Table.

In addition to the standard functionalities for saving, printing, copying content to the clipboard and navigating in the report, the following tools in the Class Analysis Table are particularly useful:



## Y Select by value

Select items in the table based on a numerical search criterion. Choose the measure (that is, column) to be used for refinement, and set the lower and/or upper limits of your search interval.



## Create spot set

Create a spot set to annotate spots of interest for use at a later stage. All currently selected spots are automatically included in the newly created spot set, which appears as a column in the table. You see a checked box for spots that belong to the set, or an empty box for spots that do not belong to the set.



## **Statistics**

Set the statistics to be used for calculating the Center and Dispersion value of each class, and consequently the Gap, Ratio and Normalized values. These settings are common to the Class Analysis Table and Class Analysis Histograms.

**Note:** The Center and Dispersion values define the interval that characterizes the protein sample of each class in a match. To characterize a class only by the central value, set the Dispersion percentage slider to 0%. This is useful if you want to calculate the difference or ratio between the central values of your classes.



Some of the above-mentioned columns may not be displayed by default. Click the **Settings** icon to show or hide columns from the table.



The displayed statistical descriptor or overlapping measure can be selected from the drop-down list in the toolbar. You can display the Center (central tendency), Dispersion, Gap, Ratio or Normalized values for the different classes.

### 8.4.5 Class analysis histograms

You can visually investigate the statistical and overlapping descriptors of populations by displaying Class Analysis Histograms. Choose Reports: Analyse Classes: Histograms.

When the Center (central tendency) value is selected for display in the dropdown list, the Class Analysis Histograms display all the individual spot values in each match, separated for each class by vertical gray lines (Figure 8-18). The classes are characterized by their central tendency (blue horizontal line) and dispersion interval (bounded by the outer red lines). The Match IDs appear below each histogram.

When displaying the Dispersion, Gap, Ratio or Normalized values (Figure 8-19), each class is represented by a single value (red bar).




Figure 8-18. Class Analysis Histograms with Center values displayed.




Figure 8-19. Class Analysis Histograms with Ratio displayed.

The histograms can be selected to highlight the corresponding spots on the gels and in any open reports.

The following tools in the Class Analysis Histograms are very useful:




To see the matches displayed on additional pages, use the slider in the toolbar. When a match is selected on an image or in another report, it is automatically highlighted in a histograms window.

## **Statistics**

Set the statistics to be used for calculating the Center and Dispersion value of each class, and consequently the Gap, Ratio and Normalized values. These settings are common to the Class Analysis Table and Class Analysis Histograms.


## **Settings**

Click the **Settings** icon to activate one of the following options:



**Figure 8-20.** Class Analysis Histograms with (a) unsorted values and (b) values sorted in ascending order.

 Adaptive Gradations, to adjust the histogram gradations according to the values in each histogram. Deselect this option to display an identical gradation in all histograms (Figure 8-21).



**Figure 8-21.** Class Analysis Histograms with adaptive gradations set (a) individually for each histogram and (b) according to the minimum and maximum values in all histograms.

• **Show Labels**, to display a table with the gel or class names. Alternatively, when you place the mouse over a letter in the histograms, a screentip displays the full image or class name and the value in that image or class.



The displayed statistical descriptor or overlapping measure needs to be selected in the drop-down list in the toolbar. You can display the Center (central tendency), Dispersion, Gap, Ratio or Normalized values for the different classes.

8 Data analysis8.4 Analyze classes

# **Annotations**

#### Introduction 9.1

Individual pixels and spots in a gel image can be labeled with annotations. Annotations are used to flag items with their specific characteristics (protein name, database accession number and so on) or to mark spots with common characteristics. They also offer the possibility of linking and associating gel objects to external guery engines or data sources of any format (text, html, spreadsheet, multimedia, database) located locally or on the Internet.




Figure 9-1. Annotations are composed of an annotation basis (square or cross), an annotation flagpole and a set of labels (flag).

An annotation is defined by its position on the gel (X and Y coordinates) and its set of labels. Each label belongs to a predefined or user-defined category. As shown in Figure 9-1, each annotation is composed of a basis, a flagpole, and a flag that consists of a set of colored labels.

### Spot and pixel annotations

You can create two types of annotations (Figure 9-1):

- **Annotations linked to pixels** are visualized with a **cross** basis. They are simply connected to a pixel and have the same coordinates as that pixel.
- Annotations linked to spots are visualized with a square basis. They are linked to a spot and have the same coordinates as that spot, that is, as the spot's center of gravity. These annotations are automatically selected when the linked spot is selected and vice versa.

## 9.2 Create annotations and labels

#### 9.2.1 Create annotations

The creation of annotations essentially consists of taking three steps.

- 1 Select the spots or pixel on which you would like to create an annotation:
  - Choose the Select tool and double-click on a spot or a pixel.
  - Select several spots, make sure to select only the images on which
    you want to create new annotations, and choose Edit: Annotations:
     Add.
- 2 Choose one of the existing categories or create a new category. The creation of new categories is described in Section 9.3.
- 3 Enter the desired label content.

### 9.2.2 Add labels to existing annotations

Labels can be added to an existing annotation.

- 1 Select the annotation to which you would like to add a label:
  - Choose the **Select** tool and double-click at the basis of the annotation to which your label should be added.
  - Select several annotations, make sure to select only the images on which you want to create new labels, and choose Edit: Annotations: Labels: Add.
- 2 Choose one of the existing categories or create a new category. The creation of new categories is described in Section 9.3.
- 3 Enter the desired label content.

**Note:** One annotation may have many labels, but it can only contain one label from each category. If one spot contains several proteins, it may need to carry different labels from the same category. You can do this by linking additional annotations to the spot.

### 9.2.3 Link annotations to spots

The software allows you to link an annotation to a spot. This is helpful when you want to link an additional annotation to a spot (to attach multiple labels of the same category), or when you missed a spot to which you wanted to link a newly created annotation.

To link an annotation with a spot, click on the annotation basis and drag it to a spot. If, for some reason, an annotation already exists within a spot but is not yet

linked, you can also select the annotation and choose Edit: Annotations: Link with Spot.

To unlink an annotation from a spot, select the annotation and drag it outside the spot or choose Edit: Annotations: Unlink from Spot.

#### 9.3 Create label categories

#### **Predefined label categories** 9.3.1

The software offers some predefined label categories:

- **Ac**: This category is provided to hold the protein's accession number (AC) taken from a user-defined database and can be the linked to ImageMaster's remote database guery engine. When such a link is defined, double-clicking on a label of this category displays the corresponding protein entry in the selected database with the default Internet browser.
- pl MW: This category contains the known isoelectric point (pl) and molecular weight (MW) information, which is subsequently used to compute approximate pl and MW values for any point in a gel. You should enter the pI value first and MW value second, separated by a space. By replacing one of the values by -1, you indicate to the program that no value is set.
- **Comment**: This category is an example of a general category where users may store their comments.

#### 9.3.2 Create new categories

When you create a new category during one of the procedures described in Section 9.2, you must enter the constraints and attributes for the new category in the Create Category window (Figure 9-2). This section explains the different options.

**Note:** User-defined categories are only available from the category list as long as there is at least one label of this category in the open gels (in any of the open sheets). To create categories that are permanently available, you must define them in the **Annotations** tab of the **Options** window (available by choosing **Tools : Options**). The same category constraints must be defined.

At any moment, you can change the category constraints by choosing **Edit**: Annotations: Categories: Edit Attibutes, or rename a category by choosing Edit: Annotations: Categories: Rename a Category.




Figure 9-2. Create Category window.

## Data type

By default, the labels can contain any character. However, to ensure coherent annotation data, the label contents can be constrained to one of the following *Data Types*:

- Text. Can contain any character.
- Number: Can only contain numerical values.
- Auto-Numbering. An incremental number (per gel) is entered automatically as the new label.
- Set. Use this data type to mark spots with common properties. The labels
  in such a category do not contain specific information. They only display
  the name of the category to which they belong. Note that labels of this
  category are displayed in the form of check boxes in tables. A checked box
  indicates that the spot belongs to the category.

#### Is unique

When you check the *Is Unique* box, you indicate that each label on a gel within the new category should be unique. The software will not accept a new label when an identical one already exists. You are asked to enter a new label.

#### External engine

ImageMaster offers the possibility to link spots on gel images to protein data in 2-DE or other databases. All you have to do is input the appropriate query format (database address and query engine) in the *External Engine* field of the new label category and enter valid database accession numbers (AC) as labels. This functionality is further described in Section 9.4.

When you subsequently double-click on a label of such a category, the software opens your default Internet browser and launches an HTTP guery that takes the form of a Web page address.

**Note:** You can define a different query engine for each label category and therefore you can link one protein spot to different database entries.

#### Display properties

User-defined categories use a gray background color by default. You can change the default color by clicking on the Color box. The new background color is used for all labels of that specific category.

You can choose to hide the category by default, by clicking the **By default** hidden box. In this case, you must choose View: Annotations: Visible Categories to make the category visible by checking its box.

#### 9.4 Connect to protein databases

#### 9.4.1 Introduction

ImageMaster can link spots on gel images to protein data in 2-DE or other databases. Such databases contain information on proteins identified on 2-DE images, such as pl and MW values, bibliographical references to protein related literature, information on protein functions, etc. If your computer has Internet access, you can remotely query and retrieve protein data related to spots on your gels.

**Note:** The ImageMaster software provides access to several databases on the Internet. It is the responsibility of the user to acquire the database licenses, if needed. In particular, the PROSITE and SWISS-2DPAGE databases are copyrighted, and all commercial users of these databases are required to purchase a database license from GeneBio. No license fee is charged to academic users for non-commercial use. For questions about obtaining a license subscription for the PROSITE and SWISS-2DPAGE databases, please contact GeneBio (www.genebio.com).

ImageMaster takes advantage of the fact that virtually all databases on the Web, and in particular those containing 2-DE and other protein data, use CGI scripts to enable data queries. A CGI (Common Gateway Interface) script is a program or script file executed on a Web server in response to a user request. The CGI script transmits information (such as a database accession number or object identifier) from the client to a database engine, receives back the results, and displays them to the client.

#### 9.4.2 Set the database

To query a remote database through the Internet, you must send HTTP queries (for instance, http://www.expasy.org/swiss-2dpage/P02649) to a CGI script on a server. The HTTP query must be composed of:

- The database HTTP address (http://www.expasy.org/).
- The database query engine (swiss-2dpage/).
- The database accession number (P02649).

In ImageMaster, the database HTTP address and query engine are entered as constraints to the annotation category (see Figure 9-3). Type them in as one string in the *External Engine* field (for example, http://www.expasy.org/swiss-2dpage/).




Figure 9-3. Setting the External Engine as a category constraint.

#### Federated 2-D PAGE databases

A list of federated 2-D PAGE databases, with the required database query formats, can be found at http://www.expasy.org/ch2d/2d-access.html, or by clicking on the *Url* button to the right of the External Engine field. Copy the desired database address and query engine from this site.

#### Other databases

It is possible to find the required query format for databases that are not federated or that do not contain 2-DE data. Directly query the database until you find a specific protein (or other) entry. Then copy the address of the corresponding Web page in your browser to the External Engine field of the Create Category window, without including the accession or identification

number of the current entry. Generally, this address consists of the database HTTP address and guery engine followed by a question mark (?) or equal sign (=).


For example, you might display the entry for the protein structure 1BMT in the Protein Data Bank (PDB): http://www.rcsb.org/pdb/explore/explore.do?structureId=1BMT. The required query format for ImageMaster is: http://www.rcsb.org/pdb/explore/explore.do?structureId=

#### 9.4.3 Query the database

The database accession number (for example P02649) is entered as a label of the particular category linked to a spot.

When you subsequently double-click on the label while the **Selection** tool is selected, ImageMaster opens your default Internet browser and launches an HTTP query that takes the form of a Web page address.

As a result, the entry for the protein with the given accession number opens in your browser (Figure 9-4). In the case of the above example, this would be the entry for Human Apo E (Gels) in the SWISS-2DPAGE database.



**Figure 9-4.** The SWISS-2DPAGE entry for Human Apo E. Entries from this database contain full protein names, bibliographic references, annotations (such as protein function or pathological variations), and the pl and MW of the related spots on the 2-DE maps. It also includes cross-references to numerous other databases.

#### 9.4.4 Connext to an executable

Analogous to an HTTP query, where the content of a label is transmitted to a CGI script on a server, ImageMaster allows you to pass on the content of a label as the first parameter to any executable. When you double-click on a label that has an executable defined in the *External Engine* field of the label category, the executable runs with the label content as a parameter.

To define an executable in the External Engine field, click on the *Exe* button and locate the .exe file.

## 9.5 Create specific links

As seen above, it is possible to link labels to remote database entries by defining an external search engine for a particular category. By using specific keywords in the labels of any category, you can also create links to Web pages, files, and text (Figure 9-5).

To create a specific link, you should add an annotation to a gel and include the following items in the label field:

- A short descriptor that will be the visible part of the label.
- A keyword indicating the type of link (http:, file:, text:).
- A **link content**, which contains the information necessary to establish the link or the content of the link (in the case of a text link).

To indicate that a label is linked, its short descriptor is followed by three dots. When you double-click on such a linked label, you do not enter the typical editing mode, but the link (Web page, file or text) is automatically opened with the appropriate software. Alternatively, to open any linked label, choose *View*:

\*\*Annotations: Linked Data\* in the menu.

**Note:** You can define links in any label category but you can only have one link per label.

### 9.5.1 Http link:

You can link spots or pixels to specific Web pages. A double-click on an http-linked label will launch your Internet browser and automatically call the corresponding Web page.

You can, for example, create a direct link to the ExPASy Proteomics Server (Figure 9-5). In this case, the label content should contain the string "http:" followed by the address of the Web page.

#### 9.5.2 File link:

You can link spots or pixels to software files. Double-clicking on a file-linked label launches the specified file with the default system application associated to the file extension.

The linked files can be placed in a specific directory, which is defined by choosing *Tools: Options* in the menu and by setting the *Annotations* folder in the *Annotations* tab. In this case, you only need to give the name and file extension when creating the link. You can create subfolders in the Annotation folder to arrange your files. The file names indicated in the labels must then contain the name of the subfolder (for example, AA composition\AA P10413.xls).

Alternatively, you can link labels with files located anywhere on your system. You should then include the complete file path when referencing the file.


For example, you can link a protein spot to an Excel file containing the amino acid composition of a protein (Figure 9-5). The label should then consist of the string "file:" followed by the file name (with its extension).

#### 9.5.3 Text link:

In some cases, you might want to associate a long text comment with a specific protein spot, but without overloading the display. The solution is to create a text link, rather than a very long annotation label. Double-clicking on the linked label is sufficient to display a window containing the entire text (Figure 9-5).

Text links are particularly useful for attaching bibliographic references to a spot, for instance, or any other comment such as the one in Figure 9-5. Please note that the string "text:" must first be inserted, followed by the text you would like to associate with the spot.

To connect general information about the gel, other tools are better adapted. Comments can be attached to projects, match sets, and classes in the workspace. Additionally, you can specify gel descriptions.




**Figure 9-5.** The annotation on the spot containing protein P10413 includes linked labels such as a link to a Web site (top), a link to a file (bottom right) and a text link (bottom left).

## 9.6 Select annotations and labels

# 9.6.1 Select

Annotations and labels can be selected with the *Select* tool. The selected labels or annotations are highlighted in green and displayed in front of the other annotations (Figure 9-6).



**Figure 9-6.** (a) Annotation hidden by some other annotations. (b) Same annotation displayed in front of all other annotations, after selection.

If the annotation is attached to a spot, the spot is also selected. Similarly, if you select a spot, the attached annotation is selected.

To select a label, click on the label. You can select more than one label by using the Shift key.

To select an annotation, click on the annotation basis. To select more than one annotation, select the first one, then hold down the Shift key and select additional annotations.

To select all annotations in a region, position the cursor at the top left position of the desired region, hold down the left mouse button, and then drag the cursor to the bottom right position. All annotations in the defined region are selected. To select annotations in more than one region, hold down the Shift key while selecting additional regions.

To deselect all annotations, click on a gel (not on an annotation).

#### 9.6.2 Select menu

You can select specific labels and annotations with the options in the **Select: Annotations** menu. Please note that this also selects hidden annotations.

- **By Content**: This feature enables the selection of labels (belonging to one or several categories that must be selected) based on their content. When the **Regular Expression** box in the window is not checked, the entered string of characters is taken literally, and the program selects all labels containing this string. By activating the Regular Expression option, regular expressions can be used in the search field (see below for details).
- **By Category**: This feature enables the selection of all labels belonging to one or several categories. Use the Shift and Ctrl keys to pick several category names at a time.
- All: This highlights all annotations in the selected gels.

• Common Labels: This option allows the retrieval of all sets of identical labels within a gel or among a series of gels, for all the categories chosen in the window.

## **Regular expressions**

Regular expressions provide a mechanism to select specific strings from a set of character strings. Regular expressions use symbols and syntax elements to describe a generalized pattern. ImageMaster invokes the standard Extended Regular Expressions to search patterns in labels, the essentials of which are summarized in Table 9-1.

Syntax	Description	Example
	Matches any one character.	<b>e.oli</b> matches:eaoli, eboli, ecoli
[]	Matches any character listed between the brackets. [a-z] indicates the range of characters between A and Z and [0-9] is any numeral from 0 to 9.	<b>P[a-d]</b> matches:Pa, Pb, Pc, Pd
[^]	Matches any character except those listed between the brackets.	<b>P[^bd]</b> matches:Pa, Pc, Pe but not Pb or Pd
?	Matches the preceding element zero or one times.	P0?1 matches:P1 and P01.
+	Matches the preceding element one or more times.	<b>P0+1</b> matches:P01, P001, P0001,
*	Matches the preceding element zero or more times.	<b>P0*1</b> matches:P1, P01, P001, P0001,
{n}	Matches the preceding element exactly n times.	<b>P0{3}1</b> matches:P0001 but not P01 or P001
{n,}	Matches the preceding element at least n times.	<b>P0{2,}1</b> matches:P001, P0001, but not P01
{n,m}	Matches the preceding element at least n times, but not more than m times.	<b>P0{1,3}1</b> matches:P01 P001, and P0001, but not P1 or P00001
()	The characters between parentheses form a subexpression.	<b>P(24)+</b> matches:P24, P2424, P242424,
	Matches the expression before or after the vertical line. Mostly used within a subexpression.	P(ab cd)1 matches:Pab1 and Pcd1

Syntax	Description	Example
۸	A circumflex outside a bracket expression anchors the element it starts with to the beginning of a string; such an element can match only a sequence starting at the first character of a string.	^(ec).* matches:ecoli, ecoli_eftu but not eftu_ecoli
\$	A dollar sign outside a bracket expression anchors the element it terminates with to the end of a string; such an element can match only a sequence ending at the last character of a string.	.*(ecoli)\$ matches:ecoli, eftu_ecoli but not ecoli_eftu

Table 9-1. Regular expressions available to search patterns in labels. Please note that the term element used in the description indicates a character or a subexpression.

The characters  $^{[\$]}$ +?{\ have a special meaning in certain contexts. If your labels contain any of these special characters, you must enter a backslash in front of them if you want to include them as normal characters in your search expression. You must also release the backslash character itself from the expression. For example, the search pattern R\\*.\* returns the result R\*3.24 but not R/2.87.

Nevertheless, bracketed expressions are an exception to the rule. Inside bracketed expressions, all special characters, including the backslash, lose their special meaning (for instance, [\*\+?{}.] matches exactly any of the characters inside the brackets).

The order of rank for the regular expressions described above is as shown in Table 9-2. For example, the regular expression abc2|3de matches either the string abc2 or the string 3de (rather than the string abc2de or abc3de) because concatenation has a higher ranking order than alternation.

Escaped characters	\ <special character:<="" th=""></special>
Bracket expression	[]
Grouping	()
Single-character duplication	* + ? {m,n}
Concatenation	
Anchoring	^ \$

Alternation	
-------------	--

Table 9-2. Ranking order (from high to low) for regular expressions.

### 9.6.3 Reports


You can select annotations by selecting their corresponding lines in the Annotation Table, available by choosing *Reports: Annotation Table*.

## 9.7 Display annotations and labels

You can change the way annotations and labels are displayed.

## 9.7.1 Annotation flag position

Sometimes you may want to move an annotation flag because you are preparing an illustration or want to see what lies underneath. To interactively change an annotation's flag position, click on one of the labels and drag the flag to the desired position (Figure 9-7).






Figure 9-7. (a) Default and (b) modified annotation flag position.

### 9.7.2 Visibility of annotations and labels

You can quickly cover entire gel images with a considerable number of annotations and labels, which are not always necessary at any given moment in the analysis. Therefore, you can hide all annotations or certain label categories in selected gels. These options are available from the *View*:

\*\*Annotations\*\* menu:

Visible categories: Sets the visibility state of the various categories on the
selected gels. To hide a category, make sure its box is unchecked. On the
other hand, select an empty check box to show the corresponding
category. Click once or twice in a grayed check box (category that takes
different visibility states in the various selected images) to hide or show the
corresponding label category in the selected gels.

- **Show All**: Makes all the annotations on the selected gels visible, including the labels that were hidden.
- *Hide All*: Hides all the annotations on the selected gels, including the cyan cross or square that remains visible when all labels, but not annotations, are hidden

When you click on an annotation that has hidden labels, all of its labels are displayed on the screen during the time it remains selected (Figure 9-8). The hidden labels disappear when the annotation is deselected.




Figure 9-8. (a) Unselected annotation with two labels. (b) Unselected annotation with 1 label hidden. (c) Unselected annotation with all labels hidden. (d) When the annotation is selected the hidden labels become visible.

#### 9.8 **Annotation table**

The Annotation Table (Figure 9-9), available by choosing **Reports**: **Annotation Table** provides specific information about annotations, and consequently labels and categories:

- Name of the image on which the annotation was created.
- SpotID, if the annotation is attached to a spot.
- Coordinates of the annotation.
- Calculated pl and MW values, if pl\_MW annotations were defined on the image, or a matched image.
- A column for each label category, with the label content in the corresponding cells.




Figure 9-9. Annotation Table.

## 9.9 Edit annotations and labels

You can add and modify annotations and labels in different ways. The **Select** tool is helpful if you want to add or modify a single or a few labels. The menu options are more adapted to the creation of a large number of annotations or labels simultaneously. Reports are useful for editing existing annotations.

# 9.9.1 Select

When double-clicking on a label while the **Select** tool is activated, the Edit Label window is displayed. Change the text in this box to modify your label content.

The **Select** tool should also be used to change the position of an annotation. In this case, simply select an annotation and drag its basis to the new location.

#### 9.9.2 Edit menu

Various options for editing labels or annotations are available from the *Edit*: *Annotations* menu. All of these features can be used to edit several labels or annotations at a time.



Two possibilities are available for deleting selected labels or annotations:

- Edit: Annotations: Delete deletes the selected annotations.
- Edit: Annotations: Labels: Delete deletes the selected labels.

To delete labels from a specific category only, first select all labels from the desired category with **Select: Annotations: By Category**.



You can change the content of selected labels (belonging to a single category) by choosing *Edit : Annotations : Labels : Delete* and entering the desired modifications in the Edit Labels window.

#### Copy / paste annotations

You can select annotations in a given gel image and copy them to the corresponding locations in other gels. This is a means of creating a set of similar annotations in a series of gels. Subsequently, you may need to adjust the annotation positions.

Select annotations on one image and choose Edit: Annotations: Labels: Copy. Then select the gels into which you want to paste the annotations and choose Edit: Annotations: Labels: Paste. Adjust the annotation positions using the Select tool.

#### Copy / paste labels

Instead of copying entire annotations, you can also copy distinct labels from one gel to selected spots or annotations in a series of gels.

Select the labels you would like to copy and choose *Edit*: *Annotations*: *Labels* : Copy. Then select spots or existing annotations into which you want to paste the labels and choose Edit: Annotations: Labels: Paste.

#### Propagate to matched

When gels have been matched, labels selected in one gel can be copied to their corresponding spots in other gels by choosing *Edit*: *Annotations*: *Labels*: **Propagate to Matched**. This is particularly useful when you have annotated one image during analysis, and want to propagate the labels to all matched gels.

#### **Duplicate labels**

You can copy selected labels to another category by choosing *Edit*: **Annotations: Labels: Duplicate.** Since the selected labels may belong to different categories, this option can be used to merge several categories into a new one. However, only one label per annotation can be duplicated at a time.

#### Annotation table 9.9.3

You can edit labels via the Annotation table, available from Reports: Annotation Table. Make sure you have displayed the desired categories for editing, by clicking the Settings icon.

#### Create labels and label categories

Click on the Add Label icon in the Annotation Table to create new labels for selected spots. During this process, you can create a new label category, which will be inserted as an extra column in the table.

#### Add or modify labels

You can directly add new labels to the appropriate cells of the Annotation Table, or edit existing ones. Double-click in a cell to start typing your label, or modifying an existing label. When finished, a single click in any cell guits the editing mode. Please note that categories using the Data Type Set are displayed in the form of check boxes in the table. These boxes show if the corresponding item belongs to the category (checked box) or not.

# 9.9.4 A Import annotations

You can import annotations from a file into open gels. Choose *Edit*: *Annotations: Import* to import annotations from an Annotation Report or a tab-delimited text file containing the required columns: SpotID, X, Y, and a column for each category to be imported.

If the Spot ID is not known, use -1 in this field (or remove this column) and the software will position the label in the corresponding X and Y positions of the gel. If X and Y positions are not known, use -1 in these fields and the software will position the label on the spot with the corresponding Spot ID.

# 10 Data integration

#### Introduction 101

Projects that were created with earlier versions of ImageMaster 2D Platinum can be imported to version 7.0. It is also possible to import images from TWAIN compatible scanners. After analysis, data can be exported as reports or to spot excision robots

#### Convert projects from earlier software versions 10.2

ImageMaster can Donvert Projects created with versions 4, 5 and 6 of ImageMaster 2D Platinum. Images analyzed with older versions must be added to a project in order to be imported into this new version and for spots, annotations and match data to be recovered. You can do this batch conversion for many projects at a time.

Choose File: Import: ImageMaster 2D Platinum or Melanie Data, and indicate the folder where some or all of your projects (.prj files) are saved. The software searches and displays all project files found in this folder, and allows you to select projects for conversion (Figure 10-1). Then you must indicate the location where the converted projects should be saved. After conversion, the projects are automatically inserted in the workspace.



Figure 10-1. Select projects from previous ImageMaster versions for Batch File Conversion

#### Acquire images from Twain compatible scanners 10.3

ImageMaster can also acquire images directly from TWAIN compatible scanners.

# 10.3.1 Select source

You must indicate the scanning source by choosing *File:Import:Twain:Select Source*. All TWAIN compatible scanners attached to your PC are automatically recognized by ImageMaster. This selection only needs to be done once (unless you want to change to a new image capture device).

# 10.3.2 Acquire

Launch the scan with *File: Import: Twain: Acquire*. The scanner software opens, giving you the opportunity to change the necessary settings, and subsequently initiate the scanning process. Once this is done, the image is saved in ImageMaster file format and added to the Image Pool.

## 10.4 Export data

#### 10.4.1 Reports

Data obtained in the analysis can be exported for use in other applications by saving the various tabular reports. Tables can be saved in text format (.txt), as a Microsoft Excel Workbook (.xls), or in XML (.xml) format. See below for more details on XML.

Graphical reports can be saved in PNG, TIFF or BMP formats.

#### 10.4.2 XML

XML stands for eXtensible Markup Language and was created as a crossplatform, software, and hardware independent tool to structure, store, and exchange information. It allows the creation of customized tags, enabling the definition, transmission, validation, and interpretation of data between applications and organizations.

XML files can be viewed in the latest versions of Web browsers such as Internet Explorer, Netscape and Mozilla Firefox. However, as XML was designed to describe data and not to display data, it does not look like a Web page. An XML document contains color-coded root and child elements. A plus (+) or minus (-) sign to the left of the elements can be clicked to expand or collapse the element structure. If you want to view the raw XML source, you must select *View: Source* from the browser menu.

XML does not use predefined tags, as is the case for HTML. Therefore, the browser does not understand the meaning of the tags and does not know how to display the XML document. Therefore, XSL (eXtensible Stylesheet Language) stylesheets must be used in addition to the XML document to transform the XML into the sort of document that is recognized by the browser. This is the case when tabular reports or the history are printed from ImageMaster. The software

uses the XSL stylesheets located in the Template folder of the ImageMaster installation directory to print attractive documents. If you are familiar with XML and XSL, you can even create personalized templates for printing. Note that XSL stylesheets can also be used to convert an XML file into another XML file.

The main interest in XML format is that external applications can easily extract necessary data. Moreover, the files can be converted to other user-defined formats.

Note: Because the XSL stylesheets are specific to the browser you use, you will find that different versions (both for printing reports and the history) are installed with the software (in the Template\Report and Template\Script folders of the installation directory). ImageMaster will therefore ask you to choose the appropriate XSL template each time you print a report or history. Look at which template works with your browser and delete the other one. Next time, the software automatically opens the remaining file and does not ask you to make a choice

#### 10.4.3 Gel and report identifiers

ImageMaster allocates a unique identifier (ID) to each project, match set, gel image, etc. This is useful to assure data consistency, allow reliable identification of objects across a computer network and enable database integration.

The IDs in ImageMaster are UUIDs (Universal Unique IDentifiers), which are 128-bit numbers that are guaranteed to be unique through combinations of hardware addresses, time stamps and random seeds. These IDs allow the ImageMaster objects to be uniquely recognized. The software can therefore detect, for example, if you delete a gel and replace it with another one with the same name.

**Note:** You must be very careful when manipulating project, match set, gel and other ImageMaster files outside the software, in order not to corrupt the data.

# 10.5 Export to spot excision robots

You can export spots to an excision robot. In addition, it may be useful to save them as part of a spot set, or annotate them. This allows you to easily select them later on to add experimental data such as mass spectrometry information.

#### 10.5.1 GE Healthcare Ettan spot picker

To use the Ettan Spot Picker, two adhesive markers should be placed on the gel before scanning. These markers are used for the calibration of the coordinates, that is, for determining the correspondence between the X and Y positions of the

analyzed gel image and the coordinates of the actual gel located on the spot picker (Figure 10-2).

Once the gel has been digitized and analyzed with ImageMaster, the software can generate a pick list. This list contains the location, in pixels, of the center of each spot you wish to pick, as well as the pixel coordinates for the centers of the two reference markers.


To export a file with spot coordinates for use by the spot picker, you first have to open the image files and perform image analysis (spot detection is mandatory). You should then annotate the reference markers on the marker spot or on the pixel in the center of the marker. The basis of an annotation on the marker spot is displayed as a small square, while a pixel annotation has a crossed basis.

**Note:** The two annotation options can be used in a single gel (one marker with a spot annotation and the other with a pixel annotation), as long as the annotations are centered on the markers.

#### To create a pick list:

- 1 Identify the two reference markers on your image.
- 2 Zoom the image to better see the left reference marker.
- 3 Check if the marker is detected as a nice, round spot. If it is not, select the spot(s) on the marker and delete by choosing *Edit: Spots: Delete*.
- 4 Click on the **Select** tool.
- 5 Double-click on the marker spot, or if such a spot is not present, on the pixel that is in the middle of the marker.
- 6 Select the **Comment** category.
- 7 Enter **IR1** as the label text.
- 8 If an annotation attached to a pixel is not in the center of the reference marker, you can move it to the appropriate position. Do this by clicking on its basis (cross) and holding down the left mouse button while dragging the annotation to its new position.
- 9 Move your gel to see the second (right) reference marker.
- 10 Repeat the procedure, but enter the label text **IR2** this time.
- 11 Once the two markers have been annotated, select spots to pick (Figure 10-2).
- 12 Choose File: Export: Spots to Picker: GE Healthcare Ettan.

13 For each gel, you will be asked to save a pick list in text or XML format (only the text file can be read by the Ettan Spot Picker).



**Figure 10-2.** Reference markers IR1 and IR2. IR1 is attached to a pixel (cross basis). IR2 is attached to a spot. Both options can be used in a single gel. Spots to be picked are selected (highlighted in green).

## 10.5.2 Spot pickers from other manufacturers

The software can also export spots to spot pickers from Bruker Daltonics, Genetix and Genomic Solutions.

- 10 Data integration10.5 Export to spot excision robots

# 11 Undo, redo and history

#### 111 Undo, redo

The software allows you to cancel any unwanted modifications to the images or layout of the current sheet by choosing Edit: Undo in the menu (Figure 11-1).




Figure 11-1. Undo window.

Each action in the Undo/Redo list is registered along with the time at which it was carried out.

By default, all actions are displayed. If you specifically want to track permanent changes that have been made to the image data (essentially the commands available under the Edit menu), choose Only Edit from the Filter drop-down list.

All operations performed can be reversed except for modifications to the Workspace and functionalities linked to opening and saving files. This includes image rotation, flipping, cropping, and the inversion of gray levels.

## To undo specific actions:

- Choose Edit: Undo.
- Select a prior action to be undone.
- Click **OK**. The selected action and all following actions are undone automatically.

If you are not satisfied with your latest undo or you canceled too many operations, you can restore the actions.

### To redo specific actions:

- 1 Choose *Edit: Redo*.
- 2 Select the action to be redone.
- 3 Click **OK**. The selected action and any previous actions are redone.

## 11.2 History

You can display a history of the modifications that have been carried out during the present work session on the images or layout of the current sheet. Choose *Edit: History: Show* to open the History window (Figure 11-2).




Figure 11-2. History window.

Please note that the action names displayed in the History window are the same as those used in Undo/Redo. As for Undo/Redo, you can choose *Only Edit* from the *Filter* drop-down list to restrict the displayed actions to those that bring about permanent modifications to the image data.

Some actions are preceded by a + or - node allowing the item to be expanded or collapsed. Once an action is expanded (by clicking on the plus sign), you see parameters and values that further describe the action.

Modifications to the Workspace and functionalities linked to opening and saving files (including image rotation, flipping, cropping, and inversion of gray levels) are not included in the History.

#### Insert marker

You can place a marker in the History by choosing *Edit: History: Insert Marker*.

#### Clear

Clear the list of actions in the History by choosing *Edit: History: Clear*.

#### Refresh

The content of the History is not automatically updated when you work on the images while the History window is still open. Click the *Refresh* icon to update the list.

## Save, print, copy

By using the options in the Save icon drop-down list in the History window, the content of the History can be:

- Saved in an XML type file with the extension .hst.
- Printed. It will first display in your Web browser using the XSL stylesheet located in the Template\Script folder of the ImageMaster installation directory. Use the print option in your browser to get a paper copy.
- Copied to the clipboard.

11 Undo, redo and history 11.2 History

# Appendix A Shortcuts

# A.1 Shortcut keys

Several menu commands can be activated by keyboard shortcuts. They are indicated to the right of the corresponding command in the menus. Please note the logic behind the key combinations:

Ctrl for gels.

**Shift** for spots.

**Alt** for annotations.

Ctrl + Shift for matches

Some exceptions do exist. The most important ones are the following two shortcuts used for undoing and redoing actions carried out on a gel.

Shortcut	Menu Command
Ctrl+Z	Edit : Undo
Ctrl +Shift+Z	Edit : Redo

# A.2 Tool shortcuts

Shortcut	Tool
Ctrl+1	Move
Ctrl+2	Zoom
Ctrl+3	Region
Ctrl+4	Selection
Ctrl+5	Measure
Ctrl+6	Landmark

# A.3 Gel shortcuts

Shortcut	Menu Command
Ctrl+F	View : Sheet : Navigation : Switch
Ctrl+I	View : Global : Grid Lines : Show
Ctrl+J	View : Global : Grid Lines : Edit
Ctrl+O	File : Open
Ctrl+P	View : Global : Show Profile
Ctrl+S	File : Save
Ctrl+W	File : Close : Sheet
Ctrl+ <down:< td=""><td>View : Gels : Navigation : Move : Down</td></down:<>	View : Gels : Navigation : Move : Down
Ctrl+ <left:< td=""><td>View : Gels : Navigation : Move : Left</td></left:<>	View : Gels : Navigation : Move : Left
Ctrl+ <right:< td=""><td>View : Gels : Navigation : Move : Right</td></right:<>	View : Gels : Navigation : Move : Right
Ctrl+ <up:< td=""><td>View : Gels : Navigation : Move : Up</td></up:<>	View : Gels : Navigation : Move : Up
F2	View : Show All
F3	View : Hide All
F5	Select : Unselect All
Page Down	View : Sheet : Navigation : Previous Page
Page Up	View : Sheet : Navigation : Next Page
Shift+ <down:< td=""><td>View : Gels : Navigation : Zoom : Out</td></down:<>	View : Gels : Navigation : Zoom : Out
Shift+ <up:< td=""><td>View : Gels : Navigation : Zoom : In</td></up:<>	View : Gels : Navigation : Zoom : In

# A.4 Spot shortcuts

Shortcut	Menu Command
Shift+1	View : Spots : Crossed
Shift+2	View : Spots : Outlined
Shift+3	View : Spots : None
Shift+A	Select : Spots : All
Shift+E	Edit : Spots : Edit Enabled
Shift+N	Select : Spots : Inverse Selection
Shift+X	Edit : Spots : Delete

# A.5 Annotation shortcuts

Shortcut	Menu Command
Alt+A	Select : Annotations : All
Alt+C	Edit : Annotations : Labels : Copy
Alt+D	Edit : Annotations : Delete
Alt+E	Edit : Annotations : Labels : Edit
Alt+F	Edit : Annotations : Add
Alt+H	View : Annotations : Hide All
Alt+J	View : Annotations : Visible Categories
Alt+L	Edit : Annotations : Link with Spot
Alt+U	Edit : Annotations : Unlink from Spot
Alt+V	Edit : Annotations : Labels : Paste
Alt+X	Edit : Annotations : Labels : Delete
Alt+Y	View : Annotations : Show All

# A.6 Match shortcuts

Shortcut	Menu Command
Ctrl+Shift+A	Select : Matches : All
Ctrl+Shift+G	Edit : Matches : Add Match

# A Shortcuts A.6 Match shortcuts

Shortcut	Menu Command
Ctrl+Shift+J	View : Matches : Show ID
Ctrl+Shift+K	View : Matches : Hide All ID
Ctrl+Shift+N	Select : Matches : Inverse Selection
Ctrl+Shift+U	Edit : Matches : Delete Match
Ctrl+Shift+Y	View : Matches : Show Vectors

#### Appendix B References

#### **B.1 Software**

Appel RD, Hochstrasser DF, Roch C, Funk M, Muller AF and Pellegrini C (1988) Automatic classification of two-dimensional gel electrophoresis pictures by heuristic clustering analysis: A step toward machine learning. Electrophoresis 9: 136-142.

Appel RD, Hochstrasser DF, Funk M, Vargas JR, Pellegrini C, Muller AF and Scherrer J-R (1991) The MELANIE project - from a biopsy to automatic protein map interpretation by computer. Electrophoresis 12: 722-735.

Appel RD, Palagi PM, Walther D, Vargas JR, Sanchez J-C, Ravier F, Pasquali C and Hochstrasser DF (1997) Melanie II - a third generation software package for analysis of two-dimensional electrophoresis images: I. Features and user interface. Electrophoresis 18: 2724-2734.

Appel RD, Vargas JR, Palagi PM, Walther D and Hochstrasser DF (1997) Melanie II - a third generation software package for analysis of two-dimensional electrophoresis images: II. Algorithms. Electrophoresis 18: 2735-2748.

Appel RD and Hochstrasser DF (1998) Computer analysis of 2-D images. In: Link AJ (ed) Methods in Molecular Biology, Vol 112: 2-D Protocols for Proteome Analysis, pp 363-381. Totowa NJ: Humana Press.

Miller MJ, Olson AD and Thorgeirsson SS (1984) Computer analysis of twodimensional gels: automatic matching. Electrophoresis 5: 297-303.

Pun T, Hochstrasser DF, Appel RD, Funk M, Villars-Augsburger V and Pellegrini C (1988) Computerized classification of two-dimensional gel electrophoretograms by correspondence analysis and ascendant hierarchical clustering. Applied and Theoretical Electrophoresis 1: 3-9.

Vargas RJ (1996) Two-dimensional gel electrophoresis computer analysis systems: from image acquisition to protein identification. Ph.D. thesis, Faculty of Science, Geneva University.

Wilkins MC, Hochstrasser DF, Sanchez J-C, Bairoch A and Appel RD (1996) Integrating two-dimensional gel databases using the Melanie II software. Trends in Biochemical Sciences 21: 496-497.

#### **B.2** Statistical methods

Armitage P and Berry G (1987) Statistical methods and medical research. Oxford, London: Blackwell Scientific Publications.

Kim JO and Mueller CW (1978) *Introduction to factor analysis: What is and how to do it?* Newbury Park: Sage Publications.

Tabachnick B and Fidell LS (1996) *Using multivariate statistics* (3rd edition). New York: Harper Collins College Publishers.

# B.3 Further reading

Appel RD, Bairoch A, Sanchez J-C, Vargas JR, Golaz O, Pasquali C and Hochstrasser DH (1996) Federated 2-DE database: a simple means of publishing 2-DE data. *Electrophoresis* 17: 540-546.

Appel RD, Sanchez J-C, Bairoch A, Golaz O, Ravier F, Pasquali C, Hughes GJ and Hochstrasser DF (1996) The SWISS-2DPAGE database of two-dimensional polyacrylamide gel electrophoresis. *Nucleic Acids Research* 22: 3581-3582.

Binz PA, Mueller M, Walther D, Bienvenut WV, Gras R, Hoogland C, Bouchet G, Gasteiger E, Fabbretti R, Gay S, Palagi P, Wilkins MR, Rouge V, Tonella L, Paesano S, Rossellat G, Karmime A, Bairoch A, Sanchez JC, Appel RD and Hochstrasser DF (1999) A molecular scanner to automate proteomic research and to display proteome images. *Analytical Chemistry* 71: 4981-4988.

Binz PA, Wilkins MR, Gasteiger E, Bairoch A, Appel RD and Hochstrasser DF (1999) Internet resources for protein identification and characterization. In: Kellner R, Lottspeich F, Meyer HE (eds) *Microcharacterization of Proteins*, 2nd ed., pp. 277-300. Weinheim: Wiley-VCH.

Hoogland C, Baujard V, Sanchez J-C, Hochstrasser DF and Appel RD (1997) Make2ddb: a simple package to set up a 2-DE database on the WWW. *Electrophoresis* 18: 2755-2758.

Hoogland C, Sanchez J-C, Bairoch A, Hochstrasser DF and Appel RD (1999) The SWISS-2DPAGE database: what has changed during the last year. *Nucleic Acids Research* 27: 289-291.

Link AJ (ed) (1998) Methods in molecular biology, Vol 112: 2-D Protocols for Proteome Analysis. Totowa NJ: Humana Press.

Lopez MF (2000) Better approaches to finding the needle in a haystack: optimizing proteome analysis through automation. *Electrophoresis* 21: 1082-1093

Sanchez J-C, Wilkins M, Appel RD and Hochstrasser DF (1997) Identifying proteins for proteome studies. In: Creighton ET (ed) *Protein Function: a practical approach*, 2nd ed., pp. 1–27. IRL Press.

Wilkins MR, Williams KL, Appel RD and Hochstrasser DF (eds) (1997) *Proteome* research: new frontiers in functional genomics. Berlin Heidelberg: Springer Verlag.

Unlu M, Morgan ME, and Minden JS (1997) Difference gel electrophoresis: a single gel method for detecting changes in protein extracts. *Electrophoresis* 18: 2071-2077.

Tonge R, Shaw J, Middleton B, Rowlinson R, Rayner S, Young J, Pognan F, Hawkins E, Currie I and Davison M (2001) Validation and development of fluorescence two-dimensional differential gel electrophoresis proteomics technology. Proteomics 1: 377-396.

- B References
- B.3 Further reading

## Index

**Symbols** 

```
% 74
%vol 89
.bmp 85
.gda 64
.mda 64
.png 85
.prj 64
.tif 85
Numerics
3D view 75
     animate 78
     auto rotate 77
     color palette 78
     display options 78
     next/previous 78
     rotation 77
     set reference 78
     show 77
     spot shape 78
     stack 78
     tools 77
     translation 77
     transparency settings 78
     transparent mode 78
     visibility 78
     zoom/contrast 77
access code 11
adaptive gradations 124, 144
add
     files to project 57, 65
     labels 163
     matches 113
     spots 96
adjust
     contrast 71
     region contrast 75
algorithm 89
     matching 105
align images 80
analysis
     class 117
     gel 117
analyze
     classes 28, 117, 135
     gels 117, 120
animate 78
```

```
annotations 147
    copy 163
    create 147
    create categories 163
    delete 162
    display 160
    flag position 160
    hide 160
    import 164
    MW and pl 100
    options 34
    paste 163
    report 160, 161, 163
    select 156
    shortcuts 177
    show 160
    table 28
    tool 156
apply
    calibration 48
    region to all gels 70
    zoom to all gels 67
area 89, 134
auto rotate 77
automatic matching 110
autonumbering 150
backup 65
BMP 85
bookmarks
    create 71
    delete 71
    load 71
Bruker Proteineer SP spot picker 167
С
calibration 37
    apply 48
    control 48
    copy to clipboard 47
    create 44
    fitting table 47
    images 44
    MW and pl 100
    open 47
    print 47
    remove 48
    save 47
    table 47
    tablet file 44
categories 163
center 141
central tendency 118, 123, 141
```

```
CGI script 151
choose matchset 114
classes
    analysis 135
    analysis histograms 143
    analysis table 141
    create 63
    folder 55
co-detection 91
coef. variation 123
coefficient variation 123
collect eLicense 12
color palette 78
colors 75
combine spot sets 94
command prompt 11
common labels 157
compare images 80
composite spots 98
concurrent license 11
contact
    information 14
    us 14
contextual menus 27
contrast 74, 77
control calibration 48
сору
    annotations 163
    formula to clipboard 120
    labels 163
    to clipboard 29
Corr 120
correlation coefficient 120
create
    annotation categories 163
    annotations 147
    calibration 44
    gel descriptions category 49
    label categories 149
    labels 147
    match hierarchy 58
    pick list 168
    projects 56
    specific link 154
    spot set 124, 142
    spot sets 94
    spots 97
    toolbar 21
crop 41
    area 41
    export area 41
    import area 41
current sheet 27
cursor information 79, 102
```

```
customize
    reports 30
    toolbars 21
data
    analysis 117
    type 150
database options 34
DD 44
DeCyder 2D algorithm 87
default mapping 73
define landmarks 107
delete
    annotations 162
    labels 162
    matches 113
    spot sets 94
    spots 97
describe images 49
descriptions 49
desktop shortcut 13
DIGE 9, 87
    co-detect spots 91
    exclude spots 93
    file naming convention 37
    histogram 133
    spot detection parameter 92
    spot quantification 92
    value 117
disable spots 96
dispersion 119, 123, 141
display
    annotations 160
    calibration information 44
    gel descriptions 49
    labels 160
    match hierarchy 106
    matches 111
    options 34, 78
    projects 63
    properties 151
    spots 97
    zone 17, 22
displayed
    items 129
    value 142, 144
dockable windows 22
    floating 33
    minimize 33
    pinned 33
    tabbed groups 33
    un-pinned 33
double detection 91
```

```
dual color 81
duplicate labels 163
Ε
EDF 140
edit 17
     annotations 163
     enabled 96
     gel descriptions 49
     grid lines 83
     matches 113
     region 70
     spots 97
     table cells 32
eLicense 10
     collect 12
     file 12
     place 12
     server 11
     test 13
empirical distribution
     function 140
     plot 140
enabled spots 96
ExPASy 154
export 85, 166
     crop area 41
     match set 62
     spot coordinate file 167
     spot excision robot 167
external engine 151
external protein databases
     query 153
     set 152
factor analysis 127, 135
     comments 132
     interpret 130
     table 124, 129
factor projection
     plot 129
     table 129
files 17
    format 37
    link 154
     names 37
first time software launched 56
fitting table 47
flag position 160
FLEXIm license finder 13
flip 41
floating 33
floating license 11
```

```
free 22, 23
    rotate 41
further reading 180
G
gap 135, 141
Gaussian distribution 137
GE Healthcare Ettan spot picker 167
gels 67
     analysis histograms 124
     analysis table 123
     description options 34
     descriptions 49
     export 85
    ID 167
     print 85
     reports 84
     save 85
     shortcuts 176
    table 28, 49, 84
general settings 117
good images 37
graphical user interface 17
gray level mapping 72
gray+saturation 75
arid lines
     display 83
     edit 83
    show 83
grow spots 97
GUI 17
half-range size 119
help 17
hide
    annotations 160
     ID 112
    labels 160
histograms
    class analysis 143
     DIGE 133
     gel analysis 124
history 172
    clear 172
     copy 172
    insert marker 172
     print 172
     refresh 172
     save 172
http
     link 154
     queries 151
```

```
image pool 22, 39
    display 39
    hide 39
    properties 39
    remove 39
    sheet 39
images 67
    calibrate 44
    calibration 37
    compare 80
    depth 37
    describe 49
    editing 37
    export 85
    file format 37
    file names 37
    good 37
    manipulate 67
    open 38
    preview 38
    print 85, 86
    process 41
    resolution 37
    save 85
    selection 24
    to clipboard 86
    to file 86
import 164, 166
    annotations 164
    crop area 41
    labels 164
    match set 62
install
    eLicense server 11
    software 10
installation cd-rom 10
installer 10
intensity 74, 89
interface 17
interpret a factor analysis 130
invert 75
    gray levels 41
ipconfig 11
is unique 150
keyboard shortcuts 18
Kolmogorov-Smirnov 140
labels 164
    add 163
    copy 163
```

```
create 147
    create categories 149
    delete 162
    display 160
    duplicate 163
    hide 160
    modify 163
    paste 163
    report 160, 161
    select 156
    show 160
landmarks
    define 107
    tool 107
launch
    software 13
    software first time 56
license
    concurrent 11
    floating 11
    information 14
    machine 10
    node-locked 10
LMTOOLS 11
load spots 96
log volume ratio 133
machine license 10
MAD 119
manage projects 65
manipulate images 67
Mann-Whitney 138
manual editing 97
mapping
    default 73
    gray levels 72
match
    count 110, 123, 141
    ID 112
match sets
    existing 61
    export 62
    import 62
    merge 59
matches
    add 113
    create hierarchy 58
    delete 113
    display 111
    display hierarchy 106
    displayed 129
    edit 113
    gels tool 110
```

```
hierarchy 53, 105
    multiple 113
    reference 105
    report 114
    review 113
    shortcuts 177
    statistics table 114
    vectors 111
matching 105
    automatic 110
    reference 59
max 123, 141
    intensity 134
    slope 134
    volume 134
mean 118
    absolute deviation 119
    squared deviation 119
measure 70, 133
measures overlapping 135
median 118
menu
    bar 17
    select 157
merge
    match set 59
    spots 97
midrange 118
min area 88
minimize 33
modify labels 163
move
    alternatives 67
    double-click 67
    projects 63
    tool 67
MSD 119
multiple matches 113
MW 100
navigator 25
    contextual menus 27
    icons 27
next in selection 29
next/previous 78
node-locked license 10
non-DIGE 9, 87
    detect spots 87
    spot detection parameters 88
    spot quantification 89
normalization 124
normalized 124, 126, 135, 141
number 150
```

```
0
OD 44
one
    column 22, 23
    row 22, 23
one-way ANOVA 137
open images 38
options 34, 117
    print 86
overlapping
    measures 135
    spots 82
overview option 68
page setup 86
panes 22
    close 23
    free 23
    layout 23
    one column 23
    one row 23
    selection 23
    stacked 23
    tiled 23
paste
    annotations 163
    labels 163
physical address 11
pl 100
pick list 168
pinned 33
pixel distances 70
place eLicense 12
plots
    empirical distribution 140
    scatter 120
PNG 85
populations 55
preview images 38
previous in selection 29
print 29, 85
    images 86
    options 86
    sheets 86
process images 41
product information 14
profile 78
projectname.prj 64
projects
    add 65
    backup 65
    convert data to new version 165
    create 56
```

```
display 63
    folder 64
    import older versions 165
    manage 65
    move 63
    properties 63
    remove 63, 65
    restore 65
    save 64
propagate spots 99
properties 63
purchase
    software 9
quantification
    options 34
    value 117
query external protein databases 153
range ratio 123
ratio 124, 126, 135, 141
raw images 64
recycle bin 64
redo 171
reference 118
    match 59
    set 80
    sheet 80
references 180
    software 179
    statistical methods 179
region
    adjust contrast 75
    apply to all gels 70
    edit 70
    tool 69
regular expressions 157
relative 124, 126
remove
    calibration 48
    projects 63
    spots 96
rename category 163
reports 17, 163
    3D view 28, 75
    analyze classes 28
    analyze gels 28
    annotation 161
    annotations 160
    columns 31
    customize 30
    edit labels 163
```

```
export 166
    factor projection 127
    gel 84
    gel analysis histograms 124
    ID 167
    import 166
    labels 160, 161
    match 114
    spot 102
    toolbar 29
resolution 37
restore 65
review matches 113
rotate 41
rotation 77
saliency 88
save 29, 85
    projects 64
    spots 96
scatter
    plots 120
    table 120
scrollbars 24
select 17, 94
    all 157
    annotations 156
    by category 157
    by content 157
    by value 29, 93, 124, 142
    common labels 157
    enabled spots 96
    labels 156
    matches 110
    menu 157
    spot sets 96
    tool 110, 156
selected region 75
separability 123
set 150
    external protein databases 152
    reference 59, 78, 80
    spots 96
settings 29, 124, 142, 144
setup wizard 10
sheets 22
    align images 80
    close 22
    current 27
    free 22
    image pool 39
    layout 22
    one column 22
```

```
one row 22
    print 86
    reference 80
    selection 22
    stacked 22
    tiled 22
    to clipboard 86
    to file 86
shortcuts 175
    annotations 177
    gels 176
    matches 177
    redo 175
    spots 177
    tools 176
    undo 175
show
    annotations 160
    default 3D view 77
    dual color 81
    ID 112
    labels 124, 160
    overview 68
    profile 78
    vectors 111
shrink spots 97
signal intensity 71
single detection 91
slider 120, 124, 133, 144
smooth 88
software references 179
sorted values 124, 144
specific links
    create 154
    file 154
    http 154
    text 154
specify classes 135
spot
    coordinate file 167
    excision robots 167
    ID 102
spot pickers 167
    Bruker Proteineer SP 167
    GE Healthcare Ettan 167
spots 87
    add 96
    co-detection parameter 92
    color 97
    composite 98
    create 97
    delete 97
    detection algorithm 87
    DIGE co-detection 91
```

```
disable 96
    display 97
    edit 97
    enabled 96
    frequency 133
    grow 97
    load 96
    manual editing 97
    merge 97
    non-DIGE detection 87
    overlapped 82
    propagate 99
    remove 96
    report 102
    save 96
    select 94
    sets 94
    shape 78, 97
    shortcuts 177
    shrink 97
    table 28, 102
stack 78
stacked 22, 23
statistical methods references 179
statistics 117, 124, 142, 144
    tests 137
status bar 17, 22
step tablets 44
student's t-test 137
suspend synchronization 29
switch order 24
system
    information 14
    requirements 9
tabbed groups 33
tables
    annotation 28
    calibration 47
    class analysis 141
    edit cells 32
    factor analysis 127
    fitting 47
    gel 28, 84
    gel analysis 123
    match statistics 114
    scatter 120, 122
    spot 28, 102
tests
    ANOVA 137
    Kolmogorov-Smirnov 137, 140
    Mann-Whitney 137, 138
    student's t 137
```

```
Wilcoxon 138
text 150
    link 154
TIFF 85
tiled 22, 23
toolbar 17, 20
    create 21
    format 21
    position 21
tools 17
    add files to project 57
    adjust contrast 71
    annotation 156
    cursor information 79
    landmark 107
    match gels 110
    measure 70
    move 67
    profile 78
    region 69
    select 94, 110, 156
    shortcuts 176
    show dual color 81
    zoom 67
translation 77
transparency settings 78
transparent mode 78
trimmed
    mean 118
    midrange 118
triple detection 91
t-test 137
Twain compatible scanners 165
undo 171
unit 74
un-pinned 33
update gels 164
UUID 167
value 117, 124, 126
vectors 111
view 17
    signal intensity 71
visibility 78
vol 89, 92
    ratio 92
whole image 75
Wilcoxon 138
workspace 22
```

```
class folder 63
classes 63
match folder 58
match hierarchy 58
navigator 25
projects 56
toolbar 25

X

XML format 166

Z

zoom 67, 77
alternatives 68
apply to all gels 67
overview option 68
```

For contact information for your local office, please visit www.gelifesciences.com/contact

GE Healthcare Bio-Sciences AB Björkgatan 30 751 84 Uppsala Sweden

www.gelifesciences.com

 $\ensuremath{\mathsf{GE}},$  imagination at work and  $\ensuremath{\mathsf{GE}}$  monogram are trademarks of General Electric Company.

Cy, CyDye, DeCyder, Ettan, ImageMaster, Immobiline, LabScan and ImageScanner are trademarks of GE Healthcare companies.

Any use of this software is subject to GE Healthcare Standard Software End-User License Agreement for Life Sciences Software Products. A copy of this Standard Software End-User License Agreement is available on request.

This version of ImageMaster has been developed by the Swiss Institute of Bioinformatics in collaboration with GeneBio and GE Healthcare. All intellectual property rights on the ImageMaster Platinum Software belong to the Swiss Institute of Bioinformatics.

 $\hbox{@ 2008\,Swiss}$  Institute of Bioinformatics - All rights reserved. Reproduced by permission of the owner.

Swiss Institute of Bioinformatics CMU-Rue Michel-Servet 1, CH-1211 Geneva, Switzerland

ImageMaster 2D Platinum uses the DeCyder co-detection algorithm. © 2005 General Electric Company – All rights reserved.

ImageMaster 2D Platinum uses the TIFF library. © 1988-1999 Sam Leffler and 1991-1999 Silicon Graphics, Inc - All rights reserved.

ImageMaster 2D Platinum uses software developed by the Apache Software Foundation (www.apache.org).
© 1999-2003 The Apache Software Foundation - All rights reserved.

All third party trademarks are the property of their respective owners.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

GE Healthcare UK Ltd Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK

GE Healthcare Bio-Sciences Corp 800 Centennial Avenue, P.O. Box 1327, Piscataway, NJ 08855-1327, USA

GE Healthcare Europe GmbH Munzinger Strasse 5, D-79111 Freiburg, Germany

GE Healthcare Bio-Sciences KK Sanken Bldg. 3-25-1, Hyakunincho, Shinjuku-ku, Tokyo 169-0073, Japan

